




Cephalosporin C

Antibiotics



C₁₆H₂₁N₃O₈S
 MW 415.4
 Zinc salt, dihydrat

A Commitment to Diagnostics

Order Code BAG2-C3

Generic Names None. Cephalosporin C is the basic structure of many antibiotics such as Cefazolin and Cefamandole.

Instruction for use Use a fresh allergen vial each time a new stimulation is performed
 Dissolve the allergen prior to use as follows:
 Add 250 µl of Stimulation Buffer to vial, vortex briefly
 Use allergen solution according to the protocol

Recommendation Additional to the normal concentration of the allergen given, we recommend to stimulate the cells with 1:5 diluted allergen solution
 (Mix 40 µl of allergen solution with 160 µl of stimulation buffer)

n Stimulation Up to 4 stimulations could be run with 1 vial

Storage Store unopened allergens at 2-8°C until exp. date

	CAST®-2000	Flow-CAST®	Flow2 CAST®
Content per vial			
25 µg			
Conc. in vial			
100 µg/ml			
Conc. in stimulation	20 µg/ml	50 µg/ml	22.7 µg/ml
Cut-off	≥ 40 pg/ml	≥5% CD63 SI* (≥2)	≥5% CD63 SI* (≥2)

CAST®-Allergens are intended for in vitro leucocyte stimulation in the CAST®-Assays
 *SI (Stimulation Index) is defined as the gross stimulation divided by the negative control value



BÜHLMANN Laboratories AG
 Baselstrasse 55
 CH-4124 Schönenbuch/Basel
 Switzerland

Phone +41 61 487 12 12
 Fax orders +41 61 487 12 99
 info@buhlmannlabs.ch
 www.buhlmannlabs.ch



ADS014-02E_C3

CAST® is a registered trademark of BÜHLMANN.