# Flow CAST<sup>®</sup> highsens

**Basophil Activation Test (BAT)  
Flow Cytometry**

FK-HSBAT 100 tests

Revision date: 2011-08-17

## ENGLISH

### INTENDED USE

The **Flow CAST® highsens** kit is a basophil activation test (BAT) which is intended to be used for the *in vitro* detection of immediate type allergic reactions and hypersensitivities against suspected allergens.

The test has been designed to determine the relative amount of CD63 and CD203c positive basophils from the total number of basophils (CCR3 positive) in whole blood by flow cytometry upon antigen stimulation.

### PRINCIPLE OF THE ASSAY

Stimulation Buffer and allergen are added to EDTA whole blood from suspected allergic/hypersensitive patients. The allergen mimics the *in vivo* reaction where specific IgE molecules bound to the cellular surface are bridged by the culprit allergen. This activates an intracellular signaling cascade leading to basophil activation. As a consequence, intracellular compounds bearing the trans-membrane protein CD63 are fused with the cellular membrane and therefore exposed to the extracellular matrix. In addition, the surface marker CD203c is up-regulated.

A highly specific monoclonal antibody binding to the high affinity IgE binding receptor (FcεRI) or the unspecific cell activator fMLP are used as positive controls.

During the stimulation reaction, Staining Reagent is added. The Staining Reagent contains a mixture of monoclonal antibodies to human CD63 and CD203c both labeled with the tandem dye PE-DY647, and a monoclonal antibody to CCR3 (CD193) labeled with phycoerythrin (PE). *In vivo*, CCR3 is constitutively expressed on eosinophils and basophils (Ref. 6,7,10) and therefore it is used as selection marker for basophilic cells.

Erythrocytes are removed by a lysing reaction and after a short centrifugation step the remaining cells are resuspended in Wash Buffer and analyzed by flow cytometry (*cf.* FLOW CYTOMETRIC DATA ACQUISITION on page 4).

### REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
<b>Stimulation Buffer</b> contains calcium and heparin	1 vial lyoph.	B-CCR-CSB	Reconstitute with 50 ml of water <sup>1)</sup>
<b>Stimulation Control</b> anti-FcεRI mAb	2 vials lyoph.	B-CCR-STCON	Reconstitute with 1.5 ml of B-CCR-CSB
<b>Stimulation Control</b> fMLP <sup>2)</sup>	2 vials lyoph.	B-CCR-FMLP	Reconstitute with 1.5ml of B-CCR-CSB
<b>Staining Reagent</b> Mix of anti-CD63-PE-DY647, anti-CD203c-PE-DY647 und anti-CCR3-PE mAb	1 vial 4.4 ml	B-HSBAT-SR	Ready to use
<b>Lysing Reagent</b> 10x concentrated	1 vial 25 ml	B-CCR-LYR	Dilute with 225 ml of deionized water <sup>1)</sup>
<b>Wash Buffer</b>	1 vial 100 ml	B-CCR-WB	Ready to use

Table 1

<sup>1)</sup> For required water quality, see Chapter Procedural Notes

<sup>2)</sup> N-formyl-methionyl-leucyl-phenylalanine

### STORAGE AND SHELF LIFE OF REAGENTS

Unopened reagents	
Store at 2-8°C. Do not use past kit expiration date.	
Opened / reconstituted reagents	
Stimulation Buffer	Stable at -20°C for 6 months. Store in aliquots, if repeated use is expected.
Stimulation Control	
Stimulation Control fMLP	
Lysing Reagent	Stable at 2-8°C for 6 months.
Staining Reagent	Stable at 2-8°C until expiration date.
Wash Buffer	

Table 2

### ALLERGENS REQUIRED BUT NOT PROVIDED

Validated Allergens for the analysis in CAST® assays are offered by BÜHLMANN. Refer to the Allergen list on the webpage to obtain the respective order codes ([www.buhmannlabs.ch](http://www.buhmannlabs.ch)).

– **Protein Allergens** are shipped as concentrated liquids (1µl/vial) and must be diluted before use.

– **Drug and Chemical Allergens** are shipped lyophilized and must be reconstituted before use.

Refer to the BÜHLMANN Allergen Booklet and **Allergen Data Sheets** available on the website [www.buhmannlabs.ch](http://www.buhmannlabs.ch).

### ALLERGENS FROM OTHER SOURCES

Allergens from other sources might be used in the Flow CAST® assay with the following limitations:

- No matrix-bound allergens (solid or liquid phase).
- No allergen preparations containing cytotoxic compounds (stabilizers, preservatives) such as glycerol, phenol, sodium azide or merthiolate (thimerosal).

If you need assistance with the preparation allergens not provided by BÜHLMANN, please contact us or ask your local distributor.

### WARNINGS AND PRECAUTIONS

#### Reagents Containing Human Source Material:

Stimulation Buffer (B-CCR-CSB) contains components of human origin. Each serum used in the preparation was tested by an approved method and found negative for HBV surface antigen, HCV and HIV1/2 antibodies. Although those methods are highly accurate there is no guarantee that this material cannot transfer Hepatitis or AIDS. *Therefore all patient specimens and kit components should be handled as if capable of transmitting infections and reasonable precautions should be taken.*

### MATERIALS REQUIRED BUT NOT PROVIDED

- K-EDTA venipuncture tubes.
- Centrifuge at 500 x g.
- Disposable, pyrogen-free polystyrene test tubes and appropriate test tube racks for the stimulation  
NOTE: Polystyrene tubes should fit with the Flow Cytometer used (e.g. 12 x 75 mm FALCON tubes from Becton Dickinson; order code: 352052).
- Vortex Mixer.
- Precision pipettes with disposable, pyrogen-free tips: 10-100 µl, 100-1000 µl, 1-5 ml adjustable pipette and a 40-200 µl adjustable dispenser.

- Cylinder (50 ml) for preparing the Stimulation Buffer.
- Sterile, ultrapure and apyrogenic water for preparing the cell stimulation reagents (cf. Chapter Procedural Notes).
- Water bath set at 37°C.
- Distilled or deionized water as well as beaker or cylinder for the preparation of Lysing Reagent.
- Bottle-top dispensers for Lysing Reagent and Wash Buffer, respectively.
- Flow Cytometer instrument with 488 nm (blue) or 532 nm (green) excitation wavelength, equipped to detect Forward Scatter (FSC), Side Scatter (SSC) and the two fluorochromes PE and PE-DY647 (analogue to PC5) including appropriate software (cf. chapter FLOW CYTOMETRIC DATA ACQUISITION).

### SPECIMEN COLLECTION AND STORAGE

Collect sufficient blood in **K-EDTA venipuncture tubes**. Fill the venipuncture tubes up to the specified volume. In tubes filled < 75 %, the EDTA concentration in the sample is higher and thus false negative results may be obtained.

1 ml of whole blood is sufficient for about 9 tests.

Perform the cell stimulation immediately or store the blood sample refrigerated (2-8°C) for up to 48 hours. In order to be able to detect responses to drugs, use blood samples only for up to 24 hours after collection. **Do not centrifuge or freeze blood samples.**

### PROCEDURAL NOTES

- **RECOMMENDED WATER QUALITY TO BE USED FOR FLOW CAST<sup>®</sup> highsens.** The use of sterile, ultrapure and apyrogenic water for reconstituting Stimulation Buffer (B-CCR-CSB) is essential for good and reproducible basophil stimulation. The following sources of water may be used: Cell culture grade water, infusion grade water or deionized, double distilled water that is filtrated in a periodically sanitized 10 kDa ultra filter.

The Lysing Reagent (B-CCR-LYR) must be reconstituted with deionized, double distilled water or the same water quality that is used for the cell stimulation reagents.

- **PRECAUTIONS TO AVOID ALLERGEN CONTAMINATION DURING CELL STIMULATION.** Aeroallergens in the laboratory may contaminate open blood samples and cell suspensions from patients potentially causing an elevated background release. Therefore, care must be taken to cover blood samples and cell stimulation tubes. Avoid dust mites, pollinating plants, latex gloves or equipment potentially containing latex and open windows in the laboratory where the cell stimulation is performed. Therefore, we recommend carrying out the cell preparation and stimulation steps in a laminar flow hood.

### ASSAY PROCEDURE

1. Mix the anti-coagulated blood sample by inverting the venipuncture tube several times.
2. Prepare fresh and pyrogen-free 3.5 ml polystyrene tubes suitable for Flow Cytometry measurements.
3. For each patient, label the tubes e.g.:  
PB = patient background  
PC1 = stimulation control with anti-FcεRI Ab  
PC2 = stimulation control with fMLP  
A1-1 for allergen 1 with dilution 1  
A1-2 for allergen 1 with dilution 2  
etc.

### Stimulation and Staining

4. Pipet 100 µl of the following “stimulators” into each tube  
PB tube: 100 µl of **Stimulation Buffer (background)**  
PC1 tube: 100 µl of **Stimulation Control**  
PC2 tube: 100 µl of **Stimulation Control fMLP**  
Ax-y tube: 100 µl of **Allergen**

Repeat the sequence for each patient.

5. Add 200 µl of Stimulation Buffer to each tube.
6. Add 100 µl of patient’s whole blood to each tube. Be sure that the side wall and top of the tube are free of blood.
7. Mix gently.
8. Add 40 µl of Staining Reagent to each tube.
9. Mix gently but thoroughly, cover the tubes and incubate for 15 minutes at 37°C in a **water bath**. (using an incubator will take about 10 minutes longer incubation time due to less efficient heat transfer).

### Lysing

10. Add 2 ml Lysing Reagent, adjusted to 18-28 ° C, into each tube, mix gently.
11. Incubate for 5-10 minutes at 18-28°C.
12. Centrifuge the tubes for 5 minutes at 500 x g.
13. Decant the supernatant by using blotting paper.
14. Resuspend the cell pellet with 200 µl of Wash Buffer.

**Note:** Depending on Flow cytometry instrumentation, it might be necessary to increase the Wash Buffer volume (e.g. 800 µl).

15. Vortex gently.
16. Acquire the data on the flow cytometer within the same day. If the samples are stored for several hours before the analysis (<16 h) they should be kept protected from light at 2-8°C.

## FLOW CYTOMETRIC DATA ACQUISITION

Flow cytometric acquisition can be performed on any flow cytometer working with a 488 nm argon laser diode (blue excitation light) or a 532 nm green laser diode.

The flow cytometer must be equipped with Forward Scatter (FSC), Side Scatter (SSC), and suitable detection channels for the fluorochromes PE and PE-DY647 (analogue to PC5).

Ensure that the flow cytometer is properly aligned and colour compensation is set.

During data acquisition of samples, make sure that the leukocyte population is separated into three discrete populations on the FSC/SSC histogram. Adjust the amplification (gain) of FSC and SSC signals to obtain a cell distribution the one shown in Figure 1. Refer to the flow cytometer manual for further information.


Figure 1: Three discrete populations for lymphocytes (A), monocytes (B) and granulocytes (C) in FSC/SSC histogram.

Typically, the acquisition can be stopped, after having acquired 800-1000 basophilic cells (gated as shown in figure 2). At least 300 basophilic cells must be analyzed to obtain reliable results. This requires a total amount of 100'000-200'000 leukocytes per sample to be acquired.

The analysis of allergies to drugs is more difficult, because the activation of basophils in these patients is lower than in patients suffering from other allergies (e.g. to pollen allergens). Thus the lower limit of basophilic cells to be analyzed should be set to 800. We recommend each laboratory define its own confidence limits.

## DATA ANALYSIS

The analysis of acquired data can be performed with any flow cytometry analysis software like FlowJo, FloMax, CellQuest or others.

The analysis is done in two steps:

1. Include the entire basophil population  $CCR3^{pos}$  by setting the gate R1 with the low Side Scatter  $SSC^{low}$ . Eosinophils located at  $SSC^{high}$  position (upper right side) will be excluded (refer to Figure 2).
2. Display the cells selected in R1 in a graph with CCR-PE on the x-axis and CD63+CD203c-PE-DY647 on the y-axis. Set a quadrant gate so that the ratio of cells in the upper right quadrant for the base activation (PB-tube) is between 2.0 and 2.5% of the total number of basophilic cells gated in R1 (Figure 3).

**Note:** Set a new quadrant threshold for each patient. Once the quadrant has been defined by measuring the patient background (PB), use the same quadrant settings for

the analysis of the other tests to be carried out with the respective patient sample (positive controls and allergens. Refer to Figure 4.




Figure 2: Selection of basophilic cells (R1)  $CCR3^{pos} / SSC^{low}$

Gated Region	Count (n=)	%
Total	152483	
R1	1560	100.0
Q2 (CD63+CD203c <sup>pos</sup> )	33	2.1
Q3 (CD63+CD203c <sup>neg</sup> )	1527	97.9

Figure 3: Patient Background (PB) with CSB only

Gated Region	Count (n=)	%
Total	156543	
R1	1504	100.0
Q2 (CD63+CD203c <sup>pos</sup> )	1221	81.2
Q3 (CD63+CD203c <sup>neg</sup> )	283	18.8

Figure 4: Stimulation Control (STCON)

## LIMITATION

- Flow cytometry may produce false results if the cytometer has not been aligned properly, if the fluorescence emission has not been correctly compensated and if the gating regions have not been carefully positioned.
- Verify the preparations by eye to assess the efficacy of lysis. The erythrocytes may be incompletely lysed and appear on a light diffraction histogram in the same location as the lymphocytes.

## QUALITY CONTROL

For the appropriate evaluation of results, different values should be taken into account:

- Typically, three distinct **leukocyte populations** lymphocytes, monocytes and granulocytes appear in the FSC/SSC plot. Their occurrence can be regarded as a criterion for the quality of the blood sample (time of collection, storage).
- Negative control** (buffer control): A patient individual base value of 2.0 to 2.5 % of the activated basophils should be regarded as negative. This was also basis for the determination of allergen specific cut-off's.
- Positive control** (stimulation control). Two different positive controls are included in the kit. **Anti-FcεRI mAb** mimics the bridging of the receptor *in vivo* caused by the allergen. **fMLP** is a tripeptide causing basophil activation in a non-immunologic way. If one of those two stimulators shows activation of **>10%** basophils, the sample is evaluable.
- Non-responders** are persons with a low reactivity (<10% CD63+CD203c positive cells) to fMLP and anti-FcεRI antibody. Evaluation of data from normal blood donors and patients showed that 6.0% (out of n=116) were non-responders to anti-FcεRI and 2.9% (out of n=70) to fMLP.

## INTERPRETATION OF THE RESULTS

To obtain an optimal sensitivity and specificity, slightly different cut-off values should be applied for different groups of allergens. Based on numerous studies and evaluations BÜHLMANN recommends using the following cut-off values:

Inhalant allergens:	≥15%	
Food allergens:	≥15%	
Hymenoptera venoms:	≥10%	
Betalactams*:	≥5%	SI ≥2
Analgesics*:	≥5%	SI ≥2
Food additives*	≥5%	SI ≥2

\* Drugs and other chemical allergens usually results in lower activation percentages than other allergens. Therefore, a lower cut-off should be taken, but the stimulation index (SI = allergen stimulation divided by negative control) must be equal or higher than 2 in order to consider the result as positive.

BÜHLMANN allergens have been validated. Individual cut-offs have been established in studies or internal evaluations. For details refer to the allergen data sheet on the BÜHLMANN website, [www.buhmannlabs.ch](http://www.buhmannlabs.ch).

## ASSAY PERFORMANCE

**Specificity:** The anti-CCR3 mAb is a highly specific antibody (Refs 6 and 7). CCR3 is constitutively expressed on eosinophilic and basophilic leukocytes (see Fig. 2) and in a smaller part on CD3<sup>+</sup> cells (lymphocytes). Samples from eight normal blood donors were double stained twice with anti-CCR3-PE and anti-CD3-AF647. The relative

amount (mean) of CD3<sup>+</sup> cells within the gated Basophil population was 3.9% (cf Table 3).

n	16
Mean	3.85%
95 <sup>th</sup> percentile	1.93 – 5.41%

Table 3: Specificity The relative amount (mean) of CD3<sup>+</sup> cells within the gated Basophil population.

**Basophil Yield:** After the lysing step, cells from 43 normal blood donors were resuspended in 250 µl of wash-buffer and counted from a 100 µl aliquot of this cell suspension by using Partec PAS Flow Cytometer connected to Partec RobbyWell automate. The mean Basophil yield was 958 cells (95<sup>th</sup> percentile 310 to 2336 basophils).

**Precision:** A precision profile was calculated from differently stimulated blood samples (n=20). Each sample was incubated 12 times with stimulation buffer or stimulation control and consecutively analyzed by flow cytometry. The results are expressed as coefficient of variation (%CV) and plotted against % activation (CD63+CD203c<sup>POS</sup>). The precision profile is shown in Figure 5. At basophil activation above 10% the CV is below 15%.


Figure 5: Precision profile

**Inter Technician Variation: 6.1 - 11.7 % CV.** Two blood samples from asymptomatic blood donors were tested by six different technicians in two different labs within the same day. Both positive controls included in the kit, STCON and fMLP were used. The results are expressed in Table 4.

	Sample1 (%CD63+CD203c <sup>POS</sup> )		Sample2 (%CD63+CD203c <sup>POS</sup> )	
	STCON	FMLP	STCON	FMLP
Mean	81.1	58.0	83.5	69.5
SD	5.57	6.80	5.88	4.24
<b>%CV</b>	<b>6.9%</b>	<b>11.7%</b>	<b>7.0%</b>	<b>6.1%</b>

Table 4: Inter Technician Variation

**Reference Intervals** were established with normal blood donors (age 18-60). Median value for anti-FcεRI stimulated samples is 71.1% (percentile 25<sup>th</sup>: 38.1 %, 75<sup>th</sup>: 85.8%). Median value for fMLP stimulated samples is 45.7 % (percentile 25<sup>th</sup>: 30.0 %, 75<sup>th</sup>: 59.7%). Refer to Table 5.

	n	%CD63+CD203c <sup>POS</sup>		
		Median	2.5 <sup>th</sup> – 97.5 <sup>th</sup> percentile	25 <sup>th</sup> – 75 <sup>th</sup> percentile
STCON	116	71.1	3.4 – 97.9	38.1 – 85.8
fMLP	70	45.7	8.9 – 89.4	30.0 – 59.7

Table 5: Reference intervals

**Interferences:** It is recommended patients avoiding systemically administered antiallergenic drugs such as corticosteroids, chromoglycic acid (DSCG) for at least 24 hours prior to blood sampling.

EDTA and **Antihistamine drugs** listed in Table 6 were evaluated according the approved CLSI guideline for interference testing EP7-A2 (Ref. 8). No interference has been observed for antihistamines up to the listed concentrations (Ref. 9). High concentration of **EDTA** inhibits the expression of CD63 and CD203c thus may lead to lower activation values. Not completely filled venipuncture blood collection tubes leads to higher concentration of EDTA. Filling grade less than 50% leads to inhibition of about 40% of the maximal possible activation. Refer to Table 7.

Drug	Active component	No interference up to	
Telfast®	Fexofenadin	1.9	µg/mL
Aerius®	Des-Loratadin	4.0	µg/mL
Claritin®	Loratadin	0.017	µg/mL
Cetirizin-Teva®	Cetirizin	1.5	µg/mL
Xyzal®	Levo-Cetirizin	1.4	µg/mL

Table 6: Interference of antihistamine drugs

Tube filling-level	vol (mL)	Conc. EDTA (mg/mL)	Recovery
1/1	9	1.8	<b>100%</b>
3/4	6.75	2.4	<b>100%</b>
1/2	4.5	3.6	<b>61%</b>
1/4	2.25	7.2	<b>11%</b>
1/9	1	16.2	<b>9%</b>

Table 7: Interference of EDTA simulation of a “short draw”

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
1. Sainte-Laudy, J, et al. [Analysis of membrane expression of the CD63 human basophil activation marker. Applications to allergologic diagnosis]. *Allerg Immunol (Paris)* **26**, 211-4. (1994).
2. Sabbah, A and Sainte-Laudy, J. Flow Cytometry applied to the analysis of Lymphocyte and Basophil activation. *ACI International* **8**, 116-9 (1996).
3. Sanz, ML, et al. Flow cytometric basophil activation test by detection of CD63 expression in patients with immediate-type reactions to betalactam antibiotics. *Clin Exp Allergy* **32**, 277-86. (2002).
4. DeWeck, AL and Sanz, ML. Flow cytometric cellular allergen stimulation Test (FAST/Flow-CAST): technical and clinical evaluation of a new diagnostic test in allergy and pseudo-allergy. *ACI International* **14**, 204-215 (2002).
5. Gamboa, P et al. The flow-cytometric determination of basophil activation induced aspirin and other non-steroidal anti-inflammatory drugs (NSAIDs) is useful for in vitro diagnosis of the NSAIDs hypersensitivity syndrome. *Clin Exp Allergy* **34**, 1448-57 (2004)
6. Ugucioni, M., C. R. Mackay, et al.. "High expression of the chemokine receptor CCR3 in human blood basophils. Role in activation by eotaxin, MCP-4, and other chemokines." *J Clin Invest* **100**(5): 1137-43 (1997)
7. Ducrest, S., F. Meier, et al. (2005). "Flowcytometric analysis of basophil counts in human blood and inaccuracy of hematology analyzers." *Allergy* **60**(11): 1446-50.
8. Clinical Laboratory Standards Institute (CLSI): *Interference Testing in Clinical Chemistry*; approved guideline (EP7-A2), 2005.
9. Wolanczyk-Medrala A, et al. (2009). "A new variant of the basophil activation test for allergen-induced basophil CD63 upregulation. The effect of cetirizine." *J Investig Allergol Clin Immunol* **19**(6): 465-73.
10. Hausmann OV., et al. (2010) „Robust expression of CCR3 as a single basophil selection marker in flow cytometry” *Allergy* **66**(1): 85-91

**Flow CAST® highsens**

**TIME TO RESULT: ~ 1 HOUR**

\* **Note:** Depending on Flow cytometry instrumentation it might be necessary to increase the wash buffer volume (e.g. 800 µl).

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
<b>REF</b>	Catalogue number Bestellnummer Réf�rence du catalogue Numero di catalogo N�mero de cat�logo
<b>LOT</b>	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
	Contains sufficient for <n> tests Ausreichend f�r „n“ Ans�tze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Temperature limitation Zul�ssiger Temperaturbereich Limites de temp�rature Limiti di temperatura Limite de temperatura

Symbol	Explanation
<b>BUF STIM</b>	Stimulation Buffer Stimmulations-Puffer Tampon de stimulation tampone di stimolazione Tamp�n de estimulaci�n
<b>CONTROL STIM</b>	Stimulation Control Stimulationskontrolle Contr�le de stimulation Controllo di stimolazione Control de estimulaci�n
<b>CONTROL FMLP</b>	Stimulation Control fMLP Stimulationskontrolle fMLP Contr�le de stimulation fMLP Controllo di stimolazione fMLP Control de estimulaci�n fMLP
<b>REAG STAIN</b>	Staining Reagent F�rbe-Reagenz R�actif de coloration Reagente di colorazione Reactivo de coloraci�n
<b>REAG LYS</b>	Lysing Reagent Lyse Reagenz R�actif de lyse Reagente di lisi Reactivo de l�sis
<b>BUF WASH</b>	Wash Buffer Wasch-Puffer Tampon de lavage tampone di lavaggio Tamp�n de lavado

Printing Date  
2011-11-09