CAST[®]-Allergens

Revision date: 2007-11-06

ENGLISH

INTENDED USE

The CAST®-Allergens from BÜHLMANN are intended for *in vitro* leukocyte stimulation in the CAST®-Assays, namely CAST®-2000 ELISA (EK-CAST) and Flow-CAST® (FK-BAT) Kit. CAST®-Allergens that consist of proteins are allergen extracts presented in a stable, concentrated liquid form without covalent modifications. CAST®-Allergens such as chemicals and drugs are provided in lyophilized form. All CAST®-Allergens are preweighted and filled in a ready-to-dilute or ready-to-dissolve form, respectively. CAST®-Allergens are available as separate reagents and are not included in the CAST®-Assay kits.

INSTRUCTION FOR USE

Dissolve (or dilute) the allergen prior to use in the CAST®-Assays as follows:

- Add 250 µl of CAST® Stimulation Buffer to the vial and vortex briefly.
- Use allergen solution according to the CAST®-Assay protocols (see corresponding instructions for use).

Additionally to the highest concentration after the allergen dilution described above we recommend to stimulate the cells also with diluted allergen solution (further dilution of e.g. 1:5 with Stimulation Buffer). Some patients individually show a positive stimulation either at high or at low allergen concentrations only, whereas others are reacting positive over a broad range of allergen concentrations.

STORAGE AND SHELF LIFE

Store the unopened CAST®-Allergens at 2-8°C until the expiration date marked on the vial.

IMPORTANT: Use a fresh vial of allergen each time a new cell stimulation is performed!

POINT OF TIME FOR PATIENT EVALUATION

The optimal time for diagnostic evaluation of a patient with the CAST®-Assays is between 3 and 12 weeks after an adverse reaction to suspected allergens, although specific basophil reactivity in many cases persists much longer. The CAST®-Assays can also detect a potentially adverse reaction to a specific allergen even without history of previous clinical reaction to that allergen. In some circumstances, the CAST®-Assays may therefore be used as a screening assay or a prognostic test for future adverse reactions within a series of related drugs (e.g. Betalactam antibiotics, non steroidal anti-inflammatory drugs).

The blood sample for the CAST®-Assays must be collected before skin testing or *in vivo* provocation of the patient with a culprit or suspected allergen. The exposition to the allergen may cause a general activation of the patients leukocytes.

LIMITATIONS

A NEGATIVE CAST®-ASSAY RESULT FOR A SPECIFIC ALLERGEN CAN NOT EXCLUDE THE POTENTIAL OCCURRENCE OF A (EVEN SEVERE) CLINICAL REACTION IN A PATIENT. Patients with a history of adverse reactions to protein or drug allergen with a negative CAST®-Assay result should therefore be followed up with an additional method, such as an *in vivo* provocation or skin prick test (where appropriate), before any drug will be administered or the patient will be exposed to the allergen.

PROTEIN ALLERGENS AND MIXES

This instruction applies for the following group of inhalant and food allergens and allergen mixtures thereof:

- GRASSES, WEEDS AND TREES
- MOLDS
- INSECT VENOMS, INSECTS, MITES AND HIGHER ANIMALS
- ALBUMINS
- FOOD

ALLERGEN CONCENTRATION AND CUT-OFF

In general, the protein allergens contain 25 ng concentrated protein extract per vial. The dilution of the allergen with 250 µl Stimulation Buffer results in a concentration of 100 ng/ml in the vial.

CAST®-2000 ELISA

For INHALATION and FOOD ALLERGENS as well as for LATEX and α-AMYLASE, we propose that individuals with a net stimulation yield **HIGHER THAN 200 pg/ml sLT** should be regarded as positive for the allergen tested.

In a recent study (1) the clinical cut-off for **BEE** (BAG2-I1) and **WASP VENOM** (BAG2-I3) was set to **270 pg/ml** with a grey zone from 200 to 270 pg/ml.

Flow-CAST®

Samples stimulated with INHALANT or FOOD ALLERGENS showing a basophil activation rate (% CD63 positive) ≥ 15 % can be regarded as positive. Samples stimulated ≥ 10 % with INSECT VENOM or LATEX can be regarded as positive

For further information refer to the allergen data sheets available on www.buhmannlabs.ch/2-1-5-Allergens.php.

DRUG ALLERGENS, CHEMICAL ALLERGENS AND FOOD ADDITIVES

This instruction applies for the following group of allergens and allergen mixtures thereof:

- ANTIBIOTICS
- ANALGESICS
- ANESTHETICS / MYORELAXANTS
- FOOD ADDITIVES
- ENVIRONMENTAL / OCCUPATIONAL.

CONCENTRATION AND CUT-OFF

All drug, chemical and food additive allergens are characterized for their inhibition and unspecific stimulatory properties. Based on those evaluations, the highest possible allergen concentration for the cell stimulation was established. For the allergen concentrations refer to the individual allergen data sheets available on www.buhmannlabs.ch.

In general, drugs, chemicals and food additives usually cause lower sLT release and basophil activation, respectively. Therefore lower cut-off values should be applied for these allergens. The proposed cut-off values for each single allergen can be found on the individual allergen data sheets available on www.buhmannlabs.ch.

CAST®-2000 ELISA

Due to the small increase of sLT release with this kind of allergens Bühlmann established for each single allergen an INDIVIDUAL TECHNICAL CUT-OFF. The values represent the mean + 3SD of net stimulation of up to 20 stimulated samples from normal blood donors. Due to inter-laboratory variation and the precision of the assay the cut-off values recommended are not set below 40 pg/ml.

Notes:

- Bühlmann recommends using a grey zone including ±20% of the individual technical cut-offs.
- To use the individual technical cut-off the assay quality control criteria described in the instruction for use of the EK-CAST must be fulfilled.

- The cut-off values should be used as guidelines only. A clinical relevant cut-off should be established by each laboratory and by further studies, respectively.

Flow-CAST®

Due to the small increase of basophil activation with these types of allergens, Bühlmann propose for the Flow-CAST® the use of a stimulation index (SI) as an additional discrimination criterion. The SI is defined as the gross stimulation divided by the negative control value.

RECOMMENDATIONS FOR THE USE OF THE CAST®-ASSAYS IN DRUG ALLERGY DIAGNOSIS

Drug allergy diagnosis is a complex domain with heterogeneous symptoms and pathophysiology (2-4). Therefore, and AS A GUIDELINE only, the following recommendations should be regarded when using the CAST®-Assays to test patients for allergy to drugs and/or chemicals.

CLINICAL CONDITIONS SUITABLE FOR EVALUATION WITH CAST®-ASSAYS

Immediate adverse reactions to drugs and chemicals may be caused *in vivo* by both, IgE and non-IgE mediated (so-called pseudo-allergies) stimulation of effector cells (e.g. basophils, mast cells and eosinophils). The CAST®-Assays could be regarded as a model of the above events, in which patients leukocytes are stimulated *in vitro*. However, as CAST®-Assays use pure allergens and a defined buffer system without plasma components (e.g. coagulation factors or complement), the method can only be considered as a simplified model of the *in vivo* situation.

The CAST®-Assays are optimized for the detection of immediate type allergic (type I) and pseudo-allergic reactions only. The CAST®-Assays are not likely to detect delayed types of allergic reactions (e.g. type IV reactions). Therefore, it should only be used for the diagnosis of immediate type reactions, particularly clinical symptoms as follows:

- Anaphylaxis and anaphylactoid reactions
- Rhinoconjunctivitis
- Asthma bronchiale
- Urticaria/Angiodema
- Gastrointestinal reactions.

CLINICAL CONDITIONS NOT SUITABLE FOR EVALUATION WITH CAST®-ASSAYS

The CAST®-Assays are less likely to yield positive results in cases of T-cell mediated reactions. Based on the known pathogenetic mechanisms, we do not recommend evaluating routinely the following group of patients with the CAST®-Assays showing:

- Maculopapular exanthema and Pustular exanthema
- Vasculitis
- Allergic contact dermatitis
- Bullous exanthemas, such as Stevens-Johnson syndrome and Lyell's disease.

CUSTOMER SPECIFIC ALLERGENS

The same rules as for the establishment of the Bühlmann CAST®-Allergens should be applied. Allergens from other sources can be used by considering to following limitations:

- No matrix-bound allergens (solid or liquid phase)
- No allergen preparations containing sulfidoleukotrienes
- No allergen preparations containing cytotoxic compounds (stabilizers, preservatives) such as phenol, glycerol, sodium azide or merthiolate (thimerosal).

For the procedure to establish customer specific allergens for the CAST®-Assays ask your local distributor or contact Bühlmann Laboratories AG. Specific information is available upon request.

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