# CALPROTECTIN

**ELISA**

EK-CAL2 2 x 96 tests

Revision date: 2012-01-20

## ENGLISH

### INTENDED USE

The BÜHLMANN Calprotectin ELISA kit is designed for the extraction and quantitative determination of human Calprotectin (MRP8/14; S100A8/S100A9) in stool samples (1-3).

### PRINCIPLE OF THE ASSAY

After a short extraction procedure using one volume of faeces and 49 volumes of Extraction Buffer, the test allows for the selective measurement of Calprotectin-antigen by sandwich ELISA. A monoclonal capture antibody (mAb) highly specific to the Calprotectin heterodimeric and polymeric complexes (4-5), respectively, is coated onto the microtiter plate. Calibrators, controls and patients extracts are incubated at room temperature for 30 minutes. After a washing step a detection antibody (Ab) conjugated to horseradish peroxidase (HRP) detects the calprotectin molecules bound to the monoclonal antibody coated onto the plate. After incubation and a further washing step, tetramethylbenzidine (TMB) will be added (blue color formation) followed by a stopping reaction (change to yellow color). The absorption is measured at 450 nm.

### REAGENTS SUPPLIED AND PREPARATION

**Standards and controls are sufficient for up to 4 runs in duplicates!**

Reagents	Quantity	Code	Reconstitution
<b>Extraction Buffer</b>	6 bottles 125 ml	B-CAL-EX	Ready to use
<b>Microtiter Plate</b> precoated with anti-Calprotectin mAb	2x12 x 8 wells	B-CAL-MP	Ready to use
<b>Plate Sealer</b>	6 pieces		
<b>Wash Buffer Concentrate (10x)</b> with preservatives	2 bottles 100 ml	B-CAL-WB	Dilute with 900 ml of deionized H <sub>2</sub> O
<b>Incubation Buffer</b> with preservatives	3 bottles 125 ml	B-CAL-IB	Ready to use
<b>Calibrators A to E</b> <sup>1)2)</sup> Calprotectin in a buffer matrix with preservatives	5 vials 1 ml	B-CAL-CASET	Ready to use
<b>Control Low / High</b> <sup>3)</sup> human serum with preservatives	2 vials 1 ml	B-CAL-CONSET	Ready to use
<b>Enzyme Label</b> Anti-Calprotectin Ab conjugated to HRP	2 vials 12 ml	B-CAL-EL	Ready to use
<b>TMB Substrate</b> TMB in citrate buffer	2 vials 12 ml	B-TMB12	Ready to use
<b>Stop Solution</b> 0.25 M sulfuric acid	2 vials 12 ml	B-ST512	Ready to use <b>Corrosive agent</b>

Table 1

- <sup>1)</sup> The actual Calprotectin concentration of the standards A to E are 4, 12, 40, 120 and 240 ng/ml, respectively. During extraction a 1:50 sample dilution occurs followed by an additional 1:50 dilution of the extracts for the measurement in the ELISA. To take these dilution steps into account for the final calculations the calibrators A to E the following concentrations have to be used for the lower range ELISA procedure (refer to page 4): **10, 30, 100, 300 and 600 µg/g of Calprotectin**.
- <sup>2)</sup> If you choose the extended range ELISA procedure (refer to page 5) the following calibrator concentrations have to be used in the respective ELISA protocol: **30, 90, 300, 900 and 1800 µg/g of Calprotectin**
- <sup>3)</sup> The controls contain lot specific amounts of native human Calprotectin. Refer to the additional QC data sheet for actual concentrations.

### REAGENTS SUPPLIED UPON REQUEST

#### Fecal Extraction Devices

<b>Smart-Prep</b>	50 tubes, spatulas, and base caps	B-CAL-RD
<b>Schebo® Quick-Prep™</b>	50 tubes consisting of tube, cone & dosing tip	B-CAL-SOFI
	1.3 ml, ready to use	

Table 2

### STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Store at 2-8°C. Do not use past kit expiration date printed on the labels.	
Opened / Reconstituted Reagents	
Extraction Buffer	Store at 2-8°C until expiration date.
Microtiter Plate	Return unused strips immediately to the foil pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store until expiration date at 2-8°C.
Diluted Wash Buffer	Store for up to 6 months at 2-8°C.
Incubation Buffer	Store at 2-8°C until expiration date.
Calibrators	
Controls	
Enzyme Label	
TMB-Substrate	
Stop Solution	

Table 3

### WARNINGS AND PRECAUTIONS

The Microtiter Strips, Calibrators and Controls of this kit contain components of human origin. Each serum donor unit used in the preparation of the kit components was tested by an FDA approved method and found negative for HBV surface antigen, so as for HCV and HIV1/2 antibodies. Although these methods are highly accurate, there is no guarantee that this material cannot transmit Hepatitis or AIDS. *Therefore, all patient specimens and kit components should be handled as if capable of transmitting infections.* All products containing human source material should be handled in accordance with good laboratory practice using appropriate precautions.

**Stop Solution:** The Stop Solution (B-ST5) contains sulfuric acid. The reagent is irritating to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothing. Wear suitable protective clothing, gloves and eye protection. After contact with eyes or skin, wash immediately with plenty of water.

### MATERIALS REQUIRED BUT NOT PROVIDED

#### Extraction Procedure

- 10 µl disposable inoculation loops (Sarstedt: # 86.1562.010)
- 15 ml polypropylene tubes with screw caps (Sarstedt: #62.554.502) required for standard extraction procedure; extraction devices (see above).
- Laminar flow work station
- Multi tube vortex mixer
- Precision balance (10-150 mg)
- Micro centrifuge (≥3000 g)

#### ELISA Procedure

- 10, 100 and 1000 µl precision pipettes with disposable tips.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 ml cylinder for the dilution of the Wash Buffer Conc.

- Microtiter plate washer (see Procedural Notes) or squeeze bottle for Wash Buffer.
- Microtiter plate rotator (see Procedural Notes).
- Blotting paper.
- Microtiter plate reader for measurement of absorbance at 450 nm.

### SPECIMEN COLLECTION AND STORAGE

50 to 100 mg of native stool sample are needed for the extraction procedure.

Collect stool samples into plain tubes and store them refrigerated at 2-8°C for up to 6 days.

Freezing of samples may result in slightly increased Calprotectin concentrations due to Neutrophils present in the sample. Therefore, it is of advantage for longer storage, to keep the extracts at -20°C. The extracts are stable for at least 4 months.

**Important:** The sample must be collected without any chemical or biological additions in the collection device.

### PROCEDURAL NOTES

- **Residues in the microtiter plate wells** result from the production process. They are removed in the washing step (Assay procedure step 3) and do not affect the results.

#### Extraction:

- To receive quantitative results it is important to homogenate the entire weighted stool sample in the extraction buffer. Avoid contamination at the top of the tube - insoluble (undigested) components can still be in the tube after extraction.
- The short centrifugation step of 5 minutes during the extraction can cause a turbid solution. Turbidity can be avoided by longer centrifugation but it shows no influence on the quantitative determination in the ELISA.

#### ELISA Procedure:

- In the ELISA procedure the **washing steps are essential** to guarantee reproducible results. A **minimal incubation time** of the Wash Buffer in the wells of **at least 20 seconds** must be ensured each time.
- When using **automated washer**, BÜHLMANN strongly recommends to use a "plate mode" i.e. each process step (dispense/aspiration) is performed on all of the strips sequentially, before processing to the next process step. Thus, the minimal incubation time is guaranteed.
- The indicated no. of **washing cycles is mandatory** to ensure reproducible results.
- To ensure a complete antigen/antibody interaction, the **incubation time in step 5** must be at least 30 minutes. Moderately longer incubation time (up to 5 minutes) has no influence to the final outcome.
- **Plate rotator** (shaker) must at 400 - 600 rpm (<10 Hz). Higher rotation frequency may cause poor dilution linearity at values between 300/900 and 600/1800 µg/g. Orbital rotation instead of reciprocal shaking should be used.
- The Enzyme Label is inactivated by oxygen and is highly sensitive to sodium azide, thimerosal, hypochlorous acid and aromatic chlorohydrocarbons often found in laboratory water supplies. Therefore, use only deionized high quality water.
- It is recommended to assay each control and specimen in **duplicate** each time a test is performed. Since conditions vary from assay to assay, a new standard curve must be generated each time a new assay is performed. Vertical alignment is recommended.
- If the initial concentration of an unknown sample reads higher than the top Calibrator (Calibrator E), the sample

must be further diluted with Incubation Buffer and assayed again according to the assay procedure. The resulting dilution factor must be accounted for the calculation of results.

### ASSAY PROCEDURE

#### Extraction

**Standard extraction procedure** (refer to page 31)

1. Label and weigh (tare) the empty polypropylene tube together with the inoculation loop.
2. Take out 50 to 100 mg of the stool sample by means of the inoculation loop and place it into the pre-weighted tube.
3. Estimate the net amount of sample, break off the inoculation loop and leave the lower part of the loop in the tube.
4. Add Extraction Buffer (49 times the weight volume) to the tube and close the tube:

Weight [mg] Stool	Volume [ml] Extraction Buffer
50	2.5
55	2.7
60	2.9
65	3.2
70	3.4
75	3.7
80	3.9
85	4.2
90	4.4
95	4.7
100	4.9

5. Homogenize the sample on a multi tube vortexer by vigorous shaking (at highest speed) for 30 minutes.
6. Transfer the homogenate into a 2 ml Eppendorf tube and centrifuge in a microcentrifuge for 5 minutes at 3'000 x g.
7. Take the supernatant into a fresh, labeled tube and continue with the ELISA procedure or store the extracts at ≤-20°C for at least 4 months.

#### Extraction procedures using fecal extraction devices:

The extraction procedure is described and illustrated in the instruction for use delivered with the respective extraction device.

1. Fecal Extraction Device Roche (Code 10745804 322) or BÜHLMANN **Smart-Prep** (Code: B-CAL-RD): The extraction time (vortexing) can be reduced to 1 minute.
2. **ScheBo® Quick-Prep™** (Code B-CAL-SOFI): The extraction tubes are prefilled with extraction buffer. The extraction time is about 10 minutes (vortexing).

After extraction, centrifuge the tubes for 5 minutes at 3'000 x g. Alternatively, transfer the homogenate into a 2 ml Eppendorf tube and centrifuge it in a microcentrifuge for 5 minutes at 3'000 x g.

Decant the supernatant into a fresh, labeled tube and continue with the ELISA procedure or store the extracts at ≤-20°C for at least 4 months.

The respective extraction procedures are published on the website:  
<http://www.buhlmannlabs.ch/core/inflammation/calprotectin/>

---

## ELISA PROCEDURES

The assay can be performed according to the following procedures – lower or extended range ELISA procedure. Which procedure is to be chosen depends on the expected Calprotectin concentration of the samples. For samples up to 600 µg/g choose the lower range procedure using a sample dilution of 1:50 (working range 10 – 600 µg/g) (refer to page 4, 5). If the samples tend to exceed 600 µg/g, choose the extended range procedure using a sample dilution of 1:150 (working range 30 – 1800 µg/g) (refer to page 5).

---

### LOWER RANGE ELISA PROCEDURE

#### WORKING RANGE 10 – 600 µg/g

Allow the reagents to equilibrate to 18-28°C prior to use

1. Dilute the stool extracts 1:50 with Incubation Buffer (e.g. 20 µl extract and 980 µl incubation buffer) and mix well. Let the samples equilibrate for at least 5 minutes at 18-28°C prior to proceeding to step 4c.
  2. Prepare a plate with sufficient strips to test the required number of calibrators, controls and diluted samples. Remove excess strips from the holder and re-seal them in the foil pouch together with the desiccant packs **without delay**. Store refrigerated.
  3. Wash the coated wells twice using at least 300 µl of Wash Buffer per well. Empty the wells and tap the plate firmly onto blotting paper.
- Important:** For every of the three wash steps a minimal incubation time of at least 20 seconds of the Wash Buffer in the wells must be ensured (see Procedural Notes – ELISA Procedure, page 3).
- 4a. Pipet 100 µl of Incubation Buffer in duplicate into wells A1+A2 (Blank)  
Pipet 100 µl of Calibrator A in duplicate into wells B1+B2  
Pipet 100 µl of Calibrator B in duplicate into wells C1+C2. etc.
  - 4b. Pipet 100 µl of the Low and High Controls in duplicate into wells G1+G2 and H1+H2, respectively.
  - 4c. Pipet 100 µl of each diluted sample in duplicate into the subsequent wells.
  5. Cover the plate with a plate sealer, and incubate for 30 + 5 minutes on a plate rotator set at 400-600 rpm at 18-28°C (see Procedural Notes – ELISA Procedure, page 3).
  6. Remove and discard the plate sealer. Empty the wells and wash three times using at least 300 µl of Wash Buffer per well (see procedural notes, page 3). Empty the wells and tap the plate firmly onto blotting paper.
  7. Pipet 100 µl of Enzyme Label to all wells.
  8. Cover the plate with a plate sealer, and incubate for 30 ± 5 minutes on a plate rotator set at 400-600 rpm at 18-28°C.
  9. Remove and discard the Plate Sealer. Empty the wells and **wash five times** using at least 300 µl of Wash Buffer per well. Empty the wells and tap the plate firmly onto blotting paper.

**Important:** Allow the TMB Substrate Solution to equilibrate to 18-28°C.

10. Pipet 100 µl of the TMB Substrate Solution to all wells.

11. Cover the plate with a plate sealer, protect the plate from direct light and incubate for 15 ± 2 minutes on a plate rotator set at 400-600 rpm at 18-28°C.

12. Pipet 100 µl of Stop Solution to all wells. Remove air bubbles with a pipette tip. Proceed to step 13 within 30 minutes.

13. Read the absorbance at 450 nm in a microtiter plate reader.

---

## RESULTS

### WORKING RANGE 10 – 600 µg/g

**Standard Curve:** Record the absorbance at 450 nm for each calibrator and blank well. Average the duplicate values, subtract the average of the blank wells and record averages (= corrected average absorbance). Plot the absorbance (vertical axis) versus the Calprotectin concentration of the calibrators (horizontal axis) using a semi logarithmic lin/log graph paper. Draw the best fitting curve or calculate the standard curve using a four parameter logistic.

**Samples and Controls:** Record the absorbance at 450 nm for each sample and control well. Average the duplicate values, subtract the average of the blank wells and record the averages (=corrected average absorbance). Locate the corrected absorbance value of the sample on the vertical axis, draw a horizontal line intersecting the standard curve and read the Calprotectin concentration from the horizontal axis. If you choose the lower range ELISA procedure, the following calibrator concentrations have to be used in the respective ELISA protocol: **10, 30, 100, 300 and 600 µg/g of Calprotectin.**

**Additional dilution factors have to be multiplied with the results to obtain the final results.**

Refer to Table 16 and Figure 1 for typical data (results and standard curve). *These results and standard curve are provided for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.*

---

### PERFORMANCE CHARACTERISTICS

#### WORKING RANGE: 10 – 600 µg/g

**Intra-Assay Precision: 4.7 %.** The intra-assay precision was calculated from 20 pairs of values from 3 extracted stool samples assayed in a single run according to the assay procedure. The values are presented in Table 17.

**Inter-Assay Precision: <15 %.** The inter-assay precision of the ELISA was calculated from 5 extracted stool samples. The aliquots were tested according to the assay procedure in 10 different runs by three technicians using 2 kit lots in two different labs. The values are presented in Table 18.

**Analytical Sensitivity: <10 µg/g.** Twenty duplicates of Incubation Buffer were assayed in a single run. Mean and standard deviation were calculated for the absorbance values. The minimal detectable dose of Calprotectin was calculated to be clearly below Calibrator A (10 µg/g) by adding two standard deviations to the mean absorbance and intersecting this value with the standard curve obtained in a new run.

**Functional Sensitivity: <10 µg/g.** Ten stool samples with values between 5.2 and 1254 µg/g Calprotectin were assayed 20 times in duplicates in one assay. The %CV and the mean values were calculated for each sample. The functional sensitivity was observed at 15 % CV. The resulting precision profile (Figure 2) allows the precise measurement within the whole standard range between 10 and 600 µg/g.

**Dilution Linearity: 103 %.** Seven stool samples with elevated Calprotectin values were extracted according to the assay procedure. The extracts were diluted with Incubation

Buffer and subsequently assayed according to the assay procedure. The expected values were calculated from the observed value found with the first dilution. The results are presented in Table 19.

**Spiking Recovery: 100 %.** Two extracted stool samples were spiked with different amounts of diluted, Calprotectin containing human serum. The samples were measured before and after spiking according the assay procedure. The results are presented in Table 20.

**Crossreactivity: <0.1 %.** Incubation Buffer spiked with different amounts of recombinant MRP8 and MRP14 were measured according to the assay procedure. The values are presented in Table 25.

### Working Range 30 – 1800 µg/g

#### Allow the reagents to equilibrate to 18-28°C prior to use

The working range can be extended by a factor of 3, if you dilute the samples 1:150 instead of 1:50. This procedure is recommended, if high Calprotectin concentrations are to be expected. Precision and linearity of the assay allow for this extension of the kit range.

1. Dilute the stool extracts 1:150 with Incubation Buffer (e.g. 20 µl extract and 2980 µl incubation buffer) and mix well. **NB: Only dilute stool extracts. Standards and controls are ready to use.** Let the samples equilibrate for at least 5 minutes at 18-28°C prior to proceeding to step 4c.
2. Prepare a plate with sufficient strips to test the required number of calibrators, controls and diluted samples. Remove excess strips from the holder and re-seal them in the foil pouch together with the desiccant packs **without delay**. Store refrigerated.
3. Wash the coated wells twice using at least 300 µl of Wash Buffer per well. Empty the wells and tap the plate firmly onto blotting paper.

**Important:** For every of the three wash steps a minimal **incubation time of at least 20 seconds of the Wash Buffer** in the wells must be ensured (see Procedural Notes – ELISA Procedure).

- 4a. Pipet 100 µl of Incubation Buffer in duplicate into wells A1+A2 (Blank)
    - Pipet 100 µl of Calibrator A in duplicate into wells B1+B2
    - Pipet 100 µl of Calibrator B in duplicate into wells C1+C2. etc.
  - 4b. Pipet 100 µl of the Low and High Controls in duplicate into wells G1+G2 and H1+H2, respectively.
  - 4c. Pipet 100 µl of each diluted sample in duplicate into the subsequent wells.
  5. Cover the plate with a plate sealer, and incubate for 30 + 5 minutes on a plate rotator set at 400-600 rpm at 18-28°C (see Procedural Notes – ELISA Procedure, page 3).
  6. Remove and discard the plate sealer. Empty the wells and wash three times using at least 300 µl of Wash Buffer per well (see procedural notes). Empty the wells and tap the plate firmly onto blotting paper.
  7. Pipet 100 µl of Enzyme Label to all wells.
  8. Cover the plate with a plate sealer, and incubate for 30 ± 5 minutes on a plate rotator set at 400-600 rpm at 18-28°C.
  9. Remove and discard the Plate Sealer. Empty the wells and **wash five times** using at least 300 µl of Wash Buffer per well. Empty the wells and tap the plate firmly onto blotting paper.
- Important:** Allow the TMB Substrate Solution to equilibrate to 18-28°C.
10. Pipet 100 µl of the TMB Substrate Solution to all wells.
  11. Cover the plate with a plate sealer, protect the plate from direct light and incubate for 15 ± 2 minutes on a plate rotator set at 400-600 rpm at 18-28°C.

12. Pipet 100 µl of Stop Solution to all wells. Remove air bubbles with a pipette tip. Proceed to step 13 within 30 minutes.

13. Read the absorbance at 450 nm in a microtiter plate reader.

---

## RESULTS

### WORKING RANGE 30 – 1800 µg/g

**Standard Curve:** Record the absorbance at 450 nm for each calibrator and blank well. Average the duplicate values, subtract the average of the blank wells and record averages (= corrected average absorbance). Plot the absorbance (vertical axis) versus the Calprotectin concentration of the calibrators (horizontal axis) using a semi logarithmic lin/log graph paper. Draw the best fitting curve or calculate the standard curve using a four parameter logistic.

**Samples and Controls:** Record the absorbance at 450 nm for each sample and control well. Average the duplicate values, subtract the average of the blank wells and record the averages (=corrected average absorbance). Locate the corrected absorbance value of the sample on the vertical axis, draw a horizontal line intersecting the standard curve and read the Calprotectin concentration from the horizontal axis. If you choose the extended range ELISA procedure, the following calibrator concentrations have to be used in the respective ELISA protocol: 30, 90, 300, 900 and 1800 µg/g Calprotectin.

Refer to Table 21 and Figure 3 for typical data (results and standard curve). *These results and standard curve are provided for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed.*

---

## PERFORMANCE CHARACTERISTICS

### WORKING RANGE: 30 – 1800 µg/g

**Intra-Assay Precision: 4.0 %.** The intra-assay precision (mean) was calculated from the results of 20 duplicates from 3 extracted stool samples assayed in a single run according to the assay procedure. The values are presented in Table 22.

**Inter-Assay Precision: <15 %.** The inter-assay precision of the ELISA was calculated from 5 extracted stool samples. The aliquots were tested according to the assay procedure in 10 different runs by three technicians using 2 kit lots in two different labs. The values are presented in Table 23.

**Functional Sensitivity: <30 µg/g.** 18 stool samples with values between 10.8 and 2080 µg/g Calprotectin were measured 20 times in one assay. The % CV and the mean values were calculated for each sample. The functional sensitivity was observed at 15 % CV. The resulting precision profile allows the precise measurement within the whole standard range from 30 to 1800 µg/g. The results are presented in Figure 4.

**Dilution Linearity: 102 %.** Five stool samples with elevated Calprotectin concentrations were extracted according to the assay procedure. The extracts were diluted with Incubation Buffer and subsequently assayed according to the assay procedure. The expected values were calculated from the observed value found with the first dilution. The results are presented in Table 24.

---

## EXPECTED VALUES

Estimation of faecal Calprotectin is a reliable and easy way to distinguish organic from functional gastrointestinal diseases.

In a clinical study 401 symptomatic patients scheduled for colonoscopy were investigated (publication in preparation). Endoscopy examination showed 273 patients with

functional diseases whereas 128 patients had various organic diseases (colitis, Morbus Crohn, diverticulitis, polyps, adenomas, cancer, or infectious diseases).

ROC curve analysis (AUC: 0.935) resulted in an optimal clinical cut-off at 50 µg/g. Applying this cut-off, a clinical sensitivity and specificity of 84.4% and 94.5%, respectively can be reached in the differentiation between organic and functional diseases (see Table 26).

Faecal Calprotectin levels from adults and children are comparable, whereas levels of newborns can be significantly increased (8).

These data support the following recommendation for interpretation of results:

### Normal values below 50 µg/g:

Calprotectin values <50 µg/g are not indicative of inflammation in the gastrointestinal tract. Patients with low Calprotectin levels are likely not to be in need of invasive procedures to determine the inflammation cause (12).

### Elevated values between 50 and 200 µg/g:

Calprotectin values between 50 and 200 µg/g can represent mild organic disease such as inflammation caused by NSAIDs, mild diverticulitis and IBD in remission phase. The low inflammatory response shown within this range may suggest repeating the measurement and performing further investigations.

### Elevated values above 200 µg/g:

Calprotectin values > 200 µg/g are indicative of active organic disease with inflammation in the gastrointestinal tract. Appropriate further investigative and curative procedures by specialists are suggested.

The cut-off level suggested for adults (50 µg/g) can also be used for children aged from 4 to 17 years regardless of sex (9).

---

## QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

Since there is no control serum for Calprotectin commercially available, we recommend using a pool of positive stool extractions for internal quality control.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability. The confidence limits for the Controls are lot-specific and printed on the additional QC data sheet.

If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) ELISA reader settings iii) expiration dates of reagents iv) storage and incubation conditions v) TMB Substrate Solution should be colorless vi) purity of water.

---

## LIMITATIONS

- The reagents supplied with this kit are optimized to measure human Calprotectin in extracted stool samples.
- Faecal Calprotectin values should be used as supplementary data available to the physician in establishing a diagnosis.

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Figure 1: Example of a Standard Curve


Table 17: Intra-Assay Precision

Sample	Mean [µg/g]	SD [µg/g]	CV [%]
low	52.5	1.4	2.7
medium	173.8	5.6	3.2
high	408.5	33.0	8.1
<b>Mean</b>			<b>4.7</b>

Table 19: Dilution Linearity/Parallelism

Sample	Dilution	Observed [µg/g]	Expected [µg/g]	O/E [%]
S 1	1:200	405	-	--
	1:400	182	202	90
	1:800	95	101	94
	1:1600	49	51	98
	1:3200	25	25	99
	1:6400	15.6	12.7	123
	1:12800	6.6	6.3	104
<b>Mean</b>				<b>101</b>
S 2	1:50	232	-	--
	1:100	124	116	107
	1:200	61	58	106
	1:400	28	29	96
	1:800	12	15	86
<b>Mean</b>				<b>98</b>
S 3	1:400	290	-	--
	1:800	145	145	100
	1:1600	73	72	100
	1:3200	39	36	106
	1:6400	19	18	103
<b>Mean</b>				<b>103</b>
S 4	Mean			<b>100</b>
S 5	Mean			<b>124</b>
S 6	Mean			<b>88</b>
S 7	Mean			<b>106</b>
<b>Mean</b>				<b>103</b>

Table 18: Inter-Assay Precision

Sample	Mean [µg/g]	SD [µg/g]	CV [%]
Sample 1	18.1	2.5	13.5
Sample 2	44.5	6.4	14.5
Sample 3	74.3	7.9	10.7
Sample 4	227	15.0	6.6
Sample 5	520	57.8	11.1
<b>Mean</b>			<b>11.3</b>

Figure 2: Precision Profile


Table 20: Spiking Recovery

Sample	µg/g	spiked with [µg/g]	Observed [µg/g]	Expected [µg/g]	O/E [%]
S 1	5.3	10	13.6	15.3	88.9
		30	33.6	35.3	95.2
		100	101	105	96.1
		150	147	155	94.5
		300	287	305	94.1
		400	468	405	115.4
<b>Mean</b>				<b>97.4</b>	
S 2	24.6	10	32.9	34.6	95.1
		30	49.2	54.6	90.1
		100	139	125	111.5
		150	176	175	100.5
		300	358	325	110.3
		400	467	425	109.9
<b>Mean</b>				<b>103</b>	
<b>Mean</b>				<b>100</b>	

EXTENDED RANGE PROCEDURE 30-1800 µg/g

Figure 3

Example of a Standard Curve

Table 21: Example of Results

	Conc. [µg/g]	Absorb. [OD]	Calc. Conc. [µg/g]	CV Conc [%]
<b>Blank Avg.</b>		<b>0.057</b>		
Cal A	30	0.047		
Cal A	30	0.046		
<b>Cal A Avg.</b>	<b>30</b>	<b>0.047</b>		<b>0.9</b>
Cal B	90	0.138		
Cal B	90	0.140		
<b>Cal B Avg.</b>	<b>90</b>	<b>0.139</b>		<b>1.0</b>
Cal C	300	0.464		
Cal C	300	0.452		
<b>Cal C Avg.</b>	<b>300</b>	<b>0.458</b>		<b>1.9</b>
Cal D	900	1.207		
Cal D	900	1.192		
<b>Cal D Avg.</b>	<b>900</b>	<b>1.200</b>		<b>0.8</b>
Cal E	1800	1.627		
Cal E	1800	1.630		
<b>Cal E Avg.</b>	<b>1800</b>	<b>1.629</b>		<b>0.1</b>
Ctrl Low		0.147	105	
Ctrl Low		0.162	115	
<b>Ctrl Low Avg.</b>		<b>0.155</b>	<b>110</b>	<b>6.2</b>
Ctrl High		0.618	396	
Ctrl High		0.618	396	
<b>Ctrl High Avg.</b>		<b>0.618</b>	<b>396</b>	<b>0.6</b>


Table 22: Intra-Assay Precision

Sample	Mean [µg/g]	SD [µg/g]	CV [%]
Sample 1	37.7	2.5	6.7
Sample 2	713	20.1	2.8
Sample 3	1246	32.4	2.6
<b>Mean</b>			<b>4.0</b>

Figure 4:

Precision Profile


Table 23: Inter-Assay Precision

Sample	Mean [µg/g]	SD [µg/g]	CV [%]
Sample 1	75.5	9.7	12.8
Sample 2	225	20	8.8
Sample 3	788	62	7.8
Sample 4	1000	111	11.1
Sample 5	1764	221	12.6
<b>Mean</b>			<b>10.6</b>

Table 24:

Dilution Linearity/Parallelism

Sample	Dilution	Observed [µg/g]	Expected [µg/g]	O/E [%]
S 1	1:225	2466	-	100
	1:350	1296	1585	82
	1:1200	420	462	91
	1:2400	237	233	102
	1:4800	111	116	96
	Mean			<b>92.6</b>
S 2	1:100	758	-	100
	1:120	702	632	111
	1:150	552	505	109
	1:200	406	379	107
	1:400	210	190	111
	1:800	108	95	114
Mean			<b>111</b>	
S 3	Mean			<b>105</b>
S 4	Mean			<b>101</b>
S 5	Mean			<b>99</b>
<b>Mean All</b>				<b>102</b>

Table 25\* **Cross Reactivity**

Spiked with	MRP8 [ng/ml]	MRP14 [ng/ml]
100 µg/ml	26.0	38.7
10 µg/ml	8.0	3.4
1 µg/ml	<4.0	<4.0
100 ng/ml	<4.0	<4.0
10 ng/ml	<4.0	<4.0

Table 26\* **Clinical Study**

	Calprotectin (EK-CAL)	Lactoferrin
n	401	391
cut-off	50 µg/g	3 µg/ml
Sensitivity	84.4%	74.2%
Specificity	94.5%	91.0%
PPV	87.8%	79.3%
NPV	92.8%	88.4%
LR+	15.4	0.17
LR-	8.25	0.28

\* Data have been established with the lower range ELISA procedure

**Table description:** cf. "Results" (pg. 4, 6), "Performance Characteristics" (pg.4, 6) and "Interpretation of Results" (pg. Fehler! Textmarke nicht definiert.).

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
## APPENDIX III

### REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

1. Striz I and Trebichavsky I: *Calprotectin – a pleiotropic molecule in acute and chronic inflammation*. *Physiol Res.* **53**, 245-253 (2004)
2. Tibble JA et al.: *Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease*. *Gastroenterol* **123**, 450-460 (2002)
3. Fagerhol MK: *Calprotectin, a faecal marker of organic gastrointestinal abnormality*. *Lancet* **356**, 1783-4 (2000)
4. Hessian PA and Fisher L: *The heterodimeric complex of MRP-8 (S100A8) and MRP-14 (S100A9). Antibody recognition, epitope definition and the implications for structure*. *Eur J Biochem* **268**, 353-63 (2001).
5. Goebeler M et al.: *Expression and complex formation of S100-like proteins MRP8 and MRP14 by macrophages during renal allograft rejection*. *Transplantation* **58**, 355-61 (1994).
6. Tibble JA et al.: *A simple method for assessing intestinal inflammation in Crohn's disease*. *Gut* **47**, 506-513 (2000).
7. Wassell J et al.: *Faecal Calprotectin: a new marker for Crohn's disease?* *Ann Clin Biochem* **41**, 230-232 (2004)
8. Konikoff MR and Denson LA: *Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease*. *Inflamm Bowel Dis* **12**(6), 524-34 (2006)
9. Fagerberg UL et al.: *Colorectal inflammation is well predicted by fecal calprotectin in children with gastrointestinal symptoms*. *J Pediatr Gastroenterol Nutr* **40**, 450-5 (2005)
10. Johnson M W et al.: *Faecal calprotectin: a noninvasive diagnostic tool and marker of severity in pouchitis*. *Eur J Gastroenterol Hepatol* **20** :174-179, (2008)
11. Sipponen T et al. *Faecal Calprotectin, Lactoferrin, and Endoscopic Disease Activity in Monitoring Anti-TNF-alpha Therapy for Crohn's Disease*. *Aliment Pharmacol Ther* **28**, 1221–1229, (2008)
12. Tschanguizi et al. *Stool Calprotectin as a Marker of inflammation*. Poster at UEGW, 2007
13. Jahnsen J, Røseth AG, Aadland E. *Measurement of calprotectin in faeces*. *Tidsskr Nor Legeforen* **128**, 743–5 (2008)


## CALPROTECTIN EXTRACTION

### Standard Extraction Procedure

## CALPROTECTIN ELISA

### Precoated Microtiter Plate

➔ Read absorbance at 450 nm (within 30 minutes)

TIME TO RESULT: 75 MINUTES

**APPENDIX V**  
**SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SIMBOLOS**

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad
<b>REF</b>	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo
<b>LOT</b>	Batch code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i>
	Contains sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura
<b>BUFEX</b>	Extraction Buffer Extraktions-Puffer Tampon d'extraction Tampone per estrazione Tampón de extracción
<b>MP</b>	Microtiterplate Mikrotiter-Platte Microplaque Micropiastra Microplaca

Symbol	Explanation
<b>BUF WASH 10X</b>	Wash Buffer Concentrate (10x) Wasch-Puffer Konzentrat (10x) Concentré de tampon de lavage (10x) Tampone di lavaggio concentrato (10x) Tampón de lavado concentrado (x10)
<b>BUF INC</b>	Incubation Buffer Inkubations-Puffer Tampon d'incubation Tampone di incubazione Tampón de incubación
<b>CAL A - CAL E</b>	Calibrator A -E Kalibrator A -E Calibreur A -E Calibratore A - E Calibrador A - E
<b>CONTROL L</b>	Control Low Kontrolle tief Contrôle bas Controllo basso Control bajo
<b>CONTROL H</b>	Control High Kontrolle hoch Contrôle élevé Controllo alto Control alto
<b>EL</b>	Enzyme Label Enzym-Marker Marqueur enzymatique Marcato enzimatico Marcador enzimático
<b>SUBS TMB</b>	TMB Substrate TMB-Substrat Substrat TMB Substrato di TMB Substrato de TMB
<b>SOLN STOP</b>	Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de parada

Printing Date  
2012-02-01