# anti-IFN $\beta$ BAB

anti-Interferon beta Binding Antibodies  
ELISA

EK-IFNB 96 tests

Revision date: 2006-07-28

## ENGLISH

### INTENDED USE

The BÜHLMANN anti-IFN $\beta$  ELISA kit is intended for the direct and quantitative *in vitro* diagnostic determination of IgG antibodies to therapeutically administered Interferon- $\beta$  (IFN $\beta$ ) in human serum (1-10).

### PRINCIPLE OF THE ASSAY

Serum from Multiple Sclerosis (MS) patients treated with interferon- $\beta$  and suspected to contain antibodies (Ab) to the substance administered, calibrators and controls are incubated in microtiter wells coated with a mix of different IFN $\beta$  molecules (natural human IFN $\beta$ , IFN $\beta$ -1a and IFN $\beta$ -1b). After removal of unreacted material by washing, a horseradish peroxidase (HRP)-labelled antibody to human IgG is added to the wells.

After a second washing step, the TMB Substrate (tetramethylbenzidine) is added to the wells and color develops in proportion to the amount of anti-IFN $\beta$  binding antibodies bound in the initial step. The reaction is terminated by the addition of stop solution and the color absorbance is measured in a microtiter plate reader at a wavelength of 450nm.

### REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
<b>Microtiter Plate</b> Precoated with human-IFN $\beta$ , rec. IFN $\beta$ -1a & rec. IFN $\beta$ -1b	12 x 8 wells	B-IFNB-MP	Wash twice before use
<b>Plate Sealer</b>	3 pieces		
<b>Wash Buffer Concentrate (10x)</b> With preservatives	1 vial 100 ml	B-IFNB-WB	Dilute with 900 ml of deionized water
<b>Incubation Buffer</b> With preservatives	1 vial 100 ml	B-IFNB-IB	Ready to use
<b>Calibrators A to D<sup>1)</sup></b> pool of lyoph. human serum containing anti-IFN $\beta$ IgG	4 vials 1 ml	B-IFNB- CASET	Reconstitute each vial with 1 ml of Incubation Buffer; vortex
<b>Control Low / High<sup>2)</sup></b> Lyoph. human serum	2 vials 1 ml	B-IFNB- CONSET	Reconstitute each vial with 1 ml of Incubation buffer; vortex
<b>Enzyme Label</b> Anti-human IgG conjugated HRP in a protein-based buffer with preservatives	1 vial 11 ml	B-IFNB-ELG	Ready to use
<b>TMB Substrate</b> TMB in citrate buffer with H <sub>2</sub> O <sub>2</sub>	1 vial 11 ml	B-TMB	Ready to use
<b>Stop Solution</b> 0.25 M sulfuric acid	1 vial 11 ml	B-STTS	Ready to use <b>Corrosive agent</b>

Table 1

<sup>1)</sup> The Calibrators A, B, C and D contain 500, 200, 80 and 20 BTU of human anti-IFN $\beta$  IgG antibodies, respectively.

<sup>2)</sup> The Controls contain lot-specific amounts of human anti-IFN $\beta$  IgG antibodies. Refer to the additional QC Data Sheet for actual concentrations.

### STORAGE AND SHELF LIFE OF REAGENTS

Unopened Reagents	
Store at 2-8°C. Do not use past kit expiration date.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the foil pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store until expiration date at 2-8°C.
Wash Buffer	Store for up to 2 months at 2-8°C after dilution.
Calibrators	Store for up to 2 months at -20°C.
Controls	
Incubation Buffer	Store at 2-8°C until expiration date printed on the labels.
Enzyme Label	
TMB Substrate	
Stop Solution	Stable at 18-28°C until expiration date printed on the label.

Table 2

### WARNINGS AND PRECAUTIONS

The Microtiter Plate (B-IFNB-MB), Calibrators (B-IFNB-CASET) and Controls (B-IFNB-CONSET) of this kit contain components of human origin. Each serum donor unit used in the preparation of the kit components was tested by an FDA approved method and found negative for HBV surface antigen, so as for HCV and HIV1/2 antibodies. Although these methods are highly accurate, there is no guarantee that this material cannot transmit Hepatitis or AIDS. *Therefore, all patient specimens and kit components should be handled as if capable of transmitting infections.* All products containing human source material should be handled in accordance with good laboratory practice using appropriate precautions.

**Substrate and Stop Solution:** The Substrate Solution (B-TMB) contains tetramethylbenzidine (TMB), hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) and dimethylformamide. The Stop Solution (B-STTS) contains sulfuric acid. Each of those reagents is irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothing.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 20  $\mu$ l, 100  $\mu$ l and 1 ml pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 ml cylinder for the dilution of the Wash Buffer Concentrate.
- Microtiter plate washer or squeeze bottle for the Wash Buffer.
- Blotting paper.
- Refrigerator.
- Microtiter plate rotator.
- Microtiter plate reader for the measurement of absorbance at 450 nm.

### SPECIMEN COLLECTION AND STORAGE

The procedure calls for 20  $\mu$ l of serum per duplicate determination. Collect blood into plain tubes, avoid hemolysis, mix by inverting sample tube several times and leave to clot for 45 minutes at room temperature (18-28°C) protected from light. Centrifuge at 2000 x g for 15 minutes at room temperature (18-28°C) and collect the serum.

Lipemic, hemolytic and icteric samples should not be used in this assay. Lipemic samples can be avoided by asking patients to fast for at least 12 hours prior to the sample being taken.

Store the serum samples at -20°C if not assayed within the same day. Repeated freezing and thawing of samples may decrease antibody activity.

### PROCEDURAL NOTES

- The enzyme used as the label is inactivated by oxygen and is highly sensitive to sodium azide, thimerosal, hypochlorous acid and aromatic chlorohydrocarbons often found in laboratory water supplies. Therefore, use only deionized high quality water.
- The sample dilution (1:50) is considered in the concentration specification of the Calibrators. Therefore, the amount of BTU present in the unknown sample can be directly read by interpolation with the standard curve.
- If the initial concentration of an unknown sample is greater than the highest Calibrator, the serum sample should be further diluted with Incubation Buffer and assayed again according to the assay procedure. The additional dilution must be considered when calculating the actual concentration of BTU present in the unknown sample.
- If the microtiter plate reader is not capable of reading absorbance greater than 2 or greater than the absorbance

of the Calibrator A, a second reading at a wavelength of 490 or 492 nm is recommended (reference filter at 600 to 620 nm if available). In this case, proceed to construct a second standard curve with the absorbance readings of all calibrators at 490 or 492 nm. The concentration of the off-scale samples at 450 nm are then read from the new standard curve as described above. The readings at 490 or 492 nm should not replace the on-scale readings at 450 nm.

### ASSAY PROCEDURE

1. Dilute all patient samples 1:50 with Incubation Buffer (e.g. 20 µl of serum + 980 µl of Incubation Buffer). Allow diluted samples to set for 30 minutes at 18-28°C prior to pipetting in step 4d.
2. Prepare a plate with sufficient strips to test the desired number of Calibrators, Controls and samples. Remove excess strips from the holder and re-seal them in the plastic foil bag containing the desiccant bags **without delay**. Store refrigerated.
3. Wash the coated strips twice using at least 300 µl of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
- 4a. Pipet 100 µl of Incubation Buffer in duplicate into the wells A1+A2 as reagent blank.
- 4b. Pipet 100 µl of Standard A in duplicate into wells B1+B2. Pipet 100 µl of Standard B in duplicate into wells C1+C2. Pipet 100 µl of Standard C in duplicate into wells D1+D2. Pipet 100 µl of Standard D in duplicate into wells E1+E2.
- 4c. Pipet 100 µl of Low Control in duplicate into wells F1+F2. Pipet 100 µl of High Control in duplicate into wells G1+G2.
- 4d. Pipet 100 µl of each diluted sample in duplicate into the subsequent wells.
5. Cover with a Plate Sealer and incubate for 2 hours (± 5 minutes) at 2-8°C.
6. Remove and discard the Plate Sealer. Empty the wells and wash three times using at least 300 µl of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.
7. Add 100 µl of Enzyme Label to all wells.
8. Cover the plate with a Plate Sealer and incubate for 2 hours (± 5 minutes) at 2-8°C.
9. Remove and discard the Plate Sealer. Empty the wells and wash three times using at least 300 µl of Wash Buffer per well. Empty the wells and strike the plate firmly onto blotting paper.  
**Important:** Allow the TMB Substrate Solution to reach 18-28°C before use in step 10.
10. Add 100 µl of the TMB Substrate Solution to each well.
11. Cover the plate with a Plate Sealer, place the plate on a plate mixer set at 800-1000 rpm, protect the plate from direct light and incubate for 30 minutes (± 5 minutes) at 18-28°C.
12. Add 100 µl of Stop Solution to all wells. Remove air bubbles with a pipette tip. Proceed to step 13. within 30 minutes.
13. Read the absorbance at 450 nm in a microtiter plate reader.

### RESULTS AND CALCULATION

**Standard Curve:** Record the absorbance at 450 nm for each Calibrator and Blank (NSB) well. Average the duplicate values, subtract the average of the blank wells (NSB) and record the averages (= corrected average absorbance). Plot the corrected average absorbance (vertical axis) versus the Bühlmann Titer Units (BTU) of the standards (horizontal axis) using a lin/log graph paper. Draw the best fitting curve or calculate the standard curve using a smoothed spline fitting algorithm.

**Samples and Controls:** Record the absorbance at 450 nm for each sample and control well. Average the duplicate values, subtract the average of the blank wells and record the averages (= corrected average absorbance). Locate the corrected average absorbance value of the samples and controls on the vertical axis, draw a horizontal line intersecting the standard curve and read the Bühlmann Titer Units (BTU) from the horizontal axis.

**Standardization:** The standards of the anti-IFNβ BAB ELISA kit were calibrated against an internal reference

**Bühlmann Titer Units (BTU)** were established as follows:

- Normal donor samples as well as samples from MS patients before IFNβ treatment were assayed according to the anti-IFNβ BAB ELISA assay procedure (*cf.* table below).
- Serially diluted samples of the reference pool were assayed in the same run.
- The **dilution** at which the reference pool falls short of the cut-off of the control samples, mentioned in first step, corresponds to the **titer** of the reference pool, expressed in Bühlmann Titer Units (BTU).

NOTE: Serum titer values depend on the assay method and, in particular, on the specificity and the cut-off values established with a given assay method. Titer values obtained with different assay methods cannot be compared directly.

**Typical Data:** See Table 11 and Figure 1 for typical data (results and standard curve). *These results and standard curve are provided for demonstration purposes only. A standard curve must be generated for each set of samples to be assayed. The BÜHLMANN QC laboratory uses a smoothed spline curve fit.*

### QUALITY CONTROL

A thorough understanding of this instruction for use is necessary for the successful use of the product. Reliable results will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following this instruction for use.

Since there is no control serum for anti-IFNβ-BAB commercially available, we recommend using a positive serum pool for internal quality control. The confidence limits for the BÜHLMANN Controls are lot-specific and printed on the additional QC Data Sheet.

The reproducibility of standard curve parameters and control values should be within established limits of laboratory acceptability.

If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) ELISA reader settings iii) expiration dates of reagents iv) storage and incubation conditions v) TMB Substrate Solution should be colorless vi) purity of water.

### LIMITATIONS

- The reagents supplied with this kit are optimized to measure antibodies directed against injected interferon-β in human serum.
- Anti-IFNβ antibody titer value should be used as supplementary data available to the physician in monitoring IFNβ treatment.
- Serum titer values depend on the assay method and, in particular, on the specificity and the cut-off values established with an assay method. Titer values obtained with different assay methods cannot be compared directly.

## PERFORMANCE CHARACTERISTICS

**Intra-Assay Precision (Within-Run): 4.7%.** The intra-assay precision was calculated from the results of 35 and 48 pairs of values, respectively, in a single run. The values are presented in Table 12 as Bühlmann Titer Units (BTU) of anti-IFN $\beta$  BAB.

**Inter-Assay Precision (Run-to-Run): 12.9%.** The inter-assay precision was calculated from the results of 7 pairs of values in 21 different runs. The values are presented in Table 13 as Bühlmann Titer Units (BTU) of anti-IFN $\beta$  BAB.

**Dilution Linearity/Paralellism: 102.4%.** Six anti-IFN $\beta$  BAB positive serum samples containing different concentrations of IFN $\beta$  binding antibodies, were differentially diluted with Incubation Buffer and assayed according to the assay procedure. The values are presented in Table 14 as Bühlmann Titer Units (BTU) of anti-IFN $\beta$  BAB

**Analytical Sensitivity: 5 BTU:** Twenty duplicates of Incubation Buffer were assayed in a single run. Mean and standard deviation (SD) were calculated for the absorbance values. The minimum detectable dose of IFN $\beta$  BAB was calculated to be 4.92 BTU by adding two SD to the mean absorbance of the reagent blank (Incubation Buffer) and intersecting this value with the standard curve obtained in the same run.

**Specificity:** Three sets of experiments were performed to assess the specificity of the anti-IFN $\beta$  BAB ELISA:

1. Comparison to Immunodot Blot: Thirteen patient samples containing elevated titers of anti-IFN $\beta$  antibodies (106 to 784 BTU) showed in a Immunodot Blot corresponding low to high *Relative Titer Units*. Four patient samples showing results below the technical cut-off (6 to 35 BTU) were also negative in the Immunodot Blot.

2. Neutralization of anti-IFN $\beta$  BAB by pre-incubation with IFN $\beta$ : A strong, concentration-dependent inhibition of binding to the microtiter plates was observed, when anti-IFN $\beta$  BAB positive sera were pre-incubated with Incubation Buffer containing 0.1 to 10  $\mu$ g/ml of IFN $\beta$ -1a, IFN $\beta$ -1b, natural human IFN $\beta$  or mixes thereof. Different concentrations (0.1 – 10  $\mu$ g/ml) of INF $\gamma$  did not exhibit substantial effects on the binding of anti-IFN $\beta$  antibodies from treated MS patients to the microtiter plate.

3. Specificity of anti-IFN $\beta$  BAB binding: Sera with positive anti-IFN $\beta$  BAB were tested on plates coated under the same condition, but without INF $\beta$ . All results obtained were negative or significantly below the technical cut-off.

Patient sera containing different pathological autoantibodies (i.e. anti-nuclear-Ab, anti-Ganglioside-Ab and anti-MAG-Ab) as well as Circulating Immune Complex positive sera showed results (0 to 23 BTU) significantly below the technical cut-off.

**Method Comparison:** 35 serum samples of IFN $\beta$ -1a and IFN $\beta$ 1b treated patients were analyzed using the BÜHLMANN anti-IFN $\beta$  BAB ELISA and the MxA Stimulation Test, respectively. The presence of IFN $\beta$  induces the synthesis of the antiviral MxA protein in human cells. Neutralizing anti-IFN $\beta$  antibodies can decrease the effect of IFN $\beta$  and, therefore their specific activities can be determined by reduced MxA production (quantification by ELISA). All sera containing neutralizing antibodies were also positive in the anti-IFN $\beta$  BAB ELISA. Then there were no falsely negative results observed. The cut-off applied for the anti-IFN $\beta$  BAB positive/ negative discrimination was 50 BTU. The comparison results are displayed in Table 15 (K.M. Myhr, et al. 1999, Poster at Scand. Neurol. Meeting). A second study with 242 serum samples from 39 patients treated either with IFN $\beta$ -1a or IFN $\beta$ -1b were analyzed using

the BÜHLMANN anti-IFN $\beta$  BAB ELISA and the cytopathic effect assay (CPE), respectively. The presence of IFN $\beta$  protects human cells from virally mediated cell lysis. Neutralizing anti-IFN $\beta$  antibodies can decrease the effect of IFN $\beta$  and, therefore, their specific activities can be determined by increased cell lysis. The cut-off applied for the anti-IFN $\beta$  BAB positive/negative discrimination was 50 BTU. The comparison results are displayed in Table 16 (A. Bertolotto, Torino IT, pers. communication).

## EXPECTED VALUES AND CUT-OFF

The frequency of anti-IFN $\beta$  binding antibodies (anti-IFN $\beta$  BAB) in normal human serum was determined using blood samples from 200 asymptomatic volunteer blood donors (100 men and 100 women each at the age of 18-72 years) and from 44 clinically positive Multiple Sclerosis (MS) patients before the onset of IFN $\beta$  treatment, respectively. All samples were assayed according to the assay procedure and the results obtained are given in Table 17, expressed in Bühlmann Titer Units (BTU).

NOTES: These titer ranges should be used as guidelines only. It is recommended that each laboratory establishes its own expected ranges for its control and patient population, respectively.

**Proposed Technical Cut-Off:** Elimination of the four values  $>$  mean + 3 S.D. would result in a cut-off value (mean + 3 S.D.) of 43 BTU. For practical reasons, we recommend to use a theoretical **cut-off value of 50 BTU**. Similarly, the mean + 3 S.D. value calculated from 44 clinically diagnosed MS patients, but not under IFN $\beta$  therapy, would result in a theoretical cut-off value of 42 BTU which is highly corresponding to that of normal blood donors.

**Clinical Cut-Off:** Dependent on the clinical question we recommend to establish, so-called clinical cut-offs. 1) In order to a screening application of the anti IFN $\beta$  BAB ELISA in IFN $\beta$  treated MS patients a clinical cut-off of 30 BTU is the best choice to find virtually all potential NAB positive candidates (11). 2) To study a possible direct correlation between the BAB titer and its clinical consequences, an elevated cut-off is recommended (e.g. 80 BTU or higher).

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## APPENDIX I

### REFERENCES/ LITERATURREFERENZEN/ REFERENCES/ RIFERIMENTI/ REFERENCIAS

1. Hall GL, et al.: *Beta-interferon and multiple sclerosis*. Trends Neurosci **20**, 63-67 (1997).
2. Chofflon M: *Recombinant human interferon beta in relapsing-remitting multiple sclerosis: a review of the major clinical trials*. Eur J Neurol **7**, 369-380 (2000).
3. Vallbracht A, et al.: *Interferon-neutralizing antibodies in a patient treated with human fibroblast interferon*. Nature **289**, 496-497 (1981).
4. The IFNB Multiple Sclerosis Study Group: *Neutralizing antibodies during treatment of multiple sclerosis with interferon beta-1b*. Neurology **47**, 889-894 (1996).
5. Antonelli G, Dianzani F: *Development of antibodies to interferon beta in patients: technical and biological aspects*. Eur Cytokine Netw **10**, 413-422 (1999).
6. Pachner AR: *Anticytokine antibodies in beta interferon-treated MS patients and the need for testing: Plight of the practicing neurologist*. Neurology **49**, 647-650 (1997).
7. Medenica RD, et al.: *Interferon inhibitor factor predicting success of plasmapheresis in patients with multiple sclerosis*. J Clin Apheresis **9**, 216 (1994).
8. Deisenhammer F, et al. : *Bioavailability of interferon beta 1b in MS patients with and without neutralizing antibodies*. Neurology **52**, 1239-1243 (1999).
9. Paty DW, et al.: *Guidelines for physicians with patients on IFN $\beta$ -1b: The use of an assay for neutralizing antibodies (NAB)*. Neurology **47**, 865-866 (1996).
10. Havredaki M, Barona F: *Variations of interferon activators and/or inhibitors in human serum and their relationship to interferon therapy*. Jpn J Med Sci Biol **38**, 107 (1985).
11. Pachner et al.: *Antibodies to IFN- $\beta$  in multiple sclerosis patients: measurement by a commercially available kit*. Proceedings, 125<sup>th</sup> Annual Meeting, American Neurological Association, October 2000, Boston, MA, Abstract #230.

Table 11: **Example of Typical Data**

	Conc. (BTU)	Absorbance (OD)	Calc. Conc. (BTU)	CV Conc. (%)
Blank		0.085		
Blank		0.086		
<b>Blank Avg.</b>		<b>0.086</b>		
Std. A	500	1.978	494	
Std. A	500	2.006	506	
<b>Std. A Avg.</b>	<b>500</b>	<b>1.993</b>	<b>500</b>	<b>1.8</b>
Std. B	200	1.072	199	
Std. B	200	1.077	201	
<b>Std. B Avg.</b>	<b>200</b>	<b>1.075</b>	<b>200</b>	<b>0.4</b>
Std. C	80	0.560	78	
Std. C	80	0.577	82	
<b>Std. C Avg.</b>	<b>80</b>	<b>0.569</b>	<b>80</b>	<b>3.3</b>
Std. D	20	0.235	20	
Std. D	20	0.233	20	
<b>Std. D Avg.</b>	<b>20</b>	<b>0.234</b>	<b>20</b>	<b>0.8</b>
Ctrl. LOW		0.573	81	
Ctrl. LOW		0.573	81	
<b>Ctrl. L. Avg.</b>		<b>0.573</b>	<b>81</b>	<b>0.0</b>
Ctrl. HIGH		1.163	222	
Ctrl. HIGH		1.170	224	
<b>Ctrl. H. Avg.</b>		<b>1.166</b>	<b>223</b>	<b>0.6</b>
Sample 1		0.112	<20	
Sample 1		0.108	<20	
<b>Sam. 1 Avg.</b>		<b>0.110</b>	<b>&lt;20</b>	<b>1.6</b>
Sample 2		1.444	300	
Sample 2		1.467	307	
<b>Sam. 2 Avg.</b>		<b>1.455</b>	<b>303</b>	<b>1.6</b>

Table 12: **Intra-Assay Precision (Within-Run)**

Sample diluted 1:50	Pairs [n]	Mean [BTU]	SD [BTU]	CV [%]
serum 1	35	63	3.5	5.6
serum 2	48	122	4.2	3.4
serum 3	48	172	8.6	5.0
Mean				4.7

Table 13: **Inter-Assay Precision (Run-to-Run)**

Sample diluted 1:50	Mean [BTU]	SD [BTU]	CV [%]
serum 4	10	1.7	16.4
serum 5	11	1.3	11.1
serum 6	36	3.3	12.0
serum 7	57	5.9	10.8
serum 8	89	14.2	15.9
serum 9	146	16.4	12.1
serum 10	205	22.1	11.9
Mean			12.9

Figure 1: **Example of a Standard Curve**


Table 14: **Dilution Linearity/Paralellism**

Sample	Dilution Factor	Observed [BTU]	Expected [BTU]	Recovery O/E [%]
serum 11	1: 50	349	--	--
	1:100	195	175	112
	1:200	92	87	105
	1:400	42	44	96
serum 12	1: 50	306	--	--
	1:100	154	153	101
	1:200	68	77	89
	1:400	30	38	78
serum 13	1: 50	202	--	--
	1:100	100	101	99
	1:200	44	51	87
	1:400	22	25	87
serum 14	1: 50	154	--	--
	1:100	79	77	103
	1:200	35	39	91
	1:400	18	19	94
serum 15	1: 50	84	--	--
	1:100	55	42	130
	1:200	30	21	143
	1:400	18	10	168
serum 16	1: 50	65	--	--
	1:100	28	32	87
	1:200	15	16	92
	1:400	7	8	82
Mean				102.4

Table 15: **Method Comparison vs. MxA Stimulation**

		anti-IFN $\beta$ BAB ELISA	
		+	-
MxA Assay	+	5	0
	-	4	26

Table 16: **Method Comparison vs. CPE Test**

		anti-IFN $\beta$ BAB ELISA	
		+	-
CPE Assay	+	61	6
	-	73	102

Table 17: **Expected Values**

	Normal Donors	MS patients before IFN $\beta$ treatment
Total [n]	200	44
Range of Values [BTU]	3-81	3-73
S.D. [BTU]	12	11
Mean Value [BTU]	15	9
Median Value [BTU]	11	6
Mean + 3 S.D. [BTU]	52	42
n > Mean + 3 S.D.	4	1

**Table description:** cf. "Results and Calculation", "Performance Characteristics" (page 3) and "Expected Values and Cut-off" (page 4).

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