



# ACE kinetic

Procedure

KK-ACK


A Commitment to Diagnostics

## INTENDED USE

BÜHLMANN ACE kinetic is an *in vitro* diagnostic biochemical assay for the quantitative determination of angiotensin converting enzyme (ACE) activity in serum samples. The assay aids the assessment of disease activity in patients with sarcoidosis in conjunction with other clinical and laboratory findings.  
For laboratory use only.

## PRINCIPLE OF THE ASSAY

The assay is a quantitative enzymatic test which can be easily applied on clinical chemistry analyzers or run by manual method. ACE catalyzes the conversion of angiotensin I to angiotensin II. The enzyme also mediates the cleavage of the synthetic substrate FAPGG (N-[3-(2-furyl)acryloyl]-L-phenylalanyl-L-glycyl-L-glycine) into the amino acid derivative FAP and the dipeptide GG. The linear kinetic of this cleavage reaction is measured by recording the decrease in absorbance at 340 nm. The final ACE activity in U/L in the patient sample is determined using a calibration curve generated from the measured calibrator value (Ronca-Testoni, Clin Chem 1983; Bénéteau, Clin Chem 1986).



## ASSAY PROCEDURE

### Application notes / assay installation

Assay procedures for the ACE kinetic are established on several clinical chemistry analyzers. Validated application notes describing installation and analysis on specific instruments are available from BÜHLMANN upon request.

### Validated Applications

|                         |  |
|-------------------------|--|
| Roche cobas®            | c501/502<br>c701/702<br>c303 (pure)<br>c503 (pro)                |
| Abbott                  | Alinity c<br>Architect c-series                                  |
| Beckman                 | AU480/AU680<br>AU5800/DxC700AU                                   |
| Siemens                 | Atellica CH930<br>Advia 2400<br>Dimension Vista<br>Dimension EXL |
| The Binding Site        | Optilite   |
| ThermoFisher Scientific | Indiko   |
| IDS                     | iSYS   |
| Mindray                 | BS480  |

For other clinical chemistry analyzers please contact [support@buhlmannlabs.ch](mailto:support@buhlmannlabs.ch).

Manual procedure on microtiter plate is possible applying an MP Reader with 37°C incubation and plate shaking option and optical filter at 340 nm and 415 nm.

## Pre-Analytics

|                    |  |
|--------------------|--|
| Sample required:   | ~200 µL serum<br>Gel separator tubes (SST) can be used<br>Optionally, Li-Heparin and Citrate Plasma can be used<br>EDTA Plasma inhibits ACE activity |
| Sample collection: | Serum collection tubes without anti-coagulants   |
| Sample storage:    | at 2-8°C up to 30 days<br>at -20°C at least 6 months   |

## Special Equipment

Open clinical chemistry analyzer

## Kit Components

ACE kinetic is available in different package sizes.

|                      | KK-ACK    | KK-ACK2   | KK-ACK4   | KK-ACKX    |
|----------------------|-----------|-----------|-----------|------------|
| Tests                | 100       | 2 x 50    | 400       | 1200       |
| Substrate            | 1 x 26 mL | 2 x 13 mL | 4 x 26 mL | 3 x 100 mL |
| Calibrator           | 1 x 2 mL  | 2 x 2 mL  | 2 x 2 mL  | 3 x 2 mL   |
| Controls normal/high | 1 x 2 mL  | 2 x 2 mL  | 2 x 2 mL  | 3 x 2 mL   |





# ACE kinetic

Characteristics

KK-ACK

A Commitment to Diagnostics

## Assay Performance

Data obtained on the Roche cobas® c501. Refer to specific application notes for performance on other analyzers.

**Reproducibility** **6.3-9.1% CV**

3 instruments/lots x 5 days x 5 replicates (EP05-A3)

On cobas c501, c701 and AU480

**Repeatability** **0.8-3.0% CV**

**Within-Laboratory Precision** **1.7-3.7% CV**

20 days x 2 runs x 2 replicates (EP05-A3)

**Accuracy/Recovery** **92.0 – 112.8%**

6 samples spiked with 20.5 U/L (10% volume) run in 4 replicates

### Sample carry-over

No statistical significant carry-over (EP10-A2)

**Limit of Blank (LoB)** **4.3 U/L**

**Limit of Detection (LoD)** **6.8 U/L**

**Limit of Quantification (LoQ)** **11.3 U/L**

Classical approach EP17-A2; LoQ n=60; <20% CV

**Linearity range** **4.3-535 U/L**

Samples >150 U/L automatically re-run with reduced volume; acceptance ±4 U/L or ±10% (EP06-A)

**Security zone** **up to 541 U/L**

No limiting effects observed

## INTERFERING SUBSTANCES

Susceptibility to interfering substances assessed according EP07-A2. Bias exceeding 20% considered as interference.

### Oral pharmaceuticals

No interference detected:

- Aspirin 0.65 mg/mL
- Azathioprine 3.0 µg/mL
- Chlorambucil 7.2 µg/mL
- Cyclophosphamide 0.375 mg/mL
- Eprosartan 0.36 mg/mL
- Hydroxychloroquine up to 0.06 mg/mL
- Ibuprofen 0.5 mg/mL
- Losartan 0.09 mg/mL
- Methotrexate 2.0 µg/mL
- Prednisone 0.3 µg/mL

### Serum indices

Interference detected at concentrations:

- triglycerides 2.24 mg/mL
- conjugated bilirubin 0.06 mg/mL
- unconjugated bilirubin 0.047 mg/mL
- hemoglobin 1.19 mg/mL

No triglycerides interference observed after short centrifugation (10 min / 12'000 x g) and separation of lipid-containing supernatant.

This document is for information purpose only, before performing the assay please carefully refer/read the respective IFU available (<https://www.buhlmannlabs.ch/support/downloads/eifus/>).

## REFERENCE INTERVALS

### Adults

2.5<sup>th</sup> – 97.5<sup>th</sup> percentile from healthy participants in three independent studies in Switzerland (n=80, age: 20 – 70), Germany (n=159, age: 18 – 64, ref. 3) and USA (n=327, age: 16 – 77):

**20 – 70 U/L**

### Children

2.5<sup>th</sup> – 97.5<sup>th</sup> percentile from healthy pediatric participants in a single study in Germany (n=84, age: 0.5 – 18):

**33 – 112 U/L**

### Plasma samples

Samples from healthy blood donors collected into lithium-heparin and citrate tubes compared to serum samples from the same donors:

Li-Hep Plasma (n=38)  $y=0.9x + 2.5$  ;  $r=0.975$   
Mean bias: -1.1%

Citrate Plasma (n=44)  $y=0.8x + 1.7$  ;  $r=0.990$   
Mean bias: -10.8%

### Ordering codes:

- KK-ACK 100 tests
- KK-ACK2 2x50 tests
- KK-ACK4 400 tests
- KK-ACKX 1200 tests

CE 0123



BÜHLMANN Laboratories AG  
Baselstrasse 55  
4124 Schönenbuch  
Switzerland

Phone +41 61 487 12 12  
Fax orders +41 61 487 12 99  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch