Extraction Kit

for BÜHLMANN fCAL[®] tests

B-CAL-EX3
B-CAL-EX12

Release Date: 2019-12-10
Version A1

ENGLISH

INTENDED USE

The extraction buffer (B-CAL-EX) is a reagent designed for the extraction of calprotectin in human stool specimens using the BÜHLMANN Smart-Prep or manual extraction method for laboratory analysis with BÜHLMANN fCAL® tests.

REAGENTS SUPPLIED

Code	Reagents	Quantity	Preparation
B-CAL-EX3	Extraction Kit extraction buffer	3 bottles x 125 mL	Ready to use
B-CAL-EX12	Extraction Kit extraction buffer	12 bottles x 125 mL	Ready to use

Table 1

MATERIALS REQUIRED BUT NOT PROVIDED

- Smart-Prep (Code: B-CAL-RD)
- 10 µL disposable inoculation loops
- 15 mL polypropylene tubes with screw caps required for standard extraction procedure
- Multi tube vortex mixer
- Precision balance (10-150 mg)
- Micro centrifuge (≥3000 g)
- Plate shaker

REAGENT STORAGE AND STABILITY

Do not freeze reagents! Shelf life of unopened/ opened reagents at 2-8 °C: see expiry date on package label.

WARNINGS AND PRECAUTIONS

- Read the instructions carefully prior to the use of the buffer. Buffer performance will be adversely affected, if the reagent is incorrectly stored as detailed in this instruction for use.
- Components must not be used after the expiry date printed on the labels.
- Let the reagent adjust to room temperature and mix well (vortex) before use.

Extraction

- To receive quantitative results it is important to homogenize the stool sample completely with the extraction buffer.

SPECIMEN COLLECTION AND HANDLING

For the extraction procedure, less than 1 g of native stool specimen is required. Collect stool specimen into plain tubes.

Important: The specimen must be collected without any chemical or biological additives.

Specimen Transport

Stool specimens should be received for processing by the laboratory within 3 days of collection. Stool specimens may be shipped at room temperature or refrigerated.

Specimen Storage

Stool specimens should be refrigerated at 2-8 °C and extracted within 3 days of receipt at the laboratory. Do not store samples at elevated temperatures.

Extract Storage

Calprotectin in extracts obtained by manual extraction method is stable at 2-8 °C for ≤ 7 days or at -20 °C for 36 months.

STOOL SAMPLE EXTRACTION

With the BÜHLMANN Smart-Prep

Refer to the corresponding Smart-Prep (Code: B-CAL-RD) instruction for use (IFU).

With the manual extraction method ¹⁾

1. Label and weigh the empty polypropylene tube including the inoculation loop. Note the weight (tare).
2. Take out 50 to 100 mg of the stool sample by means of the inoculation loop and place it into the pre-weighted tube and weigh it again (gross weight). Avoid taking up dietary fibers present in the sample during the sampling process.
3. Calculate the net amount of sample, by subtracting tare from the gross weight, break the inoculation loop away and leave the lower part of the loop in the tube.
4. Add extraction buffer according to the formula
$$\frac{x \text{ mg stool} \times 49}{100} = y \text{ } \mu\text{L extraction buffer (e.g. } \frac{50 \text{ mg stool} \times 49}{100} = 2450 \text{ } \mu\text{L buffer)}$$
 to the tube and close the tube.
5. Extract the samples by
 - Vigorously vortex the extraction tube containing buffer and stool sample on a (multitube) vortexer (at highest speed) for 30 seconds.
 - Next, incubate the extraction tube for 25 +/-5 minutes on a plate shaker at ca. 400 rpm. The inoculation loop inside the tube serves to enhance agitation.
 - Again, vigorously vortex the extraction tube for 30 seconds.
6. Transfer 1.5 mL of the homogenate into a 2 mL Eppendorf tube and centrifuge in a microcentrifuge for 5 minutes at 3000 x g.
7. Decant the supernatant into a fresh labeled tube and continue with the assay procedure of the corresponding BÜHLMANN kit.

Important: The extraction buffer (B-CAL-EX) should only be used in combination with BÜHLMANN fecal calprotectin tests.

¹⁾ According to the IFU of the BÜHLMANN fCAL® ELISA and Quantum Blue® fCAL tests.

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