BÜHLMANN GanglioCombi™ *Light* ELISA

with enzyme labels IgG/IgM Mix, IgG and IgM

Detection of anti-ganglioside
auto-antibodies by ELISA
(GM1, GD1b, and GQ1b)

EK-GCL-S 96 wells

Revision date: 2017-09-21

ENGLISH

INTENDED USE

BÜHLMANN GanglioCombi™ *Light* ELISA, is an *in vitro* diagnostic test intended to detect antibodies against defined relevant neural antigens GM1, GD1b and GQ1b (IgG and IgM) in serum samples from patients with suspected peripheral neuropathies with an unknown etiology. It allows quantitative classification of results into titre categories and serves as an aid to diagnosis of neuropathies (ref. 1-7).

INTENDED APPLICATION

With regard to the 3 different enzyme labels, the device components allow three application options:

1. Testing with the IgG/IgM mix conjugate allows to screen for the presence of antibodies and indicate a possible autoimmune neuropathy.
2. Testing with individual IgG and/or IgM conjugates allows antibody isotype determination.
3. For laboratory work-up initial sample screening using the IgG/IgM enzyme label mix (option 1), may be followed by differentiation of mix-positive samples using individual IgG and IgM conjugates (option 2) after prior consultation with a clinician/referring neurologist.

PRINCIPLE OF THE ASSAY

BÜHLMANN GanglioCombi™ *Light* ELISA is based on the enzyme-immunometric assay technique. The wells of the provided microtiter plate are coated with gangliosides: GM1, GD1b and GQ1b.

Calibrator, controls, and patient sera are incubated in the microtiter wells and anti-ganglioside antibodies present in the samples bind to the immobilized gangliosides. After washing off unbound substances, the antibodies are detected with horseradish-peroxidase (HRP) labelled antibodies against human IgG and/or IgM. Following a second washing step in which unbound enzyme label is removed, a substrate solution containing tetramethylbenzidine (TMB) is added. A blue colour develops in proportion to the amount of antibodies bound to the immobilized gangliosides. Colour development is stopped by adding an acidic stop solution (diluted sulphuric acid) which turns the blue solution into yellow. The intensity of the colour is measured at 450 nm.

The measured absorbance is proportional to the titre of antibodies present in a given sample. The titres of antibodies are expressed as % ratios of the calibrator and can be assigned to titre categories (negative, grey zone, positive, strongly positive).

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Reconstitution
Microtiter Plate precoated with gangliosides	12 x 8 wells	B-GCL-MP	Ready to use
Plate Sealer	3 pieces		
Wash Buffer Concentrate (10X) with preservatives	1 bottle 100 mL	B-GCO-WB	Dilute with 900 mL of deionized water
Incubation Buffer with preservatives	1 bottle 100 mL	B-GCO-IB	Ready to use

Reagents	Quantity	Code	Reconstitution
Calibrator Lyophilized with preservatives	1 vial	B-GCO-CA	Add 1.5 mL of Incubation Buffer
Negative, Low and Medium Control Lyophilized with preservatives	3 vials	B-GCO-CONSET	Add 1.5 mL of Incubation Buffer
Enzyme Label IgG Anti-human IgG Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 mL	B-GCO-ELG	Ready to use
Enzyme Label IgM Anti-human IgM Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 mL	B-GCO-ELM	Ready to use
Enzyme Label Mix Anti-human IgG and IgM Ab conjugated to HRP in a protein-based buffer with preservatives	1 vial 11 mL	B-GCO-ELGM	Ready to use
TMB Substrate TMB in citrate buffer	1 vial 11 mL	B-TMB	Ready to use
Stop Solution 0.25 M sulfuric acid	1 vial 11 mL	B-ST5	Ready to use Corrosive agent

Table 1

STORAGE AND SHELF LIFE OF REAGENTS

Sealed / Unopened Reagents	
All sealed / unopened kit components are stable at 2-8 °C until the expiration date printed on the labels.	
Opened / Reconstituted Reagents	
Microtiter Plate	Return unused strips immediately to the aluminium pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store for up to 4 months at 2-8 °C.
Diluted Wash Buffer	Store for up to 4 months at 2-8 °C.
Calibrator	Store for up to 4 months at 2-8 °C.
Controls	Do not freeze!
Incubation Buffer	Store at 2-8 °C until expiration date printed on the labels.
Enzyme Labels	
TMB Substrate	
Stop Solution	Store at 18-28 °C until expiration date printed on the labels.

Table 2

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 20 µL, 100 µL and 1000 µL pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 mL cylinder for the reconstitution of the wash buffer.
- Squeeze bottle for wash buffer or automatic microtiter plate washer.
- Blotting paper.
- Orbital shaker for microtiter plates.
- Microtiter plate reader for measurement of absorbance at 450 nm.

PRECAUTIONS

Safety precautions

- Both, calibrator (B-GCO-CA) and controls (B-GCO-CONSET) of this kit contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.
- **Stop solution:** The stop solution (B-STs) contains sulfuric acid (0.25 M). The reagent is an irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothes. After contact with eyes or skin, wash immediately with plenty of water.
- **Reagents:** Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, irritation / burns can occur.
- Unused solution should be disposed of according to local state and federal regulations.

Technical precautions

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- **Residues in the microtiter plate wells** result from the production process. They are removed in the washing step (assay procedure step 3) and do not affect the results.
- **Prepare reagents before starting the assay procedure.** Reagents used in **steps 3-9** must be cold (2-8 °C) and kept cold while pipetting and washing. Put the TMB substrate at room temperature (18-28 °C).
- **Steps 3-9:** Use cold (2-8 °C) reagents for all these steps and keep them cold while pipetting. Recommendation: Prepare the wash buffer the evening before performing the assay and place it into the fridge overnight.
- **Wash steps 3, 6 and 9:** The wash steps are crucial for removing residues in the microtiter plate wells resulting from the production process (step 3) as well as any unbound antibodies (steps 6 and 9).
 - Always perform the wash steps with cold (2-8 °C) wash buffer.
 - Make sure that all wells are completely empty after the last washing cycle.
- **Step 9:** Adjust TMB substrate to room temperature (18-28 °C) before using it.
- **Step 11:** Shake the microtiter plates during the incubation with substrate. Depending on the orbital plate shaker, we recommend 400-600 rpm. The solution should move in the wells but must not spill over.
- If an automated washer is used, "plate mode" should be chosen so that dispensing is performed sequentially on all strips before aspirating.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.

- Every effort should be made to ensure that no cross contamination occurs between reagents, samples or between wells.
- Microwells cannot be re-used.

SPECIMEN COLLECTION AND STORAGE

- The procedure requires <0.1 mL of blood or <50 µL of serum, respectively.
- Refer to page 7 to learn about the interference of haemolyzed, lipemic or icteric samples.
- Collect blood into plain tubes (no anti-coagulant), avoid haemolysis, leave to clot for one hour, centrifuge for 10 minutes at approximately 1500 x g at room temperature (18-28 °C), collect the serum.
- We recommend freezing aliquots of patient samples if you need to store samples in order to avoid repeated freezing / thawing.
- Store serum samples at ≤-20 °C up to 4 months. For long-term storage we recommend -70 °C (samples are stable for >1 year). Frozen samples should be thawed and vortexed thoroughly prior to use.

ASSAY PROCEDURE

You can choose between three basic options:

- (1) Detection of IgG/IgM mix-isotypes: step 4a-4e and 7
- (2) Detection of IgG and IgM isotypes: step 4a'-4f' and 7'
- (3) Two step approach: Option 1 → autoimmune antibody positive samples → Option 2.

Note: Equilibrate TMB substrate to room temperature (18-28 °C)

1. **Dilute** all patient samples to be investigated 1:50 with incubation buffer. Use 30 µL of serum + 1470 µL (cold: 2-8 °C!) incubation buffer. Mix by vortexing and leave diluted samples as well as reconstituted calibrator and controls for 30 minutes at 2-8 °C prior to pipetting (refer to step 4a and 4b).
2. **Prepare a plate-frame** with the required number of strips to test the patient samples. Reseal the remaining strips in the foil pouch together with the desiccant packs **immediately**. Store refrigerated.

Note: Use cold reagents in steps 3 to 9.

3. Wash coated wells twice using at least 300 µL of cold! wash buffer per well. Empty wells and tap plate firmly onto blotting paper to remove remaining liquid completely.

Note: Immediately proceed to the next steps.

Option 1: Detection of IgG/IgM mix-isotypes

- 4a. Calibrator: Pipet 100 µL of calibrator into the well A1 (refer to figure 1A).
- 4b. Controls: Pipet 100 µL of medium control into well B1, of low control into well A2 and of negative control into well B2 (refer to figure 1A).

Note: If more than three strips per run are used, calibrator and controls can be tested in duplicates (see figure A).

- 4c. Patient serum: Pipet 100 µL of diluted patient serum 1 into wells C1-E1 (refer to figure 1A).

4d. Patient serum: Pipet 100 µL of diluted patient serum 2 into wells F1-H1 (refer to figure 1A)

4e. Pipet 100 µL of diluted patient sera 3-24 into subsequent wells (refer to figure 1A).

Option 2: Detection of IgG isotypes

4a'. Calibrator: Pipet 100 µL of calibrator into the well A1 (refer to figure 1B).

4b'. Controls: Pipet 100 µL of medium control into well B1, of low control into well A2 and of negative control into well B2 (refer to figure 1B).

Note: If more than three strips per isotype are used, calibrator and controls can be tested in duplicates (see figure 1B).

4c'. Patient serum: Pipet 100 µL of diluted patient serum 1 into wells C1-E1 (refer to figure 1B).

4d'. Patient serum: Pipet 100 µL of diluted patient serum 2 into wells C1-H1 (refer to figure 1B).

4e'. Pipet 100 µL of diluted patient sera 3-12 into subsequent wells.

Detection of IgM isotypes.

4f'. Repeat steps 4a'-4e' using subsequent wells or a new microtiter plate if necessary (refer to figure 1B).

For options 1 and 2: Sample incubation and washes

5. Cover the plate with a plate sealer and incubate for 2 hours ±5 minutes at 2-8 °C (do not shake the plate).

6. Remove plate sealer. Empty the wells and wash three times using at least 300 µL of cold wash buffer (2-8 °C) per well. Empty the wells and strike the plate firmly onto blotting paper in order to remove washing buffer completely.

For option 1: Detection of IgG/IgM mix-isotype

7. Add 100 µL of enzyme label IgG/IgM mix to the wells.

For option 2: Detection of IgG and IgM isotypes

7'. Add 100 µL of enzyme label IgG or IgM to the respective wells (refer to figure 1B).

For option 1 and 2: Incubation with enzyme labels, washes, detection

8. Cover the plate with a plate sealer and incubate for 2 hours ±5 minutes at 2-8 °C (do not shake the plate).

9. Remove plate sealer. Empty the wells and wash three times using at least 300 µL of cold wash buffer (2-8 °C) per well. Empty the wells and strike the plate firmly onto blotting paper.

Note: Adjust TMB substrate solution to room temperature (18-28 °C).

10. Add 100 µL of TMB substrate solution to each well.

11. Cover plate with a plate sealer, incubate plate on an orbital plate shaker at 400-600 rpm for 30 ±2 minutes at 18-28 °C. Protect the plate from direct light.

12. Add 100 µL of stop solution to all wells. Proceed to step 13 within 30 minutes.

13. Read absorbance at 450 nm in a microtiter plate reader.

QUALITY CONTROL

A good understanding of this instruction for use is necessary to obtain reliable results. These will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following the instruction for use. Since there is no control serum for anti-ganglioside antibodies commercially available, we recommend using a positive, and negative serum pool for internal quality control. A minimal OD value of 1.2 is recommended for the calibrator. All controls must be within the established expected ranges (% ratio). The expected ranges of the controls are lot-specific and indicated in the QC data sheet. Performance characteristics should be within established limits. If these characteristics are not in conformity with established limits and repetition excludes handling failures, check the following issues: i) Have all reagents, used in step 3-10, been kept at 2-8 °C? ii) accuracy of pipets, thermometers, and timers, iii) settings of ELISA washer and reader, iv) expiration date of the reagents v) storage and incubation conditions vi) colour of the TMB substrate solution (should be colourless) and vii) purity of the water.

STANDARDIZATION

The calibrator included in this kit has been calibrated against internal reference material. It has been adjusted to 100 % ratio.

RESULTS AND CALCULATION

Calculation of results:

1. Record absorbance (OD) at 450 nm for each well (calibrator, controls and patient samples).
2. Average the duplicate calibrator and control values (if available).
3. Results are expressed as ratio of absorbance of samples and the (averaged) absorbance of the calibrator.

IgG/IgM Mix isotypes

$$\% \text{ Ratio: } \frac{\text{absorbance of samples or controls}}{\text{absorbance of calibrator}} \times 200$$

IgG and IgM isotypes

$$\% \text{ Ratio: } \frac{\text{absorbance of samples or controls}}{\text{absorbance of calibrator}} \times 100$$

Programs to calculate results as % ratio are available on most microplate readers.

Note: Results presented in tables 7 and 8 are examples. Calibrator and controls must be used in each individual assay.

LIMITATIONS

1. It is recommended that grey zone and positive results obtained with enzyme label IgG/IgM mix be first discussed with a clinician/referring neurologist prior to further isotype determination with individual IgG and IgM enzyme labels.
2. Dominant auto-immune responses may be accompanied by cross-reactivity with other assayed gangliosides. The cross-reactivity will typically show high inter-assay variation and may be clinically non-relevant. The interpretation of results should therefore only be made together with an expert/specialist.
3. Due to poly-reactivity of auto-immune antibodies (ref. 6) and differences in geographical prevalence, assay results should only be used to support the clinical interpretation of the neuropathy by an expert/specialist in combination with the patient's clinical picture.
4. The BÜHLMANN GanglioCombi™ *Light* ELISA has not been validated for plasmapheresis samples.

INTERPRETATION OF RESULTS

Corresponding limitations, if applicable, are indicated by numbers in superscript.

	Isotype: IgG/IgM Mix		
Titre Category / Ratio (%)	<30	30-50	>50
GD1b	Negative	Refer to limitations ¹	Refer to limitations ¹
GQ1b			
GM1			

Table 3

	Isotype: IgG		
Titre Category / Ratio (%)	<30	30-50	>50
GD1b	Negative	Retest at a later time point	Positive
GQ1b			
GM1			

Table 4

	Isotype: IgM		
Titre Category / Ratio (%)	<30	30-50	>50
GD1b	Negative	Retest at a later time point	Positive
GQ1b			
GM1			

Table 5

REFERENCE INTERVALS AND CUT-OFF RATIOS

In co-operation with the institutions mentioned below, a cut-off of 50 % has been established. Values <30 % have been clearly classified as negative. The established titre categories are based on n = 100 blood donors¹ (adult men and women between 18 to 70 years of age) and n = 277 pathological samples². Sera were assayed for each of the 3 gangliosides (GD1b, GQ1b and GM1), according to the assay procedure. Results of normal blood donors are shown in table 16.

Guidelines for the use of cut-off and titre categories / ratios (%):

Clinical Interpretation	Titre Category / Ratio (%)
Negative	<30
Grey zone	30-50
Cut-off	50
Positive	>50-100
Strongly positive	>100

Table 6

¹) Collected at the blood donation centre, University Hospital of Basel

²) Received from Friedrich-Baur-Institute, Ludwig-Maximilians-University of Munich; Department of Neurology, University of Basel; Department of Neurology, University of Lyon.

PERFORMANCE CHARACTERISTICS

Intra-assay precision (within-run):

GD1b: 2.3-4.4 % CV; GQ1b: 2.6-5.2 % CV;

GM1: 1.4-5.7 % CV.

For each of the three gangliosides coated on the microtiter plate of the BÜHLMANN GanglioCombi™ *Light* ELISA test two anti-ganglioside positive sera were selected. The serum samples were assayed in twelve replicates in a single run for each of the enzyme labels: IgG/IgM Mix, IgG and IgM. The study was performed with one reagent lot. The results are summarized in table 9, 10 and 11.

Inter-assay precision (between-run):

GD1b: 8.0-13.2 % CV; GQ1b: 8.2-19.6 % CV;

GM1: 9.0-21.0 % CV.

For each of the three gangliosides coated on the microtiter plate of the BÜHLMANN GanglioCombi™ *Light* ELISA test two anti-ganglioside positive sera were selected. The serum samples were assayed in single replicates in 20 independent runs, with one run per day, for each of the enzyme labels: IgG/IgM Mix, IgG and IgM. The study was performed with one reagent lot. The results are summarized in table 12, 13 and 14.

Limit of Blank (LoB): ≤6.1 % ratio

The LoB was established according to the CLSI guideline EP17-A. Twelve blank replicates (incubation buffer) per ganglioside were assayed in a single run for all three enzyme labels: IgG/IgM Mix, IgG and IgM. The Limit of Blank (LoB), expressed as the % ratio to the calibrator absorbance, was calculated to be 6.1 % for enzyme label IgG/ IgM Mix, ≤3.1 % for enzyme label IgG and ≤2.6 % for enzyme label IgM. The highest LoB obtained with the three different enzyme labels was taken to determine the overall Limit of Blank (LoB). The study was performed with one reagent lot. The LoB was evaluated by using parametric analysis.

Limit of Detection (LoD): ≤8.1 % ratio

The LoD was established according to the CLSI guideline EP17-A. For each of the three gangliosides coated on the microtiter plate a single clinical sample representing low antibody concentration was measured in twelve replicates, in a single run for each of the three enzyme labels: IgG/IgM Mix, IgG and IgM. The LoD, expressed as the % ratio to the calibrator absorbance, was calculated to be ≤8.1 % for enzyme label IgG/ IgM Mix, ≤6.3 % for enzyme label IgG, and ≤3.0 % for enzyme label IgM. The highest LoD obtained with the three different enzyme labels was taken to determine the overall Limit of Detection (LoD). The study was performed with one reagent lot. The LoD was calculated using parametric analysis.

Functional sensitivity: ≤8.1 % ratio

A precision profile was generated from the results obtained in the inter-assay precision study described above, for two positive sera samples for each of the three gangliosides coated on the microtiter plate and detected with each of the enzyme labels: IgG/IgM Mix, IgG and IgM. The data points were fitted using a polynomial non-linear fit with a cubic function.

The 95 %-confidence intervals of the fit were determined and the functional sensitivity was determined as the intersection of the cubic fit with the 20 % CV acceptance criterion. The results are summarized in figure 2.

Linearity

The linear range of the BÜHLMANN GanglioCombi™ *Light* ELISA was determined according to the CLSI guideline EP06-A. Multiple sera with concentrations over the entire measuring range of the test, allowing the evaluation of the majority of ganglioside / enzyme label combinations, were used. Serum samples were diluted according to the instruction for use. In case of high positive serum samples a higher dilution of 1:2000 was applied. Subsequently dilution series from each sample were prepared in graduations of 10 % using incubation buffer as the diluent. Linearity was defined as the interval in which the relative difference between the linear and, if significant, higher order polynomial fit was below 20 %. For ratios ≤25 % an absolute difference of below 5 % ratio was allowed. A linear range covering and exceeding the BÜHLMANN GanglioCombi™ *Light* ELISA titre categories of <30 % - negative; 30-50 % - grey-zone; 50-100 % positive and >100 % strongly positive was confirmed. The results are summarized in table 15.

INTERFERING SUBSTANCES

No interference is detected with the following substances up to the following concentrations: unconjugated bilirubin (icteric sera): 40 mg/d; conjugated bilirubin (icteric sera): 60 mg/dL; haemoglobin (haemolyzed sera): 400 mg/dL; haemolyzate (haemolyzed blood): 400 mg/dL and Triglycerides (Intralipid®): 3000 mg/dL.

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Microtiter plate set-up


Figure 1A: IgG-Mix conjugate: 24 patients / kit (1 MP / kit)


Figure 1B: IgG & IgM conjugate: 2 profiles / patient, 12 patients / Kit (1 MP/Kit)

Example of results

A IgG/IgM mix label

B-GCO-ELGM	Absorbance (OD450)	Ratio [%]	Result / Suggestion
Calibrator	1.415 1.445		
Calibrator Avg.	1.430	200	
Medium Control	0.498 0.482	69 67	
Med. Control Avg.	0.490	68	
Low Control	0.195 0.191	27 26	
Low Control Avg.	0.193	27	
Negative Control	0.090 0.100	12 14	
Neg. Control Avg.	0.095	13	
Sample 1 GM1	0.544	76	Positive
Sample 1 GD1b	0.745	104	Strongly positive
Sample 1 GQ1b	0.090	13	Negative

Table 7

B IgG & IgM labels

Enzyme labels	Absorbance (OD450)		Ratio [%]		Category	
	IgG	IgM	IgG	IgM	IgG	IgM
B-GCO-ELG/ B-GCO-ELM						
Calibrator	1.789 1.833	2.576 2.527				
Calibrator Avg.	1.836	2.551	100	100		
Medium Control	1.267 1.237	1.743 1.764	69 67	68 69		
Med. Control Avg.	1.252	1.753	68	69		
Low Control	0.567 0.584	0.938 0.942	30 32	37 37		
Low Control Avg.	0.571	0.940	31	37		
Neg. Control	0.061 0.051	0.098 0.095	3 3	4 4		
Neg. Control Avg.	0.056	0.097	3	4		
Sample 1 GM1	0.171	3.814	9	150	Neg.	Strong. Pos.
Sample 1 GD1b	1.021	0.354	56	14	Pos.	Neg.
Sample 1 GQ1b	0.378	0.208	21	8	Neg.	Neg.

Table 8

Intra-assay precision (within-run)

Enzyme label IgG/ IgM						
Gangliosides	GD1b	GD1b	GQ1b	GQ1b	GM1	GM1
Sample	1	2	1	2	1	2
Mean [Ratio]	110	310	286	84	86	156
SD [Ratio]	4.0	8.0	9.2	3.8	2.4	5.0
CV [%]	3.6	2.6	3.2	4.5	2.8	3.2

Table 9

Enzyme label IgG						
Gangliosides	GD1b	GD1b	GQ1b	GQ1b	GM1	GM1
Sample	1	2	1	2	1	2
Mean [Ratio]	139	111	89	222	59	82
SD [Ratio]	4.0	2.8	4.5	5.8	1.6	4.7
CV [%]	2.9	2.5	5.1	2.6	2.7	5.7

Table 10

Enzyme label IgM						
Gangliosides	GD1b	GD1b	GQ1b	GQ1b	GM1	GM1
Sample	1	2	1	2	1	2
Mean [Ratio]	56	142	72	74	79	87
SD [Ratio]	2.5	3.3	3.8	2.6	1.8	1.2
CV [%]	4.4	2.3	5.2	3.5	2.3	1.4

Table 11

APPENDIX I

TABLES AND FIGURES / TABELLEN UND ABBILDUNGEN / TABLEUX ET FIGURES /
TABELLE E FIGURE / TABLAS Y FIGURAS / TABLAS E FIGURAS

Inter-assay precision (between-run)

Enzyme label IgG/ IgM						
Gangliosides	GD1b	GD1b	GQ1b	GQ1b	GM1	GM1
Sample	1	2	1	2	1	2
Mean [Ratio]	96	306	92	262	74	184
SD [Ratio]	9.6	26.8	7.6	26.6	7.0	27.8
CV [%]	10.1	8.8	8.2	10.1	9.5	15.2

Table 12

Enzyme label IgG						
Gangliosides	GD1b	GD1b	GQ1b	GQ1b	GM1	GM1
Sample	1	2	1	2	1	2
Mean [Ratio]	109	140	42	204	52	84
SD [Ratio]	14.2	18.5	8.3	23.6	6.4	17.6
CV [%]	13.0	13.2	19.6	11.6	12.4	21.0

Table 13

Enzyme label IgM						
Gangliosides	GD1b	GD1b	GQ1b	GQ1b	GM1	GM1
Sample	1	2	1	2	1	2
Mean [Ratio]	154	46	171	63	102	63
SD [Ratio]	19.1	3.7	16.2	9.9	10.9	5.7
CV [%]	12.4	8.0	9.5	15.6	10.7	9.0

Table 14

Functional sensitivity


Figure 1

Linearity

Ganglio-side	Enzyme label	Serum	Range measured	Linear range
GD1b	IgG	581-SP	5-99	5-99
	IgM	72 BR	14-163	14-148
	Mix	581-SP	18-264 12-198	18-176 34-166
GQ1b	IgG	MG	18-228	18-228
GM1	IgG	SP69	7-94	7-94
		P II	9-181	17-175
		JJ	6-127	6-127
	IgM	SP69	10-93	21-79
		SM1	14-148	26-114
			9-138	9-89
	Mix	P II	16-302	16-302
		SP69	14-218	22-178

Table 15

Normal values

BÜHLMANN GanglioCombi™ Light ELISA: Anti-ganglioside antibodies (%)									
Enzyme label isotype	IgG			IgM			IgG/ IgM		
	GD1b	GQ1b	GM1	GD1b	GQ1b	GM1	GD1b	GQ1b	GM1
Negative	92.0	98.0	96.0	99.0	99.0	95.0	97.0	100	96.0
Greyzone	6.0	1.0	2.0	0.0	1.0	3.0	3.0	0.0	3.0
Positive	2.0	1.0	2.0	1.0	0.0	2.0	0.0	0.0	1.0
Strongly positive	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0


Table 16:

Distribution frequency of individual anti-ganglioside auto-antibodies (%) in normal blood donors (n=100) detected with BÜHLMANN GanglioCombi™ Light ELISA assay, classified into titre categories.

REFERENCES / REFERENZEN / REFERENCES / RIFERIMENTI / REFERENCIAS / REFERÊNCIAS

1. Willison HJ and Yuki N: *Peripheral neuropathies and anti-glycolipid antibodies*. Brain 125, 2591-2625 (2002)
2. Latov N: *Diagnosis and treatment of chronic acquired demyelinating polyneuropathies*. Nature Reviews Neurology 10, 435-446 (2014)
3. Burns TM and Mauermann ML: *The Evaluation of Polyneuropathies*. Neurology Clinical Practice 76, 6-13 (2011)
4. Stork ACJ et al. *Prevalence, specificity and functionality of anti-ganglioside antibodies in neuropathy associated with IgM monoclonal gammopathy*. Journal of Neuroimmunology 268, 89-94 (2014)
5. Humbel RL and Schmit P: *Anticorps antigangliosides et neuropathies peripheriques*. Rev Med Liege 51, 368-375 (1996)
6. Bourque P R et al. *Autoimmune peripheral neuropathies*. Clin Chim Acta 449, 37-42 (2015)
7. Delmont E and Willison H: *Diagnostic Utility of Auto Antibodies in Inflammatory Nerve Disorders*; J of Neuromuscular Diseases 2, 107-112 (2015)

BÜHLMANN GanglioCombi™ Light ELISA

TIME TO RESULT: 4.5 HOURS

APPENDIX IV

NOTES / NOTIZEN / NOTES / NOTE / NOTAS

APPENDIX V

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLI / SIMBOLOS

Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad Data de expiração
REF	Catalogue number Bestellnummer Référéce du catalogue Numero di catalogo Número de catálogo Código de encomenda
LOT	Batch code Lot-Bezeichnung Code du lot Codice del lotto Codigo de lote Código lote
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i> Producto sánitario para diagnóstico <i>in vitro</i>
	Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperature Límite de temperatura
MP	Microtiterplate Mikrotiterplatte Plaque de microtitration Micropiastra Placa de microtitulación Placa de microtitulação
BUF WASH 10X	Wash Buffer Concentrate (10x) Waschpufferkonzentrat (10x) Concentré de tampon de lavage Tampón de lavado concentrado (x10) Tampón de lavado concentrado (x10) Concentrado de tampão de lavagem (x10)
BUF INC	Incubation Buffer Inkubationspuffer Tampon d'incubation Tampón de incubación Tampão de incubação

Symbol	Explanation
CONTROL -	Negative Control Negativkontrolle Contrôle négatif Control negativo Control negativo Controle negativo
CONTROL L	Low Control Kontrolle tief Contrôle faible Control bajo Control bajo Control baixo
CONTROL M	Medium Control Kontrolle medium Contrôle moyen Control medio Control medio Control medio
CAL	Calibrator Kalibrator Calibreur Calibratore Calibrador Calibrador
EL IgG	Enzyme Label IgG Enzymmarker-IgG Marqueur enzymatique IgG Marcatore enzimatico IgG Marcador enzimático IgG Marcador enzimático de IgG
EL IgM	Enzyme Label IgM Enzymmarker-IgM Marqueur enzymatique IgM Marcatore enzimatico IgM Marcador enzimático IgM Marcador enzimático de IgM
EL MIX	Enzyme Label IgG/IgM Mix Enzymmarker- IgG/IgM Mix Marqueur enzymatique IgG/IgM mix Marcatore enzimatico Mix Mezcla de Marcador enzimático Marcador enzimático combinados
SUBS TMB	TMB Substrate TMB Substrat Substrat TMB Substrato di TMB Substrato TMB Substrato TMB
SOLN STOP	Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de interrupción Solução stop