


# BÜHLMANN GanglioCombi® MAG ELISA

## Pipetting Template



## Pipetting Guidance

Precoated Microtiter Plate



wash 2x

100 µL of Calibrator, Controls or Serum Samples (1:50)



incubate 2 hours ± 5 min at 2-8°C  
wash 3x (use cold wash buffer from refrigerator)

add 100 µL of Enzyme Label(s)



incubate 2 hours ± 5 min at 2-8°C  
wash 3x (use cold wash buffer from refrigerator)

add 100 µL of TMB Substrate  
(ambient temperature)




incubate 30 min ± 2 min at 18-28°C

add 100 µL of Stop Solution

Read absorbance at 450 nm (within 30 minutes)

# Autoimmune Peripheral Neuropathy: four easy Steps to complete the diagnostic assessment in your lab



\*„MAG“ is a chemically synthesized, sulfated disaccharide which mimics a MAG (HNK-1) epitope.



BÜHLMANN Laboratories AG  
Baselstrasse 55  
4124 Schönenbuch  
Switzerland

Phone +41 61 487 12 12  
Fax +41 61 487 12 34  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch

BÜHLMANN GanglioCombi® MAG ELISA  
anti-MAG Autoantibodies ELISA  
κ/λ Reagent Set  
anti-SGPG Autoantibodies ELISA  
BÜHLMANN GanglioCombi® Light ELISA  
anti-GM1 Autoantibodies ELISA

EK-GCM  
EK-MAG  
B-KLSET  
EK-SGPG  
EK-GCL-S  
EK-GM1-GM

BÜHLMANN GanglioCombi® is a registered trademark of BÜHLMANN in Switzerland.