anti-MAG Autoantibodies ELISA

MAG = Myelin Associated Glycoprotein

EK-MAG 96 tests

Release Date: 2018-10-26
Version A1

ENGLISH

INTENDED USE

The anti-MAG Autoantibodies ELISA is intended for the quantitative *in vitro* diagnostic determination of human IgM-auto-antibodies directed against Myelin Associated Glycoprotein (MAG) in serum. It serves as an aid to diagnosis of Paraproteinemic Demyelinating Neuropathy (PDN) in conjunction with other clinical and laboratory findings (1-7).

PRINCIPLE OF THE ASSAY

The anti-MAG Autoantibodies ELISA employs the quantitative enzymatically amplified sandwich-type immunoassay technique. The microtiter plates of the test are precoated with highly purified MAG from human brain (1). Calibrator, controls, and patient sera are incubated in the microtiter wells and any anti-MAG auto-antibodies present in the samples bind to the immobilized human MAG. After washing off of unbound substances, the anti-MAG auto-antibodies are detected with horseradish-peroxidase (HRP) labelled antibodies against human IgM. Following a second washing step, in which unbound enzyme label is removed, a substrate solution containing tetramethyl-benzidine (TMB) is added. A blue color develops in proportion to the amount of antibodies bound to the immobilized human MAG. Color development is stopped by adding an acidic stop solution (diluted sulphuric acid) which turns the blue solution into yellow. The intensity of the color is measured at 450 nm.

The measured absorbance is proportional to the titre of anti-human MAG auto-antibodies present in a given sample. A set of human anti-MAG auto-antibody calibrators is used to plot a calibration curve of absorbance versus human anti-MAG auto-antibody titer units. The concentrations of anti-human MAG auto-antibodies in the samples is calculated based on the calibration curve.

REAGENTS SUPPLIED AND PREPARATION

| Reagents | Quantity | Code | Reconstitution |
|---|---------------------------------|--------------|---------------------------------------|
| Microtiter Plate 96 wells precoated with human MAG | 12 x 8-wells strips with holder | B-MAG-MP | Ready to use |
| Plate Sealer | 3 pieces | | |
| Wash Buffer Concentrate (10x) | 1 bottle 100 mL | B-MAG-WB | Dilute with 900 mL of deionized water |
| Incubation Buffer with preservatives | 1 bottle 100 mL | B-MAG-IB | Ready to use |
| Calibrators A to D¹ Human serum with preservatives | 4 vials | B-MAG-CASET | Add 1 mL of Incubation Buffer |
| Low and High Control² Human serum with preservatives | 2 vials | B-MAG-CONSET | Add 1 mL of Incubation Buffer |
| Enzyme Label IgM Anti-human IgM antibody conjugated to HRP in a protein-based buffer with preservatives | 1 vial 11 mL | B-MAG-ELM | Ready to use Blue solution |

| Reagents | Quantity | Code | Reconstitution |
|---|-----------------|--------|--|
| TMB Substrate TMB in Citrate buffer | 1 vial 11 mL | B-TMB | Ready to use |
| Stop Solution 0.25 M Sulfuric acid | 1 vial 11 mL | B-STTS | Ready to use Corrosive agent |

Table 1

¹ After reconstitution, calibrators A, B, C and D contain 70000, 15000, 3000 and 1000 BÜHLMANN Titer Units (BTU) of anti-MAG antibodies, respectively.

² Controls contain lot-specific amounts of anti-MAG antibodies. Refer to the additional QC data sheet for exact concentrations.

STORAGE AND SHELF LIFE OF REAGENTS

| Sealed / Unopened Reagents | |
|--|---|
| All sealed/unopened kit components are stable at 2-8 °C until the expiration date printed on the labels. | |
| Opened / Reconstituted Reagents | |
| Microtiter Plate | Return unused strips immediately to the foil pouch containing the desiccant packs and reseal along the entire edge of zip-seal. Store for up to 2 months at 2-8 °C. |
| Diluted Wash Buffer | Store for up to 2 months at 2-8 °C. |
| Calibrators | Store for up to 2 months at -20 °C. |
| Controls | |
| Incubation Buffer | Store at 2-8 °C until expiration date printed on the labels. |
| Enzyme Label IgM | |
| TMB Substrate | |
| Stop Solution | Store at 18-28 °C until expiration date printed on the labels. |

Table 2

MATERIALS REQUIRED BUT NOT PROVIDED

- Precision pipettes with disposable tips: 20 µL, 100 µL and 1000 µL pipettes.
- Disposable polystyrene or polypropylene tubes for the preparation of sample dilutions.
- 1000 mL cylinder for the dilution of the wash buffer.
- Squeeze bottle for wash buffer or automatic microtiter plate washer.
- Blotting paper.
- Orbital shaker for microtiter plates.
- Microtiter plate reader for measurement of absorbance at 450 nm.

PRECAUTIONS

Safety precautions

- The microtiter-plate (B-MAG-MP), calibrators (B-MAG-CASET) and controls (B-MAG-CONSET) of this kit contain components of human origin. Although tested and found negative for HBV surface antigen, HCV and HIV1/2 antibodies, the reagents should be handled as if capable of transmitting infections and should be handled in accordance with good laboratory practices using appropriate precautions.
- **Stop solution:** The stop solution (B-STTS) contains sulfuric acid (0.25 M). The reagent is an irritant to eyes, skin and mucous membranes. Avoid contact with eyes, skin and clothes. After contact with eyes or skin, wash immediately with plenty of water.

- **Reagents:** Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, irritation / burns can occur.
- Unused solution should be disposed of according to local state and federal regulations.

Technical precautions

- Read carefully the instructions prior to carrying out the test. Test performance will be adversely affected, if reagents are incorrectly diluted, modified or stored under conditions other than those as detailed in this instruction for use.
- **Residues in the microtiter plate wells** result from the production process. They are removed in the washing step (assay procedure step 3) and do not affect the results.
- **Prepare reagents before starting the assay procedure.** Reagents used in **steps 3-9** must be cold (2-8 °C) and kept cold while pipetting and washing. Put the TMB substrate at room temperature (18-28 °C).
- **Steps 3-9:** Use cold (2-8 °C) reagents for all these steps and keep them cold while pipetting. Recommendation: Prepare the wash buffer the evening before performing the assay and place it into the fridge overnight.
- **Wash steps 3, 6 and 9:** The wash steps are crucial for removing residues in the microtiter plate wells resulting from the production process (step 3) as well as any unbound antibodies (steps 6 and 9).
 - Always perform the wash steps with cold (2-8 °C) wash buffer.
 - Make sure that all wells are completely empty after the last washing cycle.
- **Step 10:** Adjust TMB substrate to room temperature (18-28 °C) before using it.
- **Step 11:** Shake the microtiter plates during the incubation with substrate. Depending on the plate shaker, we recommend 400-600 rpm. The solution should move in the wells but must not spill over.
- If an automated washer is used, “plate mode” should be chosen so that dispensing is performed sequentially on all strips before aspirating.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Every effort should be made to ensure that no cross contamination occurs between reagents, samples or between wells.
- Microwells cannot be re-used.

SPECIMEN COLLECTION AND STORAGE

- The procedure requires <0.1 mL of blood or <50 µL of serum, respectively.
- Lipemic, hemolytic and icteric samples should not be used in this assay. Lipemic samples can be avoided by asking patients to fast for at least 12 hours prior to the sample being taken.
- Collect blood into plain tubes (no anti-coagulant), avoid hemolysis, leave to clot for one hour at room temperature

(18-28 °C), centrifuge for 10-30 minutes at 1000-2000 x g, collect the serum.

- Store serum samples at ≤-20 °C. Samples are stable for ≥1 year if stored at ≤-20 °C. Avoid repeated freeze-thaw cycles. Frozen samples should be thawed and mixed thoroughly by gentle swirling or inversion prior to use.
- We recommend freezing aliquots of patient samples in order to avoid repeated freezing/thawing.

ASSAY PROCEDURE

1. Dilute all patient samples to be investigated 1:1000 with incubation buffer. Use 2 µL of serum + 2000 µL of incubation buffer. Mix thoroughly by vortexing and leave diluted samples for one hour at 18-28 °C. Put samples as well as reconstituted calibrators and controls for 10 minutes on ice prior to pipetting in step 4.
2. **Prepare a plate-frame** with the required number of strips to test the calibrators, controls and patient samples. Reseal the remaining strips in the foil pouch together with the desiccant packs **immediately**. Store refrigerated.

Note: Use cold reagents in steps 3 to 9.

3. Wash coated wells four times using at least 300 µL of cold! wash buffer per well. Empty wells and tap plate firmly onto blotting paper to remove remaining liquid completely.
- 4a. Pipet 100 µL of incubation buffer (blank) in duplicate into wells A1+A2.
 - Pipet 100 µL of calibrator A in duplicate into wells B1+B2.
 - Pipet 100 µL of calibrator B in duplicate into wells C1+C2.
 - Pipet 100 µL of calibrator C in duplicate into wells D1+D2.
 - Pipet 100 µL of calibrator D in duplicate into wells E1+E2.
- 4b. Pipet 100 µL of the low control in duplicate into wells F1+F2.
 - Pipet 100 µL of the high control in duplicate into wells G1+G2.
- 4c. Pipet 100 µL of each diluted sample in duplicate into the subsequent wells.
5. Cover the plate with a plate sealer and incubate for 2 hours (±5 min) at 2-8 °C.
6. Remove plate sealer. Empty the wells and wash four times using at least 300 µL of cold wash buffer (2-8 °C) per well. Empty the wells and strike the plate firmly onto blotting paper in order to remove wash buffer completely.
7. Add 100 µL of enzyme label IgM to all wells.
8. Cover plate with a plate sealer and incubate for 2 hours (±5 min) at 2-8 °C.
9. Remove plate sealer. Empty the wells and wash four times using at least 300 µL of cold wash buffer (2-8 °C) per well. Empty the wells and strike the plate firmly onto blotting paper.

Note: Adjust TMB substrate solution to room temperature (18-28 °C).

10. Add 100 µL of TMB substrate solution to each well.
11. Cover plate with a plate sealer, incubate plate on an orbital plate shaker at 400-600 rpm for 30 ± 2 minutes at 18-28 °C. Protect the plate from direct light.
12. Add 100 µL of **stop solution** to all wells. Proceed to step 13 within 30 minutes.

13. Read the absorbance at 450 nm in a microtiter plate reader.

QUALITY CONTROL

A good understanding of this instruction for use is necessary to obtain reliable results. These will be obtained only by using precise laboratory techniques (current GLP guidelines) and accurately following the instruction for use.

Since there is no control serum for anti-MAG antibodies commercially available, we recommend using a positive, and negative serum pool for internal quality control.

All controls must be within established confidence ranges (BTU). The confidence ranges of the controls are lot-specific and printed on the QC data sheet delivered with this kit.

Performance characteristics should be within established limits. If these characteristics are not in conformity with established limits and repetition excludes handling failures, check the following issues: i) Have all reagents, used in step 3-10, been kept at 2-8 °C? ii) accuracy of pipets, thermometers, and timers, iii) settings of ELISA washer and reader, iv) expiration date of the reagents v) storage and incubation conditions vi) color of the TMB substrate solution (should be colorless) vii) purity of the water.

STANDARDIZATION

The calibrators included in this kit have been standardized against an internal reference pool. The reference pool consists of more than ten human sera containing low, medium and high titers of anti-MAG antibodies, respectively. Titers and specificity of each serum in the reference pool were first analyzed by anti-MAG immunoblots.

BÜHLMANN Titer Units (BTU) were established as follows:

- Normal donor samples were assayed according to the anti-MAG Autoantibodies ELISA assay procedure.
- Serially diluted samples of the reference pool were assayed in the same run.
- The dilution at which the reference pool sample falls short of the cut-off corresponds to the titer of the reference pool.

Note: Serum titer values depend on the assay method and, in particular, on the specificity and the cut-off values established with an individual assay method. Therefore, titer values obtained with different methods cannot be compared directly.

RESULTS AND CALCULATION

Calibration curve

Record the absorbance at 450 nm for each well to construct a calibration curve from the average of the calibrator duplicates, subtracting the average blank OD from each calibrator. It is recommended to use a software program to calculate the calibration curve and to determine the concentration of the samples, using a 4 parameter logistic (4 PL) fit.

Samples and controls

- Record the absorbance at 450 nm for each sample and control well. Subtract the average blank OD from the average of duplicate values for each sample and control.
- Locate the absorbance values of the samples and controls on the vertical axis, draw a horizontal line

intersecting the calibration curve and read from the horizontal axis.

If the microtiter plate reader is not capable of reading absorbance greater than 2 or greater than the absorbance of the highest calibrator (calibrator A), a second reading at a wavelength of 490 or 492 nm is recommended (reference filter at 600 or 620 nm if available). In this case, proceed to construct a second calibration curve with the absorbance readings of all calibrators at 490 or 492 nm. The concentration of the off-scale samples at 450 nm are then read from the new calibration curve as described above. The readings at 490 or 492 nm should not replace the on-scale readings at 450 nm.

Note: The results presented in table 3 and figure 1 are examples. Calibrator and controls must be used in each individual assay.

Results and calibration curve are provided for demonstration purposes only. A calibration curve must be generated for each set of samples to be assayed.

LIMITATIONS

1. Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.
2. The anti-MAG Autoantibodies ELISA has not been validated for plasmapheresis samples.

REFERENCE INTERVALS AND CUT-OFF

The frequency of anti-MAG auto-antibodies in normal human sera was determined using blood samples from asymptomatic volunteer blood donors (adult men and women at the age of 18 to 70 years). 150 samples were assayed according to the assay procedure and the results listed in table 4 were obtained.

Note: These titer ranges should be used as guidelines only. It is recommended that each laboratory establishes its own expected ranges for its patient population.

Proposed cut-off Titer

The mean +3 SD results in a technical cut-off value of 729 BTU. For practical reasons, we recommend to use a cut-off value of 1000 BTU, corresponding to the lowest calibrator (=calibrator D) in the standard curve.

PERFORMANCE CHARACTERISTICS

Intra-Assay Precision (Within-Run): 6.5 %. The intra-assay precision was calculated from the results of 20 pairs of values obtained in a single run (table 5).

Inter-Assay Precision (Run-to-Run): 15.4 %. The inter-assay precision was calculated from the results of 20 pairs of values obtained in 20 different runs (table 6).

Limit of Blank (LoB): 444 BTU. 20 duplicates of incubation buffer were assayed in a single run. Mean and standard deviation were calculated for the absorbance values. The minimal detectable dose of anti-MAG antibodies was calculated to be 444 BTU by adding two standard deviations to the mean absorbance and intersecting this value with the standard curve obtained in this run.

Limit of Quantification (LoQ): 900 BTU. 20 duplicates of low titer anti-MAG auto-antibodies were assayed in a single

run. Mean, standard deviation (SD) and coefficient of variation (CV) were calculated from the absorbance values. We found a titer of 900 BTU with a CV less than 10 %.

Dilution Linearity/Parallelism: 147 %. Seven human serum samples containing high titer of anti-MAG antibodies were diluted with incubation buffer 1:1000 to 1:64000, left for one hour at 18-28 °C and subsequently assayed according to the assay procedure (table 7). It is suggested that the relatively high deviation in about 50 % of the samples is due to antibody aggregations. In general, pathological sera show strongly elevated titer of auto-antibodies therefore it has no influence to the positive/negative discrimination.

Specificity: Four sets of experiments were performed to assess the specificity of the anti-MAG Autoantibodies ELISA:

1. Neutralization of anti-MAG auto-antibodies: Five sera with high anti-MAG titers could be inhibited from binding to the microtiter plates coated with MAG in a concentration-dependent manner when preincubated for 1 hour with incubation buffer supplemented with 1 to 200 µg/mL of MAG prior to testing in the ELISA.

2. Specificity of anti-MAG auto-antibody binding: Five sera with medium to high anti-MAG auto-antibody titers and five negative sera were tested on GanglioCombi (EK-GCO) plate coated with GA1, GM1, GM2, GD1a, GD1b and GQ1b. None of the sample showed ratio higher than 10 %. Similarly, six sera with high anti-Ganglioside auto-antibody titers (see above) and five negative sera, were tested on plates coated with MAG. None of these disease state sera gave a signal higher than 300 BTU.

3. Comparison to immunoblot (Western Blot): 127 patients sera (40 female; 87 male; age 4-90 years) with neurological diseases suspectedly caused by anti-MAG auto-antibodies were tested by the ELISA and a western blot method, respectively. MAG isolated from human central nervous system (CNS) was used in both methods. All twelve sera showing medium to very high titers of anti-MAG IgM auto-antibodies in the ELISA procedure were also clearly positive when analyzed on immunoblots. 114 out of 115 sera showing anti-MAG IgM antibody titers below the ELISA cut-off value were also negative when analyzed on immunoblots, whereas one serum was slightly positive by the immunoblot method only (F. Ferracin and A.J. Steck, unpublished results).

4. Comparison to indirect immunofluorescence assay (IFA): Sera from 150 patients (unknown sex and age) diagnosed as having IgM monoclonal gammopathies were tested by the ELISA and an IFA method, respectively. In the IFA method, sera were incubated on frozen and acetone-fixed monkey sciatic nerve sections and bound anti-MAG antibodies were detected by FITC-labeled anti-human IgM antibody. 26 patients (17.3 %) were positive and 115 patients (76.7 %) were negative in both methods, whereas five sera (3.3 %) were positive in the IFA method only and four sera (2.7 %) were only positive in the ELISA method, respectively (2).

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Example of Standard Curve (OD₄₅₀)


Figure 1

Cut-off values and Standardization

| | |
|-------------------|--------|
| Total (n) | 150 |
| Range (BTU) | 0 -832 |
| Mean Value (BTU) | 173 |
| SD (BTU) | 186 |
| Mean + 3 SD (BTU) | 729 |

Table 4

Intra-Assay Precision (Within-Run)

| Sample Type | Mean [BTU] | SD [BTU] | CV [%] |
|----------------|---------------|-------------|-----------|
| Serum 1 (Low) | 3764 | 176 | 4.7 |
| Serum 2 (High) | 35384 | 2920 | 8.3 |
| Mean | | | 6.5 |

Table 5

Inter-Assay Precision (Run-to-Run)

| Sample Type | Mean [BTU] | SD [BTU] | CV [%] |
|----------------|---------------|-------------|-----------|
| Serum 3 (Low) | 4208 | 714 | 17.0 |
| Serum 4 (High) | 17112 | 2362 | 13.8 |
| Mean | | | 15.4 |

Table 6

Dilution Linearity

| Serum Samples | Range Min-Max | Mean Observed/Expected |
|---------------|------------------|------------------------|
| 5 | 78-96 % | 86 % |
| 6 | 99-119 % | 107 % |
| 7 | 105-151 % | 123 % |
| 8 | 162-229 % | 206 % |
| 9 | 135-231 % | 196 % |
| 10 | 112-292 % | 233 % |
| 11 | 59-98 % | 76 % |
| Mean | | 147 % |


Table 7

APPENDIX II

REFERENCES / LITERATURREFERENZEN / RÉFÉRENCES / RIFERIMENTI / REFERENCIAS/ REFERÊNCIAS

1. Willison, H.J. et al.: Use of antibody testing in nervous system disorders. European Handbook of Neurological Management, Vol 1, 2nd edition, edited by Gilhus N.E. et al., (2011).
2. Bourque, R. P. et al.: Autoimmune peripheral neuropathies. Clin Chim Acta 449, 37-42 (2015).
3. Nobile-Orazio E Update on neuropathies associated with monoclonal gammopathy of undetermined significance . JPNS 15, 302-306 (2010) .
4. Kuijf M et al.; Detection of anti-MAG antibodies in polyneuropathy associated with IgM monoclonal gammopathy. Neurology 73(9), 688-95 (2009).
5. Hadden R.D.M et al., Paraproteinemic demyelinating neuropathies; European Handbook of Neurological Management: Volume 1, 2nd Edition, © 2011 Blackwell Publishing Ltd. ISBN: 978-1-405-18533-2
6. Pruppers M.H.J. et al., Improving future assessment and Research in IgM anti-MAG peripheral neuropathy: a consensus collaborative effort; Neuromuscular Disorders **27**, 1065-1072, (2017)
7. Dalakas M C.; Advances in the diagnosis immunopathogenesis and therapies of IgM- anti-MAG antibody-mediated neuropathies; Ther Adv Neurol Disord **11**, 1-12 (2018)

anti-MAG Autoantibodies ELISA

TIME TO RESULT: 4.5 HOURS

APPENDIX IV

NOTES / NOTIZEN / NOTES / NOTE / NOTAS

APPENDIX V

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLI / SIMBOLOS

| Symbol | Explanation | Symbol | Explanation |
|---|---|---------------------|--|
|  | Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad Data de expiração | MP | Microtiter Plate Mikrotiterplatte Microplaque Micropiastra Microplaca Placa de microtitulação |
| REF | Catalogue Number Bestellnummer Référéncie du catalogue Numero di catalogo Número de catálogo Número de catálogo | BUF WASH 10X | Wash Buffer Concentrate (10x) Waschpuffer Konzentrat (10x) Tampon de lavage concentré (10x) Tampone di lavaggio concentrato (10x) Tampón de lavado concentrado (10x) Concentrado do tampão de lavagem (10x) |
| LOT | Batch Code Chargenbezeichnung Code du lot Codice del lotto Codigo de lote Código lote | CALA - CALD | Calibrator A - D Kalibrator A - D Calibreur A - D Calibratore A - D Calibrador A - D Calibrador A - D |
| IVD | <i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i> Producto sánitarío para diagnóstico <i>in vitro</i> | CONTROL L | Low Control Kontrolle tief Contrôle bas Controllo basso Control bajo Controle baixo |
|  | Contains Sufficient for <n> Tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos Conteúdo suficiente para <n> tests | CONTROL H | High Control Kontrolle hoch Contrôle élevé Controllo alto Control alto Controle alto |
|  | Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso Leia cuidadosamente as instruções | EL IgM | Enzyme Label IgM Enzymmarker IgM Marqueur enzymatique IgM Marcatore enzimatico IgM Marcador enzimático de IgM Marcador enzimático IgM |
|  | Temperature Limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura Límite de temperatura | SUBS TMB | TMB Substrate TMB Substrat Substrat TMB Substrato di TMB Substrato TMB Substrato TMB |
| BUF INC | Incubation Buffer Inkubationspuffer Tampon d'incubation Tampone di incubazione Tampón de incubación Tampão de incubação | SOLN STOP | Stop Solution Stopp-Lösung Solution stop Soluzione stoppante Solución de parada Solução de parada |