Quantum Blue[®] Infliximab

Kvantitativni test
lateralnog protoka

za *in vitro* dijagnostičku upotrebu

LF-TLIF25 25 testova

LF-TLIF10 10 testova

Datum izdavanja: 2023-02-20
Verzija A6

 **Proizvođač**

BÜHLMANN Laboratories AG
Baselstrasse 55
4124 Schönenbuch, Švajcarska
Tel.: +41 61 487 12 12
Fax: +41 61 487 12 34
info@buhlmannlabs.ch

NAMENA

Quantum Blue® Infliximab je *in vitro* dijagnostički imunotest lateralnog protoka za kvantitativno određivanje donjeg nivoa infliksimaba u uzorcima seruma. Test služi kao pomoć terapijskog praćenja leka kod pacijenata sa inflamatornim bolestima creva (IBD), u okviru terapije infliksimabom, zajedno sa drugim kliničkim i laboratorijskim nalazima. Quantum Blue® Infliximab se koristi sa Quantum Blue® Reader uređajem. Za laboratorijsku upotrebu.

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Test je dizajniran za selektivno određivanje infliksimaba imunotestom na principu sendviča. Rekombinantni faktor nekroze tumora alfa (TNF α) je konjugovan sa koloidom zlata. Kada uzorak se aplikuje na test kasetu konjugat se oslobađa iz jastučića u reakcionom sistemu. Infliksimab prisutan u uzorku veže se za konjugat zlata. Monoklonsko antitelo, visoko specifično za analit, imobilise se na test membrani i veže kompleks konjugata zlata i infliksimab analita, rezultujući bojenjem test linije (T). Preostali slobodan konjugat TNF α /konjugat zlata veže se za kontrolnu liniju (C). Intenzitet signala test linije (T) kontrolne linije (C) se meri kvantitativno uz pomoć Quantum Blue® Reader uređaja.

OBEZBEĐENI REAGENSI I PRIPREMA

Reagensi	Količina		Kod	Komentari
	LF-TLIF25	LF-TLIF10		
Test kasete	25 komada	10 komada	B-LFTLIF-TC	Vakuumski zatvoren u foliji
Chase pufer	1 bočica 10 mL	1 bočica 10 mL	B-LFTLIF-CB	Spreman za upotrebu
Kontrole Niska*/Visoka*	2 bočice 0.5 mL	2 bočice 0.5 mL	B-LFTLIF-CONSET	Spreman za upotrebu
RFID Chip kartica	1 komad	1 komad	B-LFTLIF-RCC	Bela plastična kartica
RFID Chip kartica	1 komad	1 komad	B-LFTLIF-RCC15	Zelena plastična kartica
Bar-kod kartica	1 komad	1 komad	B-LFTLIF-BCC	2D bar-kod plastična kartica

Tabela 1

* Koncentracija infliksimaba koju sadrže kontrole je specifična za lot. Pogledati dodatne QC liste za aktuelne koncentracije.

PROVERITE TEST KOMPLET

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
Slika 3

REFERENCE

1. Blirup-Jensen et al.: *Protein standardization V: value transfer. A practical protocol for the assignment of serum protein values from a Reference Material to a Target Material*. Clin Chem Lab Med; 46, 1470 – 9 (2008)
2. Mitrev N et al.: *Review article: consensus statements on therapeutic drug monitoring of anti-tumor necrosis factor therapy in inflammatory bowel disease*. Aliment Pharmacol Ther. 46(11-12):1037-1053 (2017)
3. Vande Castele, N.: *Trough Concentration of Infliximab Guide Dosing for Patients With Inflammatory Bowel Disease*. Journal of Gastroenterology, 148, 1320–1329 (2015)
4. Papamichael K, et al.: *Improved Long-term Outcomes of Patients With Inflammatory Bowel Disease Receiving Proactive Compared With Reactive Monitoring of Serum Concentrations of Infliximab*. Clin Gastroenterol Hepatol. 15(10):1580-1588 (2017)
5. Deora V, et al.: *Therapeutic drug monitoring was helpful in guiding the decision-making process for children receiving infliximab for inflammatory bowel disease*. Acta Paediatr Int J Paediatr. 106(11):1863-1867. (2017)
6. Mitchell RA, et al.: *The Utility of Infliximab Therapeutic Drug Monitoring among Patients with Inflammatory Bowel Disease and Concerns for Loss of Response: A Retrospective Analysis of a Real-World Experience*. Can J Gastroenterol Hepatol. 2016
7. Amiot A, et al.: *Therapeutic drug monitoring is predictive of loss of response after de-escalation of infliximab therapy in patients with inflammatory bowel disease in clinical remission*. Clin Res Hepatol Gastroenterol. 40(1):90-98. (2016)
8. Burgess C. et al.: *Utility of regular infliximab levels in pediatric Crohn's disease*. J Pediatr Gastroenterol Nutr. 63:S224 (2016).
9. Vande Castele N. et al.: *Detection of infliximab levels and anti-infliximab antibodies: a comparison of three different assays*. Aliment Pharmacol Ther. 36, 765-771 (2012)
10. Afonso J. et al.; *Therapeutic drug monitoring of CT-P13: a comparison of four different immunoassays*. Therap Adv Gastroenterol. 10(9):661-671 (2017)
11. Magro F. et al.: *The performance of Remicade®-optimized quantification assays in the assessment of Flixabi® levels*. Therap Adv Gastroenterol. 11 (2018)

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