

INTRODUCTION

BÜHLMANN GanglioCombi™ ELISA is a multiparametric test to determine anti-ganglioside antibodies in a semi-quantitative way. It allows for a targeted investigation of immune-mediated neuropathies. There are no recognized reference materials nor reference measurement procedures for anti-ganglioside antibodies. To ensure consistency of results over time, we guarantee a transparent traceability chain (Fig. 5). This is achieved by using an internal reference material (IRM) to produce standardized calibrators. Due to the similar structure of gangliosides on our ELISA, anti-GM1 antibodies can serve as a surrogate standardization.

METHODS

Sialic acid residues on gangliosides used in BÜHLMANN GanglioCombi™ ELISA determine fine specificities of anti-ganglioside antibodies. GM1 represents the antigenic core structure of four C6 sugars and one sialic acid, that is commonly shared among different gangliosides (Fig. 1).

Therefore, anti-GM1 antibodies were used exemplarily to standardize multiparametric anti-ganglioside antibodies ELISAs (Fig. 2).

RESULTS


An IRM was generated from monoclonal anti-GM1 IgG and IgM antibodies (Fig. 4). Following the protocol by Blirup-Jensen et al., 2008, the value of the IRM is assigned to a calibrator stock, which is subsequently gravimetrically diluted into calibrators (Fig. 5). Based on the anti-GM1 antibodies standardization, the combined relative uncertainty of the calibrators of the BÜHLMANN GanglioCombi™ ELISA is calculated. The IRM traceable Ganglioside-ELISAs will be compared to the current version in a validation stage. The assays will be submitted to IVDR, the new European regulatory basis for in vitro diagnostic medical devices.

FIGURES

Figure 1

Ganglioside structures


Glycolipid structures of gangliosides involved in human neuropathies with GM1 as major backbone component



Adapted from Goodfellow & Willison, 2018, Progress in Molecular Biology and Translational Science ISSN 1877-1173; Jerkovic & Kalanj-Bognar, 2010, DOI: 10.2478/v10134-010-0043-6

Figure 2

Accessibility GM1 by ganglioside binding antibodies



adapted from Horwacik et al. Molecular & Cellular Proteomics, 2015

Figure 3

Principle of the assay: sandwich ELISA




Figure 4

Reference material components

Two human, monoclonal anti-GM1 antibodies





Figure 5

Metrological Traceability & Standardization




$$U_{comb} \text{ of Calibrator at 95\% CI} = 2x \sqrt{u_{value\ transfer}^2 + u_{CASET\ dilution}^2}$$

$$2x \sqrt{14.64^2 + 0.41^2} = 29.3\%$$


The combined uncertainty of the calibrator at the 95% CI is 29.3%.

Value Transfer Protocol


	# dilutions	# runs	# replicates	# days
IRM	6	3	2	4
Calibrator Material	5 - 6	3	2	4
IRM Mastercontrol	1	3	2	4



Value Assignment Calibrator IgG



Value Assignment Calibrator IgM



CONCLUSIONS

Based on the modularity and structural similarity of various gangliosides, that are targeted by autoantibodies in neuropathies, anti-GM1 antibodies can serve as a surrogate standardization. The transparent anti-GM1 antibodies-based traceability chain of the anti-ganglioside antibodies ELISAs truly display the uncertainty of the assays. This concept not only fulfills the requirements of IVDR regulations, but also guarantees long-term robustness of these state-of-the-art commercial assays for the detection of anti-ganglioside antibodies.

