# Sleep Check

## Additional laboratory instructions for Direct Saliva Melatonin ELISA/RIA

B-SLEEPCHECK16  
16 profiles

Revision date: 2009-02-16

## ENGLISH

### INTENDED USE

The BÜHLMANN Sleep Check is intended for the accurate collection and handling of human saliva samples for the subsequent quantitative determination of the dim-light melatonin onset (DLMO) in human saliva (1-6).

Quantification of Melatonin in saliva will be done with the Direct Saliva Melatonin ELISA/RIA from BÜHLMANN. All basic information concerning quantification step will be found in the corresponding instruction for use. DLMO can be calculated with the help of this manual.

### BACKGROUND INFORMATION

DLMO is an important parameter for the diagnosis and treatment of circadian rhythm disorders. Endogenous melatonin, a hormone produced by the pineal gland during the dark phase of the day-night cycle, synchronizes circadian rhythms, especially the sleep/wake rhythm. Its secretion is controlled by an endogenous oscillator that is entrained by light. Melatonin concentration reaches the maximum saliva levels around 03:00-04:00 a.m., varying with chronotype, whereas diurnal levels are low to undetectable.

When endogenous melatonin starts to rise late, sleep/wake rhythm is delayed. Vice versa when it starts to rise early, sleep/wake rhythm is advanced. Melatonin rhythm correlates strongly with circadian rhythms. Therefore **DLMO is the most useful parameter to establish circadian rhythm disorders (4).**

The determination of Melatonin concentration in Saliva is a non-invasive method. With the Sleep Check Saliva collection device it is possible to determine a partial 5 points Melatonin profile, around the client sleep time. The first collection point has to be collected three hours before the client sleep time, preceded with a hourly Saliva collection 4 times. With the 5 points partial profile the determination of DLMO is possible in about 80% of all clients (pers. communication with Dr. M. Smits; Ede, NL).

### REAGENTS SUPPLIED

Components	Quantity	Code	Comments
Safety boxes	16	B-SB	Each box contains: 5 Salivette® 5 Salivette® Labels
Big Envelope	16	--	For shipment to client
Blue plastic Envelope	16	--	For reshipment to Lab
Client Manual	16	--	Instructions for clients
Client Data Form	16	--	

Table 1

One Sleep Check kit last for 16 individual partial profiles with 5 points of measurement. Using the BÜHLMANN Direct Saliva Melatonin ELISA (EK-DSM) for melatonin quantification the 16 individual partial profiles can be done with one kit.

### MATERIALS REQUIRED BUT NOT PROVIDED

The Sleep Check is a concept for the calculation of the individual DLMO and contains the required sample collection devices.

For quantification of the samples BÜHLMANN strongly propose to use either the BÜHLMANN Direct Saliva Melatonin ELISA or RIA:

- Direct Saliva Melatonin ELISA (order code: EK-DSM)
- Direct Saliva Melatonin RIA (order code: RK-DSM)

### PROCEDURAL NOTES FOR SPECIMEN COLLECTION

Collect saliva using the BÜHLMANN Saliva Collection Device (Salivette® tube) delivered with the Sleep Check kit. One of these cotton swabs can absorb up to 3 ml of saliva. The procedure calls for <0.5 ml of saliva.

Saliva collection is simple and can be performed in clients' home at his convenience, but following precautions are recommended:

*Please read the instructions carefully before you start to collect saliva. The more carefully you follow the instructions, the more reliable the results will be.*

Client should set aside an evening for this test avoiding sporting activities and any intense efforts as far as possible.

**Light:** Bright light can suppress melatonin production. Therefore, it is important to avoid bright light during the test. Muted lighting from a reading lamp or from the television is preferable.

**Eating:** Nothing should be eaten during the collection time. The last meal must be taken at least 30 minutes before starting the collection. Bananas and chocolate should not be eaten during the entire day before the collection.

**Drinking:** Drinks containing artificial colorants, caffeine (coffee, black or green tea, iced tea, cola) or alcohol are not allowed on the evening of the collection.

**Medicines:** On the collection day, if possible, no aspirin and medicines that contain ibuprofen (Algiofor, Brufen, Dysmenol, Dolocyl, Ecoprofen) should be taken. If your sleep or sleep-wake rhythm is treated with melatonin, this must be discontinued at least one week before the collection.

## COLLECTION SCHEME


Five saliva samples must be collected during one evening. The time point of collection is depending on individual sleep time.

The first collection has to be done three hours before the regular individual sleep time. The last collection point has to be collected one hour after the regular sleep time.

**Example:** If the regular sleep time of the client is 11 pm, the first collection point has to be at 8 pm, afterwards every hour one saliva sample has to be collected (at 9 pm, 10pm, 11 pm and 12 pm).

Example for Collection Time					
Sleep time	1 <sup>st</sup>	2 <sup>nd</sup>	3 <sup>rd</sup>	4 <sup>th</sup>	5 <sup>th</sup>
8 pm	5pm	6pm	7pm	8pm	9pm
10 pm	7pm	8pm	9pm	10pm	11pm
12 pm	9pm	10pm	11pm	12pm	1am
1 am	10pm	11pm	12pm	1am	2am

Collection points

## COLLECTION MODE

1. Before you begin the collection, label the tubes. Start with label 1. Write on the label the client ID number assigned to you, your date of birth, the day you started the saliva collection and the first planned collection time. Then attach the label to the tube and do the same with tubes 2 – 5. Fill out the client information form.

2. 15 minutes before each saliva sample, you should rinse your mouth thoroughly with water

3. Open the top of the tube (the cotton ball is in the top) and remove the top from the tube.

4. Put the cotton ball into your mouth straight from the top without touching it with your fingers.

5. Put the cotton ball between your teeth and cheek and move it around with your tongue for 3 – 5 minutes, until the cotton ball is thoroughly soaked with saliva.

6. Put the cotton ball straight into the tube from your mouth, without touching it with your fingers.

7. Using the top, push the cotton ball into the tube and put on the top.

8. Now repeat the whole of sequence 2 - 7 every hour to obtain 5 samples.

9. After having finished the collection period, put all tubes into the white plastic box and put the white box, together with the client information form, into the blue franked envelope.

10. Store the samples in the refrigerator until you dispatch them so that you can send them to the lab within 2 days (put them in the next mailbox).

**11. Please make sure you send samples to the laboratory on Monday to Thursday only.**

## SPECIMEN SHIPMENT AND STORAGE

**Shipment:** The collected saliva samples (Salivette® tubes) must be shipped to the laboratory packed into the white safety box within two days. Used Salivette® tubes must be kept in the fridge at approximately 2-8°C.

The safety box, containing the 5 labelled and used Salivette® tubes together with the filled client form must be packed in the blue plastic envelope and sent by mail to the analyzing laboratory.

**Attention:** Samples should not be sent on Fridays, Saturdays or the day before holiday.

**Please keep in mind:** If the identification number is missing or the labels are incorrectly filled, the laboratory will be unable to carry out the melatonin examination.

**Storage:** The saliva samples absorbed to the cotton swab may be stored in the saliva collection device for up to 7 days at 2-8°C. If not assayed within one week after collection, samples should be frozen and may be stored for at least 6 months at ≤ -20°C.

Repeated freeze-thaw cycles should be avoided.

## SAMPLE RECOVERY FROM SALIVETTE®

Centrifuge the Salivette® tubes sent by the client for around 5 min at 3000 rpm (~1500x g). Discard the suspended insert with the swap and store the tube at 2-8°C or -20°C.

## SAMPLE PRETREATMENT AND ASSAY PROCEDURE (LAB)

See corresponding instruction for use for Direct Saliva Melatonin ELISA or RIA

## CALCULATION OF THE DLMO

The calculation of the ELISA/RIA results is explained profoundly in the corresponding instruction for us of the Direct Saliva Melatonin ELISA/RIA.


The results of the five time points shall be integrated into a time axis. On the X-axis insert the results of the five time points and on the Y-axis insert the calculated concentration of Melatonin in pg/ml. Sign the point on the graph and make a point to point fitting. BÜHLMANN recommend to use the DLMO definition as the time at which individual salivary melatonin reaches a value of 4 pg/ml (4)

### Example:

Regular client's retiring time: 12 Midnight (12 am).

Collection time	Melatonin [pg/ml]
21:00 (9 pm)	1
22:00 (10 pm)	0.8
23:00 (11 pm)	2.5
00:00 (12 am)	9.9
01:00 (1 am)	17.9

After the integration of the Melatonin values on the graph, the DLMO time can be read out, on the Melatonin concentration of 4 pg/ml. In this example the DLMO of this client is around 23:15 pm.

The determination of some clients' DLMO is not possible with a 5 points partial melatonin profile:

- Concentrations of Melatonin do not reach 4 pg/ml, these clients are low melatonin secretors or the DLMO time of this client is later than the last collection point.
- The Melatonin base line is higher than 4 pg/ml over the five collection points. Most probably the client's DLMO time is earlier than the first collection point.

BÜHLMANN recommends in these cases, to collect samples for a 24 hours profile to determine the DLMO.

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