Quantum Blue[®] fCAL extended

**Quantitative
Lateral Flow Assay**

LF-CALE25

25 tests

Revision date: 2016-11-09

ENGLISH

INTENDED USE

The BÜHLMANN Quantum Blue® fCAL extended is an *in vitro* diagnostic test for the quantitative determination of calprotectin in human stool specimens intended as an aid in the assessment of intestinal mucosal inflammation. The assay results can be used as an aid to diagnosis in distinguishing organic, inflammatory disease of the gastrointestinal tract (inflammatory bowel disease, IBD, e.g. Crohn's disease or ulcerative colitis, UC) from functional disease (irritable bowel syndrome, IBS) (ref. 1-6), in patients with chronic abdominal pain, above the age of four (ref. 7-8), and as an aid to IBD disease monitoring (ref. 9-10).

For professional use only*.

*Canada, Taiwan: For laboratory use only.

PRINCIPLE OF THE ASSAY

The test is designed for the selective measurement of calprotectin antigen by sandwich immunoassay. A monoclonal capture antibody (mAb) being highly specific for calprotectin is coated onto the test membrane. A second monoclonal detection antibody conjugated to gold colloids is deposited onto the conjugate release pad and released into the reaction system after addition of the extracted and diluted stool sample. The calprotectin / anti-calprotectin gold conjugate binds to the anti-calprotectin antibody coated on the test membrane (Test Line; Test Band) and the remaining free anti-calprotectin gold conjugate binds to the goat anti-mouse antibody coated on the test membrane (Control Line; Control Band). The signal intensities of the Test Line and the Control Line are measured quantitatively by the BÜHLMANN Quantum Blue® Reader.

REAGENTS SUPPLIED AND PREPARATION

Reagents	Quantity	Code	Comments
Test Cassette	25 pieces	B-LFCALUS-TC	vacuum-sealed in a foil bag
Extraction Buffer	1 bottle 125 mL	B-CAL-EX	Ready to use
Controls Low* / High*	2 vials 0.5 mL	B-CALE-CONSET	Ready to use
RFID Chip Card	1 piece	B-CALE-RCC	White plastic card
RFID Chip Card	1 piece	B-CALE-RCC720	Green plastic card

Table 1

* The controls contain lot specific amounts of native human calprotectin. Refer to the additional QC data sheet for actual concentrations.

REAGENTS & MATERIAL SUPPLIED SUPPLEMENTARY

Fecal Extraction Devices

Fecal Extraction Devices described below are not delivered with the kit and either of them has to be ordered with the kit.

CALEX® Cap Device	50, 200 or 500 tubes filled with 5 mL Extraction Buffer, ready to use	B-CALEX-C50 B-CALEX-C200 B-CALEX-C500
Smart-Prep	50 tubes consisting of spatulas and base caps	B-CAL-RD
ScheBo® Quick-Prep™	50 tubes consisting of tube, cone & dosing tip. Filled with 1.3 mL Extraction Buffer, ready to use	B-CAL-SO50

Table 2

STORAGE AND SHELF LIFE OF REAGENTS

All kit components are stable at 2-8 °C until the expiration date printed on the labels.

PRECAUTIONS

Safety Precautions

- None of the reagents of this assay contains components of human origin.
- Patient specimens should be handled as if capable of transmitting infections and should be handled in accordance with Good Laboratory Practice (GLP) using appropriate precautions.
- **Reagents:** Avoid contact of reagents with the skin, eyes, or mucous membranes. If contact does occur, immediately wash with generous amounts of water; otherwise, irritation can occur.
- Unused solution should be disposed of according to local state and federal regulations.

Technical Precautions

Kit Components

- All reagents and test samples must be equilibrated to room temperature (18-28 °C) before starting the assay. Mix well (vortex) the reagents before use.
- Components must not be used after the expiry date printed on the labels.
- Do not mix different lots of reagents.
- Test Cassettes cannot be re-used.

Assay Procedure

- Read carefully the instructions prior to carrying out the assay. Assay performance will be adversely affected, if reagents are incorrectly diluted, handled or stored under conditions other than those as detailed in this instruction for use.
- The Quantum Blue® Reader must be switched on and programmed for the Quantum Blue® fCAL extended: B-CALE-RCC or B-CALE-RCC720, before starting the assay (see Quantum Blue® Reader Manual).
- Use the RFID Chip Card in order to change lot-specific test parameters.
- Patient samples that are not properly handled may cause inaccurate results.
- In order to receive reliable and quantitative results it is important to homogenize the stool sample entirely in the Extraction Device.

- Diluted samples should be used within several hours and cannot be stored for a longer time period.
- With BÜHLMANN Smart Prep and ScheBo® Quick-Prep™, it is important to centrifuge the extracts before storage. Centrifuge the tubes for 5 minutes at 3000 x g. After centrifugation the supernatant must be transferred into a fresh storage tube. With CALEX® Cap Device you can store the extracts directly with the device. For re-use / re-measurement of the extracts see step 2 under the chapter assay procedure.

MATERIALS REQUIRED BUT NOT PROVIDED

- Vortex mixer for stool extraction
- Precision pipettes with disposable tips: 10-200 µL and 1 mL
- Centrifuge
- 5 mL polypropylene or polystyrene tubes for dilution of the extracts
- Timer (optional)
- Quantum Blue® Reader available from BÜHLMANN (order code: BI-POCTR-ABS)
- Soft tissues or blotting paper

SPECIMEN COLLECTION AND STORAGE

If the Extraction Devices are used (refer to Table 2), less than 1 g of native stool will be required.

Collect stool samples into plain tubes. The samples can be stored refrigerated at 2-8 °C for at least 6 days.

The extracts are stable for at least 7 days at 2-8 °C and for at least 24 months at ≤-20 °C.

Important: The sample must be collected without any chemical or biological additives in the collection device.

ASSAY PROCEDURE

The assay procedure consists of three steps:

1. Extraction of Stool Samples:

The extraction is described in the instruction for use delivered with the respective extraction devices.

2. Sample Processing:

- Smart Prep or ScheBo® Quick Prep™: Let the stool extract settle for 10 minutes after extraction. Dilute the supernatant 1:10 with Extraction Buffer (e.g. 50 µL stool extract and 450 µL Extraction Buffer) and mix well. Let the samples equilibrate for at least 5 minutes at 18-28 °C prior to proceeding to the next step (step no. 3).
- CALEX® Cap Device: After extraction, let the stool extract settle for 10 minutes with the white head of the device down. Unscrew the blue cap. The supernatant can be used without further dilution in the lateral flow assay.

3. Lateral Flow Assay Procedure and Readout:

Two alternative methods can be loaded from the respective RFID Chip card: B-CALE-RCC or B-CALE-RCC720. Select one of the RFID chip cards before starting the experiments. Load the test method from the RFID Chip Card.

3.1. Method <B-CALE-RCC720> with internal timer

- Use the green plastic card.
- Load the Test Cassette onto the Test Cassette holder of the Reader.
- Add 60 µL of diluted stool extract onto the sample loading port of the Test Cassette.
- Close the Test Cassette holder and start the measurement by pressing the start button.
- The scan starts automatically after 12 minutes (720 seconds).
- For Low / High Controls: Repeat step 3.1 with 60 µL controls instead of diluted stool extract.

3.2. Method <B-CALE-RCC> without internal timer

- Use the white plastic card.
- Add 60 µL of diluted stool extract onto the sample loading port of the Test Cassette.
- Incubate for 12 minutes +/- 1 minute (set a timer manually).
- Load the Test Cassette onto the Test Cassette holder of the Reader.
- Scan the Test Cassette with the Quantum Blue® Reader by pressing the start (<ENTER>) button immediately.
- For Low / High Controls: Repeat step 3.2 with 60 µL controls instead of diluted stool extract.

Remark: Please refer to the Quantum Blue® Reader Manual to learn more about the basic functions and how to initialize and operate the Reader, especially how to select test methods, and how to load lot specific parameters from the RFID Chip Card in order to get the samples measured.

QUALITY CONTROL

- If the precision of the assay does not correlate with the established limits and repetition excludes errors in technique, check the following issues: i) pipetting, temperature controlling and timing devices ii) expiration dates of reagents and iii) storage and incubation conditions.
- Result of the self-test of the Quantum Blue® Reader performed at startup of the instrument has to be valid.

VALIDATION OF RESULTS

- For a valid test result, the Control Line (C) must be visible in any case (see Figures 1A and 1B). It is used as functional test control only and cannot be used for the interpretation of the Test Line (T). If the Test Line (T) is not detectable after 12 minutes of incubation time (Figure 1A), the concentration of calprotectin present in the stool sample is below the detection limit. If a Test Line (T) is detectable after 12 minutes of incubation time (Figure 1B), the calprotectin concentration present in the stool sample is calculated by the Quantum Blue® Reader.
- If only the Test Line (T) is detectable after 12 minutes of incubation time (Figure 1C), the test result is invalid and the assay has to be repeated using another Test Cassette.

- If neither the Control Line (C) nor the Test Line (T) are detectable after 12 minutes of incubation time (Figure 1D), the test result is invalid and the assay has to be repeated using another Test Cassette.
- As the Quantum Blue® Reader allows a quantitative evaluation of the Test (T) and Control (C) Lines, an additional validity check of the Control Line (C) is undertaken. If the signal intensity of the Control Line (C) is below a threshold after 12 minutes of incubation time, the test result is also invalid and the Quantum Blue® fCAL high range assay has to be repeated using another Test Cassette.

STANDARDIZATION AND INTERPRETATION OF RESULTS

- The Quantum Blue® fCAL is calibrated with the BÜHLMANN fCAL® ELISA (order code: EK-CAL).
- The BÜHLMANN Quantum Blue® Reader uses a lot-specific standard curve to calculate the calprotectin concentration. The assay range is between 30 and 1000 µg/g.
- For quantitative measurements, unknown samples reading above 1000 µg/g can be re-tested in the BÜHLMANN Quantum Blue® fCAL high range (order code: LF-CHR25).

LIMITATIONS

- Reagents delivered with the BÜHLMANN Quantum Blue® fCAL extended kit are intended for the determination of calprotectin levels in human stool samples only.
- Fecal calprotectin values can be used as supplementary data available to the physician in establishing a diagnosis.
- A calprotectin test result obtained with the BÜHLMANN Quantum Blue® fCAL extended can be confirmed by measuring other biomarkers (such as C-Reactive Protein) or by repeating the calprotectin test (with BÜHLMANN fCAL® ELISA or using a second Test Cassette) and comparing to the CDAI score and the clinical status of the patient, before any clinical and/or therapeutic decision is made.
- Patients who are taking NSAIDs regularly may have elevations in their fecal calprotectin levels.

EXPECTED VALUES

Estimation of fecal calprotectin is a reliable and simple way to distinguish organic from functional gastrointestinal diseases.

The following data were established with the BÜHLMANN fCAL® ELISA (order code: EK-CAL). In a clinical study, 344 symptomatic patients were compared with endoscopic findings. Endoscopy examination showed 264 patients with functional diseases, whereas 80 patients had various organic diseases (colitis, Crohn's, ulceration, diverticulitis, polyps, adenomas, cancer, or infectious diseases) (ref. 9). ROC curve analysis (AUC: 0.962) resulted in an optimal clinical cut-off at 50 µg/g. Applying this cut-off, a clinical sensitivity and specificity of 88.8 % and 93.6 %,

respectively can be reached in the differentiation between organic and functional diseases (Table 13). Fecal calprotectin levels from adults and children are comparable, whereas levels of newborns can be significantly increased (ref.6).

These data support the following recommendation for interpretation of results:

Values below 50 µg/g:

Calprotectin values <50 µg/g are not indicative of inflammation in the gastrointestinal tract. Patient- samples with low levels are likely not to be in need of invasive procedures to determine the inflammation cause (ref. 4).

Elevated Values between 50 and 200 µg/g:

Calprotectin values between 50 and 200 µg/g can represent mild organic disease such as inflammation caused by NSAIDs, mild diverticulitis and IBD in remission phase. The low inflammatory response shown within this range may suggest repeating the measurement and performing further investigations.

Elevated Values above 200 µg/g:

Calprotectin values >200 µg/g are indicative of active organic disease with inflammation in the gastrointestinal tract. Appropriate further investigative and curative procedures by specialists are suggested.

The cut-off level suggested for adults (50 µg/g) can be used for children aged from 4 to 17 years regardless of sex (ref.7).

PERFORMANCE CHARACTERISTICS

Limit of Blank (LoB): <20 µg/g calprotectin.

Limit of Detection (LoD): <30 µg/g calprotectin.

Limit of Quantification (LoQ):

Lower LoQ: ≤30 µg/g calprotectin.

Upper LoQ: ≥1000 µg/g calprotectin.

The LoQ has been established with 14 stool samples. The resulting precision profile is shown in Figure 2. The limit of quantification corresponds to the concentration of calprotectin with an imprecision below 25 % CV allowing a quantitative measurement within the range from 30 (lower LoQ) to 1000 µg/g (upper LoQ).

Linearity: The results showed linearity within the indicated measuring range of 30 to 1000 µg/g of the Quantum Blue® fCAL extended assay for all samples.

Recovery: The recovery varied between 90.6 %-101.8 %.

Repeatability: 20.1 % mean CV. The repeatability was calculated from 8 stool extract samples. The mean values of three different lots of Test Cassettes are presented in (Table 14). The repeatability varied between 17.9 %-23.9 % CV.

Inter-lot Precision: 20.8 % mean CV. The inter-lot precision of the Quantum Blue® fCAL extended assay was calculated from 8 stool extract samples containing between 30 and 1000 µg/g calprotectin. The inter-lot data of the samples varied from 16.1 %-24.1 % CV.

Method Comparison: $R^2 = 0.94$; $y = 1.05-3.91 \mu\text{g/g}$

50 stool samples exhibiting calprotectin concentrations within the indicated measuring range of the Quantum Blue® fCAL extended assay were analyzed according to the assay procedure and compared with the values obtained by the BÜHLMANN fCAL® ELISA (order code: EK-CAL). The correlation data are illustrated in Figure 3.

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
Figure 3

REFERENCES / REFERENZEN / RÉFÉRENCES / RIFERIMENTI / REFERENCIAS / REFERÊNCIAS

1. Tibble JA et al.: *Use of surrogate markers of inflammation and Rome criteria to distinguish organic from nonorganic intestinal disease.* Gastroenterol **123**, 450-460 (2002)
2. Fagerhol MK: *Calprotectin, a faecal marker of organic gastrointestinal abnormality.* Lancet **356**, 1783-4 (2000)
3. Tibble JA et al.: *A simple method for assessing intestinal inflammation in Crohn's disease.* Gut **47**,506-513 (2000).
4. Tschanguizi et al.: *Stool Calprotectin as a Marker of inflammation.* Poster at UEGW, 2007
5. Hessian PA and Fisher L: *The heterodimeric complex of MRP-8 (S100A8) and MRP-14 (S100A9). Antibody recognition, epitope definition and the implications for structure.* Eur J Biochem **268**, 353-63 (2001).
6. Konikoff MR and Denson LA: *Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease.* Inflamm Bowel Dis **12(6)**, 524-34 (2006)
7. Fagerberg UL et al.: *Colorectal inflammation is well predicted by fecal calprotectin in children with gastrointestinal symptoms.* J Pediatr Gastroenterol Nutr **40**, 450-5 (2005)
8. Jahnsen J, Røseth AG, Aadland E. *Measurement of calprotectin in faeces.* Tidsskr Nor Lægeforen **128**, 743–5 (2008)
9. Manz M et al. *Value of fecal calprotectin in the evaluation of patients with abdominal discomfort: an observational study.* BMC Gastroenterology **12**, 5 (2012).
10. Naismith GD et al. *A prospective evaluation of the predictive value of faecal calprotectin in quiescent Crohn's disease.* J Crohns Colitis. **8**, 1022-9, 2014

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLI / SIMBOLOS

Symbol	Explanation	Symbol	Explanation
	Use By Verwendbar bis Utiliser jusqu'au Utilizzare entro Fecha de caducidad Validade até		Temperature limitation Zulässiger Temperaturbereich Limites de température Limiti di temperatura Limite de temperatura Limite de temperatura
REF	Catalogue number Bestellnummer Référence du catalogue Numero di catalogo Número de catálogo Número de catálogo	TC	Test Cassette Test Kassette Cassette de testi Card di rilevazione Cartucho de prueba Cassete de teste
LOT	Batch code Chargenbezeichnung No. du lot No. del lotto No. de lote No. do lote	BUF EX	Extraction Buffer Extraktionspuffer Tampon d'extraction Tampone di estrazione Tampón de extracción Tampão de extração
IVD	<i>In Vitro</i> Diagnostic Medical Device <i>In Vitro</i> Diagnostikum Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> Producto sanitario para diagnóstico <i>in vitro</i> Produto para Diagnóstico <i>in vitro</i>	CONTROL L	Control Low Kontrolle Tief Contrôle Bas Controllo Basso Control Bajo Controle Baixo
	Content sufficient for <n> tests Ausreichend für "n" Ansätze Contenu suffisant pour „n“ tests Contenuto sufficiente per „n“ saggi Contenido suficiente para <n> ensayos Conteúdo suficiente para <n> testes	CONTROL H	Control High Kontrolle Hoch Contrôle Élevé Controllo Alto Control Alto Controle Alto
	Consult Instructions for Use- Gebrauchsanweisung beachten Consulter le mode d'emploi Consultare le istruzioni per l'uso Consulte las instrucciones de uso Consultar as intruções de uso	RCC	RFID Chip card RFID Chipkarte RFID Carte à puce RFID Carta chip RFID Tarjeta chip RFID Cartão Chip