# Flow CAST<sup>®</sup>

## Basophil Activation Test (BAT) Flow Cytometry

FK-CCR 100 tests

Revision date: 2012-08-17

## ENGLISH

### INTENDED USE

The Flow CAST® kit is a basophil activation test (BAT) which can be used for the *in vitro* detection of immediate type allergic reactions and hypersensitivities against suspected allergens.

The test is intended for the *in vitro* diagnostic determination of expression of CD63 surface marker on basophils in whole blood by flow cytometry upon antigen stimulation.

### PRINCIPLE OF THE ASSAY

The assay is based on the method first described by Sainte-Laudy *et al.* 1994 and 1996 (1,2) where basophil activation by allergens or controls is detected by flow cytometry measured by the increase of the CD63 (gp53) at the cellular surface. IgE and non-IgE mediated reactions can be detected (3-5).

Stimulation Buffer and Allergen is added to EDTA whole blood from suspected allergic/hypersensitive patients. The allergen mimics the *in vivo* reaction where specific IgE bound to the cellular surface are bridged by the culprit allergen and activates an intracellular signaling cascade leading to the activation of the basophil. As a consequence, intracellular compounds bearing the transmembrane protein CD63 are fused to the cellular membrane and therefore exposed to the extracellular matrix.

As positive control, highly specific monoclonal antibody binding to the high affinity IgE binding receptor (FcεRI) or the unspecific cell activator fMLP is used.

Together with the cellular stimulation, Staining Reagent is added containing a mixture of monoclonal antibodies to human CD63 labeled with fluorescein isothiocyanate (anti-CD63-FITC) and to human chemokine receptor CCR3 labeled with phycoerythrin (anti-CCR3-PE). CCR3 is constitutively expressed on eosinophils and basophils (6,7) Erythrocytes are removed by a lysing reaction and after a short centrifugation step the cells are resuspended in Wash Buffer and analyzed by flow cytometry (*cf.* Flow cytometric Data Acquisition on page 4).

### REAGENTS SUPPLIED AND PREPARATION

| Reagents  | Quantity      | Code        | Reconstitution                                 |
|---|---------------|-------------|--|
| <b>Stimulation Buffer</b><br>containing calcium, heparin and IL-3     | 1 vial lyoph. | B-CCR-STB   | Reconstitute with 50 ml of water <sup>1)</sup> |
| <b>Stimulation Control</b><br>anti-FcεRI mAb                          | 1 vial lyoph. | B-CCR-STCON | Reconstitute with 1.5 ml of B-CCR-STB          |
| <b>Stimulation Control</b><br>fMLP <sup>2)</sup>                      | 1 vial lyoph. | B-CCR-FMLP  | Reconstitute with 1.5ml of B-CCR-STB           |
| <b>Staining Reagent</b><br>Mix of anti-CD63-FITC and anti-CCR3-PE mAb | 1 vial 2.2 ml | B-CCR-SR    | Ready to use                                   |
| <b>Lysing Reagent</b> <sup>3)</sup><br>10x concentrated               | 1 vial 25 ml  | B-CCR-LYR   | Dilute with 225 ml of deionized water          |
| <b>Wash Buffer</b>  | 1 vial 100 ml | B-CCR-WB    | Ready to use                                   |

Table 1

<sup>1)</sup> For required water quality, see Chapter Procedural Notes

<sup>2)</sup> N-formyl-methionyl-leucyl-phenylalanine

<sup>3)</sup> Crystals may be formed during storage at 2-8°C and should be dissolved at 18-28°C prior to dilution.

### STORAGE AND SHELF LIFE OF REAGENTS

| Unopened reagents                                    |  |
|--|--|
| Store at 2-8°C. Do not use past kit expiration date. |  |
| Opened / reconstituted reagents                      |  |
| Stimulation Buffer                                   | Stable at -20°C for 6 months. Aliquot if repeated use is expected. |
| Stimulation Control                                  |  |
| Stimulation Control fMLP                             | Stable at -20°C for 6 months. Aliquot if repeated use is expected. |
| Lysing Reagent                                       | Stable at 2-8°C for 6 months.                                      |
| Staining Reagent                                     | Stable at 2-8°C until expiration date.                             |
| Wash Buffer  |  |

Table 2

### ALLERGENS AND REAGENTS TO BE ORDERED ADDITIONALLY

Validated Allergens for the analysis in CAST® assays are offered by BÜHLMANN. Refer to the Allergen list on the webpage to obtain the respective order codes ([www.buhmannlabs.ch](http://www.buhmannlabs.ch)).

- Protein Allergens are shipped as concentrated liquids (1µl/vial) and must be diluted before use.
- Drug and Chemical Allergens are shipped lyophilized and must be reconstituted before use.

Refer to the BÜHLMANN Allergen Booklet and **Allergen Data Sheets** available on the website [www.buhmannlabs.ch](http://www.buhmannlabs.ch).

### ALLERGEN REAGENTS FROM OTHER SOURCES

Allergens from other sources might be used in the Flow CAST® assay with the following limitations:

- No matrix-bound allergens (solid or liquid phase).
- No allergen preparations containing cytotoxic compounds (stabilizers, preservatives) such as glycerol, phenol, sodium azide or merthiolate (thimerosal).

For the procedure to establish customer specific allergens for the CAST®-Assays ask your local distributor or contact BÜHLMANN.

### MATERIALS REQUIRED BUT NOT PROVIDED

- K-EDTA venipuncture tubes.
- Centrifuge for centrifugation at 500 x g.
- Disposable, pyrogen-free polypropylene or polystyrene test tubes and appropriate test tube racks for the stimulation  
NOTE: Polystyrene tubes should fit with the Flow Cytometer used (e.g. 12 x 75 mm FALCON tubes from Becton Dickinson; order code: 352052).
- Vortex Mixer.
- Precision pipettes with disposable, pyrogen-free tips: 10-100 µl, 100-1000 µl, 1-5 ml adjustable pipette and a 10-50 µl adjustable dispenser.
- Cylinder for preparing the Stimulation Buffer.
- Sterile, ultrapure and apyrogenic water for preparing the cell stimulation reagents (*cf.* Chapter Procedural Notes).
- Water bath set at 37°C.
- Distilled or deionized water as well as beaker or cylinder for the preparation of Lysing Reagent.
- Bottle-top dispensers for Lysing Reagent and Wash Buffer, respectively.
- Flow Cytometer with 488 nm (blue) excitation wavelength including appropriate software (*cf.* chapter Flow cytometric Data Acquisition).

## PRECAUTIONS

### SAFETY PRECAUTIONS

- Stimulation Buffer (B-CCR-CSB) contains components of human origin. As material from human may contain infectious agents e.g. Hepatitis B virus and HIV. Therefore *patient samples and kit components should be handled as if capable of transmitting infections.*

### TECHNICAL PRECAUTIONS

- Recommended Water Quality for the Flow CAST®. The use of sterile, ultrapure and apyrogenic water for reconstituting Stimulation Buffer (B-CCR-STB) is essential for good and reproducible basophil stimulation. The following sources of water may be used: Cell culture grade water, infusion grade water or deionized, double distilled water that is ultra filtrated in a periodically sanitized 10 kDa ultra filter.
- The Lysing Reagent (B-CCR-LYR) can be reconstituted with deionized, double distilled water or the same water quality that is used for the cell stimulation reagents.
- Precautions to avoid allergen contamination during cell stimulation: Aeroallergens in the laboratory may contaminate open blood samples and cell suspensions from patients potentially causing an elevated background release. Therefore, care must be taken to cover blood samples and cell stimulation tubes. Avoid dust mites, pollinating plants, latex gloves or equipment potentially containing latex and open windows in the laboratory where the cell stimulation is performed. Therefore, we recommend carrying out the cell preparation and stimulation steps in a laminar flow hood.
- For cell stimulation and labeling reaction, the use of tissue culture grade microtiter plates is possible.
- Components must not be used after the expiry date printed on the labels. Do not mix different lots of reagents.
- Avoid contamination of reagents.
- It is important reading through the Instructions for Use prior to commencing the test. Reliable results will be obtained only if this Instruction for Use is followed accurately.
- Samples that are not properly handled may cause inaccurate results.
- Test results should be interpreted in conjunction with information available from clinical assessment of the patient and other diagnostic procedures.

### SPECIMEN COLLECTION AND STORAGE

It is recommended that patients should avoid systemically administered antiallergenic drugs such as corticosteroids, chromoglycic acid (DSCG) for at least 24 hours prior to blood sampling.

Collect sufficient blood into **K-EDTA venipuncture tubes**. Fill the venipuncture tubes up to the dedicated volume with blood. Not sufficiently filled tubes (filling grade < 50%) lead to higher EDTA concentration in the sample and thus may give false negative results.

1 ml of whole blood is sufficient for about 18 test tubes.

Perform the cell stimulation immediately or store the blood sample refrigerated (2-8°C) for up to 48 hours. For detecting responses to drugs store the blood sample only up to 24 hours. **Do not centrifuge or freeze blood samples.**

## ASSAY PROCEDURE

**Important:** The following procedure is optimized for whole blood specimen collected with EDTA as anticoagulant.

1. Mix the anti-coagulated blood sample by inverting the venipuncture tube several times.
2. Prepare fresh and pyrogen-free 3.5 ml polypropylene or polystyrene tubes suitable for Flow Cytometry measurements.
3. For each patient, label the tubes e.g.:  
PB = patient background  
PC1 = stimulation control with anti-FcεRI Ab  
PC2 = stimulation control with fMLP  
A1-1 for allergen 1 with dilution 1  
A1-2 for allergen 1 with dilution 2  
etc.

### Stimulation and Staining

4. Pipet 50 µl of the corresponding stimulus to each tube for each patient  
PB tube: 50 µl of **Stimulation Buffer (background)**  
PC1 tube: 50 µl of **Stimulation Control**  
PC2 tube: 50 µl of **Stimulation Control fMLP**  
Ax-y tube: 50 µl of **Allergen**
5. Add 100 µl of Stimulation Buffer to each tube.
6. Add 50 µl of patient's whole blood to each tube. Be sure that the side wall and top of the tube are free of blood.
7. Mix gently.
8. Add 20 µl Staining Reagent to each tube.
9. Mix gently, cover the tubes and incubate for 15 minutes at 37°C in a **water bath**.  
(using an incubator will take about 10 minutes longer incubation time due to less efficient heat transfer).

### Lysing

10. Add 2 ml pre-warmed (18-28°C) Lysing Reagent to each tube, mix gently.
  11. Incubate for 5 -10 minutes at 18-28°C.
  12. Centrifuge the tubes for 5 minutes at 500 x g.
  13. Decant the supernatant by using blotting paper.
  14. Resuspend the cell pellet with 300 µl of Wash Buffer.
- Note:** Depending on Flow cytometry instrumentation more wash buffer (e.g. 800 µl) might be necessary.
15. Vortex gently.
  16. Acquire the data on the flow cytometer within the same day. If the samples are stored for several hours it should be kept protected from light at 2-8°C.

**Note:** Samples stored, protected from light at 2-8°C over night are still analysable. A slight decrease of fluorescence intensity and a lower basophil recovery can be observed.

## FLOW CYTOMETRIC DATA ACQUISITION

Flow cytometric acquisition can be performed on any flow cytometer working with a 488 nm argon laser diode (blue-green excitation light).

The flow cytometer must be equipped to detect Forward Scatter (FSC), Side Scatter (SSC) and the two fluorochromes FITC and PE.

Ensure that the flow cytometer is properly aligned and colour compensation is set.

During acquisition of the samples, make sure that on a FSC/SSC histogram the leukocyte population is separated into three discrete populations. Adjust the amplification (gain) of FSC and SSC signals to obtain a distribution that is shown in Figure 1. Refer to the flow cytometer product manuals for instructions.


Figure 1: Three discrete populations (lymphocytes, monocytes and granulocytes) in FSC/SSC histogram.

Typically, after 500-600 basophilic cells (gated as shown in figure 2) the acquisition can be stopped. At least 200 basophilic cells must be analyzed, requiring a total amount of 50'000-100'000 leukocytes to be acquired per sample. Because of the lower activation percentage in drug allergies each laboratory has to define its own confidence limits (i.e. in drug allergies the limit of basophilic cells analyzed should be set to 300 or more).

## DATA ANALYSIS

The analysis of the acquired data can be performed with any flow cytometry analysis software e.g. FlowJo, FloMax, CellQuest or others.

The analysis is based on two steps:

1. Set a gate 1 (R1) by including the entire basophil population  $CCR3^{pos}$  with low Side Scatter  $SSC^{low}$  (see Figure 2). eosinophils located on the high right side will be excluded due to their  $SSC^{high}$  position.
2. Calculate the percentage of CD63 positive cells (brightly fluorescent FITC; Q2) compared to the total amount of basophilic cells gated in R1 (see Figure 3 and Figure 4).




Figure 2: Selection of basophilic cells  $CCR3^{pos} / SSC^{low}$

| Gated Region              | Count (n=) | %     |
|---------------------------|------------|-------|
| Total                     | 78251      | 100.0 |
| R1                        | 655        | 0.8   |
| Q2 (CD63 <sup>pos</sup> ) | 11         | 1.7   |
| Q4 (CD63 <sup>neg</sup> ) | 644        | 98.3  |

Figure 3: Patient Background (PB) with STB only

| Gated Region              | Count (n=) | %     |
|---------------------------|------------|-------|
| Total                     | 72916      | 100.0 |
| R1                        | 606        | 0.8   |
| Q2 (CD63 <sup>pos</sup> ) | 467        | 77.1  |
| Q4 (CD63 <sup>neg</sup> ) | 139        | 22.9  |

Figure 4: Stimulation Control (STCON)

## LIMITATIONS

- Flow cytometry may produce false results if the cytometer has not been aligned perfectly, if the fluorescence emission has not been correctly compensated and if the regions have not been carefully positioned.
- Verify the preparations by eye to assess the efficacy of lysis. The erythrocytes may be incompletely lysed and appear on a light diffraction histogram in the same location as the leucocytes.

## QUALITY CONTROL

For the appropriate evaluation of results, different values should be taken into account:

- Typically, three distinct **leukocyte populations** lymphocytes, monocytes and granulocytes appear in the FSC/SSC plot. Their occurrence can be regarded as a criterion for the quality of the blood sample (time of collection, storage).
- Negative control** (buffer control): A patient individual base value of 2.0 to 2.5 % of the activated basophils should be regarded as negative. This was also basis for the determination of allergen specific cut-offs.
- Positive control** (stimulation control). Two different positive controls are included in the kit. **Anti-Fc $\epsilon$ RI mAb** mimics the bridging of the receptor *in vivo* caused by the allergen. **fMLP** is a tripeptide causing basophil activation in a non-immunologic way. If one of those two stimulators shows activation of **>10 %** basophils, the sample is evaluable.
- Non-responders** are persons with a low reactivity (<10% CD63 positive cells) to fMLP and anti-Fc $\epsilon$ RI antibody. Evaluation of data from normal blood donors and patients showed that 6.1% (out of n=98) were non-responders to anti-Fc $\epsilon$ RI and 4.9 % (out of n=61) to fMLP.

## INTERPRETATION OF RESULTS

To obtain an optimal sensitivity and specificity, slightly different cut-off values should be applied for different groups of allergens. Based on numerous studies and evaluations done with different basophil activation tests Bühlmann recommends using the following cut-off:

|                     |             |             |
|---------------------|-------------|-------------|
| Inhalant allergens: | $\geq 15$ % |             |
| Food allergens:     | $\geq 15$ % |             |
| Hymenoptera venoms: | $\geq 10$ % |             |
| Betalactams*:       | $\geq 5$ %  | SI $\geq 2$ |
| Analgesics*:        | $\geq 5$ %  | SI $\geq 2$ |
| Food additives*     | $\geq 5$ %  | SI $\geq 2$ |

\* Drugs and other chemical allergens usually results in lower activation percentages than other allergens. Therefore, a lower cut-off should be taken, but the stimulation index (SI = allergen stimulation divided by negative control) must be equal or higher than 2 in order to consider the result as positive.

## ASSAY PERFORMANCE

**Specificity:** The anti-CCR3 mAb is a highly specific antibody described in 6 and 7. CCR3 is constitutively expressed on eosinophilic and basophilic leukocytes (see Fig. 2) and in a smaller part on CD3+ cells (lymphocytes). Samples from eight normal blood donors were double stained twice with anti-CCR3-PE and anti-CD3-AF647. The relative amount (mean) of CD3+ cells within the gated Basophil population was 3.9% (cf .Table 11)

Depending on the gating strategy used the relative amount of CD63+ basophils in Patient background (negative control) is very low. We analyzed 98 samples from normal blood donors and allergic patients. The median patient background was 1.5 % (95 % CI: 1.20 - 1.77 %). 6/98 samples showed a patient background >5 %.

**Basophil Recovery:** >500 basophils/stimulation tube. We analyzed 102 samples from normal blood donors and allergic patients. The median Basophil recovery was 526 cells (95% CI 481-578 basophils) stained with CCR3-PE.

**Precision (Patient background):** 16.2 % CV. The precision shows the reproducibility of the patient background within the same blood sample incubated 20 times with Stimulation Buffer and consecutively analyzed by flow cytometry. The results are expressed in Table 13 as percentage basophils activated and as mean fluorescence intensity (MFI) of CD63-FITC.

**Precision (Positive control):** 5.4 % CV. The precision shows the reproducibility of the stimulation within the same blood sample incubated 20 times with positive control (STCON) and consecutively analyzed by flow cytometry. The results are expressed in Table 13 as percentage basophils activated and as mean fluorescence intensity (MFI) of CD63-FITC.

**Inter-Technician Variation (Positive control):** 3.7-8.1 % CV. Two blood samples from normal blood donors were tested with the Flow CAST® by five different technicians within the same day. The two positive controls included into the kit, STCON and FMLP were used in duplicates.

The results are expressed in Table 14 as mean percentage basophils activated (CD63+) from double activation per technician.

**Reference Intervals** were established with normal blood donors (age 18-60). Median value for anti-Fc $\epsilon$ RI stimulated samples is 75.8 % (percentile 25th: 53.3 %, 75th: 85.6 %). Median value for fMLP stimulated samples is 48.1 % (percentile 25th: 35.4 %, 75th: 65.5%). Refer to Table 12.

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Table 13 Precision

|            | Patient Background |                       | STCON              |                       |
|------------|--------------------|-----------------------|--------------------|-----------------------|
|            | %CD63 <sup>+</sup> | MFI CD63 <sup>+</sup> | %CD63 <sup>+</sup> | MFI CD63 <sup>+</sup> |
| Mean       | 2.4%               | 10.1                  | 35.5%              | 82.1                  |
| SD         | 0.4%               | 0.8                   | 1.9%               | 7.2                   |
| <b>%CV</b> | <b>16.2%</b>       | <b>7.7%</b>           | <b>5.4%</b>        | <b>8.8%</b>           |

Table 12 Reference Intervals

|       | n  | %CD63 <sup>pos</sup> |   |  |
|-------|----|----------------------|---|--|
|       |    | Median               | 2.5 <sup>th</sup> – 97.5 <sup>th</sup> percentile | 25 <sup>th</sup> – 75 <sup>th</sup> percentile |
| STCON | 98 | 75.8                 | 3.4 – 92.0  | 53.3 – 85.6                                    |
| fMLP  | 61 | 48.1                 | 7.6 – 88.9  | 35.4 – 65.5                                    |

Table 14 Inter Technician Variation


|            | S1 (%CD63 <sup>+</sup> ) |             | S2 (%CD63 <sup>+</sup> ) |             |
|------------|--------------------------|-------------|--------------------------|-------------|
|            | STCON                    | FMLP        | STCON                    | FMLP        |
| Mean       | 64.8                     | 42.2        | 69.6                     | 48.1        |
| SD         | 3.6                      | 1.6         | 2.6                      | 3.9         |
| <b>%CV</b> | <b>5.6%</b>              | <b>3.7%</b> | <b>3.7%</b>              | <b>8.1%</b> |

**APPENDIX II**  
**REFERENCES/ LITERATURREFERENZEN/ RÉFÉRENCES/ RIFERIMENTI/ REFERENCIAS**

- Sainte-Laudy, J, et al. [Analysis of membrane expression of the CD63 human basophil activation marker. Applications to allergologic diagnosis]. *Allerg Immunol (Paris)* **26**, 211-4. (1994).
- Sabbah, A and Sainte-Laudy, J. Flow Cytometry applied to the analysis of Lymphocyte and Basophil activation. *ACI International* **8**, 116-9 (1996).
- Sanz, ML, et al. Flow cytometric basophil activation test by detection of CD63 expression in patients with immediate-type reactions to betalactam antibiotics. *Clin Exp Allergy* **32**, 277-86. (2002).
- DeWeck, AL and Sanz, ML. Flow cytometric cellular allergen stimulation Test (FAST/Flow-CAST): technical and clinical evaluation of a new diagnostic test in allergy and pseudo-allergy. *ACI International* **14**, 204-215 (2002).
- Gamboa, P et al. The flow-cytometric determination of basophil activation induced aspirin and other non-steroidal anti-inflammatory drugs (NSAIDs) is useful for in vitro diagnosis of the NSAIDs hypersensitivity syndrome. *Clin Exp Allergy* **34**, 1448-57 (2004)
- Uguccioni, M., C. R. Mackay, et al.. "High expression of the chemokine receptor CCR3 in human blood basophils. Role in activation by eotaxin, MCP-4, and other chemokines." *J Clin Invest* **100**(5): 1137-43 (1997)
- Ducrest, S., F. Meier, et al. (2005). "Flowcytometric analysis of basophil counts in human blood and inaccuracy of hematology analyzers." *Allergy* **60**(11): 1446-50.

**Flow CAST®**

**TIME TO RESULT: ~ 1 HOUR**

\* **Note:** Depending on Flow cytometry instrumentation more wash buffer (e.g. 800 µl) might be necessary.

**APPENDIX II**  
**SYMBOLS/SYMBOLE/ SYMBOLES/SIMBOLI/ SIMBOLOS**

| Symbol  | Explanation  |
|---|--|
|  | Use By<br>Verwendbar bis<br>Utiliser jusqu'au<br>Utilizzare entro<br>Fecha de caducidad  |
| <b>REF</b>  | Catalogue number<br>Bestellnummer<br>Réf rence du catalogue<br>Numero di catalogo<br>N mero de cat logo  |
| <b>LOT</b>  | Batch code<br>Chargenbezeichnung<br>Code du lot<br>Codice del lotto<br>Codigo de lote  |
|  | Contains sufficient for <n> tests<br>Ausreichend f r „n“ Ans tze<br>Contenu suffisant pour „n“ tests<br>Contenuto sufficiente per „n“ saggi<br>Contenido suficiente para <n> ensayos |
|  | Consult Instructions for Use-<br>Gebrauchsanweisung beachten<br>Consulter le mode d'emploi<br>Consultare le istruzioni per l'uso<br>Consulte las instrucciones de uso                |
|  | Temperature limitation<br>Zul ssiger Temperaturbereich<br>Limites de temp rature<br>Limiti di temperatura<br>Limite de temperatura   |

| Symbol              | Explanation  |
|---------------------|--|
| <b>BUF STIM</b>     | Stimulation Buffer<br>Stimulations-Puffer<br>Tampon de stimulation<br>tampone di stimolazione<br>Tamp n de estimulaci n                                  |
| <b>CONTROL STIM</b> | Stimulation Control<br>Stimulationskontrolle<br>Contr le de stimulation<br>Controllo di stimolazione<br>Control de estimulaci n                          |
| <b>CONTROL FMLP</b> | Stimulation Control fMLP<br>Stimulationskontrolle fMLP<br>Contr le de stimulation fMLP<br>Controllo di stimolazione fMLP<br>Control de estimulaci n fMLP |
| <b>REAG STAIN</b>   | Staining Reagent<br>F rbe-Reagenz<br>R actif de coloration<br>Reagente di colorazione<br>Reactivo de coloraci n  |
| <b>REAG LYS</b>     | Lysing Reagent<br>Lyse Reagenz<br>R actif de lyse<br>Reagente di lisi<br>Reactivo de l sis   |
| <b>BUF WASH</b>     | Wash Buffer<br>Wasch-Puffer<br>Tampon de lavage<br>tampone di lavaggio<br>Tamp n de lavado   |

Printing Date  
2012-11-30