

# Breast Cancer – The Next Stage

A molecular diagnostic test for quantitative determination of the

**4**

key biomarkers used in the sub typing of breast cancer

ER

PR

HER2

Ki-67





## Molecular information drives treatment choices in breast cancer

Estrogen receptor (ER), progesterone receptor (PR), human epidermal growth factor receptor 2 (HER2) and the marker of proliferation Ki-67 are key biomarkers in the evaluation of breast cancer tumours

**The combination of the biomarker results allows the assessment of the different St. Gallen Breast Cancer subtypes, which are a key parameter for treatment decisions**

Ki-67 is a prognostic and predictive marker. Analytical challenges like high observer variability hinder its standardized and reproducible determination

### Definition of Breast Cancer Surrogate Subtypes (St Gallen 2013)

Breast Cancer Subtypes	ER	PR	HER2	Ki-67
Luminal A-like	Pos	Pos	Neg	Neg
Luminal B-like (HER2 negative)	Pos	Pos/Neg*	Neg	Pos/Neg*
Luminal B-like (HER2 positive)	Pos	Pos/Neg	Pos	Pos/Neg
HER2 positive (non-luminal)	Neg	Neg	Pos	Pos/Neg
Triple negative (ductal)	Neg	Neg	Neg	Pos/Neg

\*with the exception of the combination PR pos and Ki-67 neg = Luminal A-like

### MammaTyper®'s RT-qPCR technology has the following accepted advantages:

- Standardized performance and fast turn-around time
- Minimized inter- and intra-laboratory variability
- Quantitative results with wide dynamic range



## MammaTyper<sup>®</sup> is an easy-to-use test delivering precise results within 6 hours


1



### Sample preparation

10 µm FFPE tissue section (tumor cell content > 20 %).

2



### RNA extraction

Use of RNXtract<sup>®</sup> or validated commercial RNA extraction systems is recommended.

3



### MammaTyper<sup>®</sup> test set up

Preparation of mastermixes and distribution on 96 well plate.  
Analysis of up to 8 patient samples per run.

4



### RT-qPCR analysis

Validated on the following qPCR instruments:

- Roche cobas z<sup>®</sup> 480 Analyzer
- Roche LightCycler<sup>®</sup> 480 II
- Applied Biosystems<sup>®</sup> 7500 Fast (Dx)
- Siemens Versant<sup>®</sup> kPCR Cycler
- Bio-Rad CFX96<sup>®</sup> (IVD, non-deep well)
- Agilent Technologies Mx3000P

5



### Data processing and reporting


Export of mRNA expression data.  
Calculation and assessment of results.  
Results provided within 6 hours.



## MammaTyper® delivers consistent treatment guidance supported by accurate MKI67 determination

Precise Ki-67 evaluation provides prognostic value for patient outcomes

Comparison of Ki-67 expression determined by MammaTyper® and IHC



■ MKI67 negative/ MammaTyper®  
■ MKI67 positive/ MammaTyper®  
■ Ki-67 negative/ IHC  
■ Ki-67 positive/ IHC

Favorable DDFS is independently associated with low MKI67 mRNA expression determined by MammaTyper®

Determination of MKI67 by MammaTyper® delivers extended information on patient's risk of developing distant metastases based on validated cut-off

DDFS = Distant disease free survival

**Clinical outcome proves that Ki-67 determination by MammaTyper® is superior to IHC\***


\*see page 11

Wirtz RM et al., Breast Cancer Res Treat 2016; 157(3), 437-446

## MammaTyper® enhances breast cancer biomarker assessment

MammaTyper® results provide a quantitative measure for each biomarker

Example: Scatter plot shows quantitative MammaTyper® PGR values vs. quantitative IHC PR values



Scatter plot based on in-house data

MammaTyper® mRNA expression results correlate very well to quantitative IHC values.

MammaTyper® cut-offs classify biomarker values into positive or negative results. Dichotomized biomarker results enable stratification into subtypes.

### Precise and reproducible biomarker assessment with MammaTyper®

- MammaTyper® cut-offs are validated based on clinical outcome
- MammaTyper® results express a wider dynamic range than IHC
- Precise and highly reproducible biomarker results through standardized assessment

Wirtz RM et al., Poster San Antonio Breast Cancer Symposium December 9-13, 2014  
 Laible M et al., BMC Cancer 2016; DOI: 10.1186/s12885-016-2476-x





## MammaTyper® provides guidance for treatment decisions of breast cancer patients

### Quantitative gene expression report



Results enable accurate molecular subtyping of tumor tissue according to St Gallen guidelines

MammaTyper® provides quantitative results reflecting the degree of mRNA expression per biomarker

Validated cut-offs are used to determine the positivity or negativity of the biomarker


Values in the ranges colored in light blue/orange have been observed in a cohort of 752 breast cancer samples. Values in the ranges colored in dark blue/orange have not yet been observed but are still valid

## Varying tumor cell content has minimal influence on MammaTyper® performance\*

Fluctuating tumor cell content (TCC) could possibly affect the validity of quantitative assessment of ER, PR, HER2 and Ki-67. Therefore, the performance of MammaTyper® was investigated under different scenarios of TCC.

MammaTyper® results of paired tumor samples with high TCC (> 80 %) derived by macrodissection and low TCC (5 – 50 %) including varying DCIS\* content (10 – 70 %)

Gene nomenclature ERBB2 = HER2, ESR1 = ER, PGR = PR, MKI67 = Ki-67



Despite varying tumor cell content, MammaTyper® results for HER2, ER, PR and Ki-67 revealed a high level of concordance between low TCC and high TCC samples

## Macrodissection of FFPE samples with a tumor cell content > 20 % is not required for gene expression analysis using MammaTyper®

\* shown for HER2 negative tissue samples with TCC > 20 %; DCIS = ductal carcinoma in situ

Wirtz RM et al., Poster San Antonio Breast Cancer Symposium December 9-13, 2014  
Laible M et al., BMC Cancer 2016; DOI: 10.1186/s12885-016-2476-x



## MammaTyper® International Multicenter Study demonstrates excellent reproducibility

### MammaTyper® International Multicenter Reproducibility Study

Evaluated the inter- and intra-site reproducibility of the quantitative detection of ERBB2, ESR1, PGR and MKI67 mRNA expression in clinical samples

- 10 different sites in Europe, North America and Asia
- Standardized application training
- Locally and centrally extracted total RNA of 24 clinical FFPE samples
- Assessment of the precision of the test under different conditions: Laboratories, operators, instruments, days and lots

### MammaTyper® Reproducibility Study results

Excellent inter-site agreement of binary single biomarker classification (positive/negative) is represented by high Kappa values:

Binary biomarker results	ERBB2	ESR1	PGR	MKI67
Kappa values	1.00	0.91	0.94	0.94

Low inter-site variance of quantitative single biomarker results is demonstrated by excellent K values  $\geq 0.98$ :

Quantitative biomarker results	ERBB2	ESR1	PGR	MKI67
Kappa values	0.987	0.992	0.998	0.980

Determination of breast cancer subtypes shows a high level of agreement across all sites (Kappa = 0.90)

**MammaTyper® is highly reproducible – reducing inter- and intra- laboratory variations and allows e.g. standardization of Ki-67 assessment**

## MammaTyper® allows accurate biomarker assessment

### Clinical performance evaluation: MammaTyper® FinHer study

Concordance of MammaTyper® results with IHC/CISH based standard diagnostic methods was evaluated using 769 tissue samples obtained within the FinHer trial.

#### Patients:

Node positive or high risk node negative invasive breast cancer.

#### FinHer trial design:

The trial evaluated the efficacy of combining FEC with Docetaxel vs. Vinorelbine. Patients with HER2 positive tumors were also assigned to receive or not receive Trastuzumab.

### Concordance between MammaTyper® and IHC/CISH-based biomarker assessments

	ESR1 (ER)	PGR (PR)	ERBB2 (HER2)	MKI67 (Ki-67)
<b>Concordance</b>	<b>91.8%</b> 660/719	<b>82.5%</b> 593/719	<b>91.8%</b> 660/719	<b>75.0%</b> 516/688
PPA	<b>95.9%</b> 490/511	<b>93.2%</b> 368/395	<b>85.9%</b> 140/163	<b>89.1%</b> 369/414
NPA	<b>81.7%</b> 170/208	<b>69.4%</b> 225/324	<b>93.5%</b> 520/556	<b>53.7%</b> 147/274
Kappa statistic	0.80 0.75–0.85 p < 0.0001	0.64 0.58–0.70 p < 0.0001	0.77 0.72–0.83 p < 0.0001	0.45 0.38–0.52 p < 0.0001

PPA Positive percent agreement, NPA Negative percent agreement

Assessment of ER, PR and HER2 by MammaTyper® correlated well with results obtained by IHC and CISH

As expected, Ki-67 shows moderate concordance between the IHC and MammaTyper® results due to known technical limitation in standardization of Ki-67 IHC assessment

**MammaTyper® results show a high degree of concordance with IHC/CISH for ER, PR and HER2**






# MammaTyper® results ensure precise determination of breast cancer subtypes correlated with clinical outcome

Accurate MKI67 assessment by MammaTyper® has substantial impact on distinction between Luminal A- and B-like breast cancers.

MammaTyper® luminal subtypes correlated with FinHer\* clinical outcome data, thus proving the accuracy of MammaTyper® results

## Concordance of molecular surrogate subtypes defined by MammaTyper® and IHC (FinHer\* study)



**By precise subtyping, MammaTyper® results support selection of appropriate treatment strategy for each patient**


\* FinHer results reclassified according to St Gallen classification  
Wirtz RM. et al., Breast Cancer Res Treat 2016; 157(3), 437-446

Hypothesis-generating data:

# MammaTyper® predicts benefit from Taxane treatment\* by accurate distinction of molecular subtypes

Patients classified as Luminal B-like (HER2 negative) by MammaTyper® show benefit from Taxane-based chemotherapy.

## Comparison of OS and different treatment regimens of Luminal A-like and Luminal B-like (HER2 negative) patients stratified by MammaTyper® or IHC



Shown for Luminal B-like (HER2 negative) patients  
FinHer results reclassified according to St Gallen classification, OS = Overall survival

Patient stratification based on MammaTyper® shows improved OS of the Luminal B-like (HER2 negative) patients receiving Docetaxel-based treatment compared to those receiving Vinorelbine

By accurate stratification of patients into Luminal A-like and Luminal B-like (HER2 negative) subtypes using MammaTyper®, patient groups that might benefit from Docetaxel-based treatment can be revealed

Classification of the patients in Luminal B-like (HER2 negative) subtype by IHC did not reveal a clear separation in responders and non-responders to a Docetaxel-based treatment

**MammaTyper® opens up new opportunities of providing predictive information about the benefit of adjuvant Taxane-based treatment**

Wirtz RM et al., Breast Cancer Res Treat 2016; 157(3), 437-446

# MammaTyper® precisely determines mRNA expression of ER, PR, HER2 and Ki-67



## High-performance test

### Quantitative RT-qPCR assay (CE marked IVD)

- Highly reproducible biomarker assessment
- Reliable results through standardized biomarker detection
- Accurate stratification of breast cancer tumors into St Gallen subtypes



## Promising clinical utility

### Clinical value validated in numerous performance evaluation studies

- Outperforms IHC by accurate Ki-67 determination
- Provides information on patient's prognosis
- Accurate subtyping supports treatment decisions



## Easy-to-use

- Reliable method for any molecular pathology laboratory
- Validated for multiple qPCR instruments\*
- From resection or core needle biopsy FFPE sample to result within 6 hours

FFPE = formalin fixed paraffin embedded

# Innovation for your breast cancer diagnostics

## High-performance

Precise quantitative results for mRNA expression of HER2, ER, PR and Ki-67

Outperforms IHC by accurate Ki-67 determination

Accurate and reliable test to stratify breast cancer into surrogate subtypes acc. to St Gallen

## Reproducible

Highly reproducible results with extremely low inter-/intra-laboratory variation

## Validated

Clinical value has been validated in numerous performance evaluation studies

## Prognostic

Extended information on patient's prognosis

## Predictive

Accurate subtyping supports treatment decisions

## Optimized

Reliable method for any pathology laboratory

## Easy-to-use

Easy-to-use test which allows results within 6 hours

\* Validated RT-qPCR instruments: Roche cobas z® 480 Analyzer, Roche LightCycler® 480 II, Applied Biosystems® 7500 Fast (Dx), Siemens Versant® kPCR Cycler, Bio-Rad CFX96®, Agilent Technologies Mx3000P

Laible M. et al., BMC Cancer 2016; DOI: 10.1186/s12885-016-2476-x  
 Wirtz RM. et al., Breast Cancer Res Treat 2016; 157(3), 437-446  
 Varga Z. et al., Breast Cancer Research 2017; 19:55. DOI 10.1186/s13058-017-0848-z  
 Sinn HP. et al., BMC Cancer 2017; 17:124

All scientific and technical information in this brochure is based on the MammaTyper® Instructions for Use Rev. 3.1 and Wirtz et al., Breast Cancer Res Treat. 2016; 157(3) 437-446; Laible et al., BMC Cancer 2016; DOI: 10.1186/s12885-016-2476-x; Varga Z et al., Breast Cancer Research 2017; 19:55. DOI 10.1186/s13058-017-0848-z; Sinn et al., BMC Cancer 2017; 17:124



To order MammaTyper<sup>®</sup> or for further information, please contact your local distributor: