

# Quantum Blue® Adalimumab trough level rapid test

adapt therapy  
without delay



- First quantitative rapid test for Adalimumab monitoring
- Time to result 15 minutes
- Measuring Range 1.3-35 µg/mL
- High agreement to reference ELISAs


# Quantum Blue® Adalimumab rapid test

## Key features

- The first quantitative rapid test for adalimumab
- Measuring range 1.3 -35 µg/mL
- Covers the range of the relevant 5 to 12 µg/mL for trough levels\*
- High correlation and agreement to the reference ELISA assays

\* Mitrev et al. 2017

## Correlations to ELISA assay:



— Passing-Bablok fit

		RIDASCREEN ADM MONITORING (ELISA, µg/mL)		
		low	optimal	high
Quantum Blue® Adalimumab (µg/mL)	low	20.0%	5.0%	0.0%
	optimal	0.0%	27.5%	0.0%
	high	0.0%	12.5%	35.0%

total agreement 82.5%



BÜHLMANN Laboratories AG  
Baselstrasse 55  
4124 Schönenbuch  
Switzerland

Phone +41 61 487 12 12  
Fax orders +41 61 487 12 99  
info@buhlmannlabs.ch  
www.buhlmannlabs.ch

Quantum Blue® is a registered trademark of BÜHLMANN in many countries.

**Ordering codes:**  
LF-TLAD10 10 cassettes  
LF-TLAD25 25 cassettes

**Related products:**  
LF-TLIF10 10 cassettes  
LF-TLIF25 25 cassettes  
LF-CAL25 25 cassettes  
LF-CALE25 25 cassettes  
LF-CHR25 25 cassettes

CE-marked products